

**Frederick National Laboratory
for Cancer Research**

sponsored by the National Cancer Institute

Vaccine, Immunity and Cancer Directorate
Standard Operating Procedure

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 1 of 16

Supersedes

2.0

Effective Date: 12Aug21

Author Name	Title	Signature/Date

Approver Name	Title	Signature/Date

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 2 of 16

Supersedes

2.0

Effective Date: 12Aug21

1. PURPOSE

- 1.1. The purpose of this procedure is to describe the proper use and maintenance of the pH Meter.

2. SCOPE

- 2.1. This procedure applies to all pH Meters.

3. REFERENCES

- 3.1. Mettler Toledo SevenCompact pH Meter S220 User Manual
- 3.2. Mettler Toledo Question and Answer on pH Electrodes:
https://www.mt.com/us/en/home/products/Laboratory_Analytics_Browse/pH-meter/sensor/pH-sensor.html#sensors_question_1
- 3.3. 10007: Non-Routine Equipment Maintenance
- 3.4. 10009: General Record Review
- 3.5. 26016: Operation, Use and Maintenance of the Water Purification Systems
- 3.6. HSL_GL_001: Waste Disposal at the Advanced Technology Research Facility

4. RESPONSIBILITIES

- 4.1. The Research Associate, hereafter referred to as Analyst, is responsible for reviewing and following this procedure, and documenting performance of equipment maintenance.
- 4.2. The Scientific Manager or designee is responsible for training personnel in this procedure and reviewing associated documentation.
- 4.3. The Quality Assurance Specialist is responsible for quality oversight and approval of this procedure.
- 4.4. Trained personnel perform equipment maintenance record review per "10009: General Record Review."

5. DEFINITIONS

- 5.1. As Needed Maintenance – maintenance that is performed outside of routine maintenance but is not performed in response to equipment malfunction.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 3 of 16

Supersedes

2.0

Effective Date: 12Aug21

- 5.2. Non-Routine Maintenance – maintenance that is performed in response to equipment malfunction or failure.
- 5.3. Routine Maintenance – maintenance that is performed at planned intervals to identify and prevent problems before they result in equipment failure.
- 5.4. TOC – Total Oxidizable Carbon
- 5.5. Type II Water – Pure/Analytical Grade, used for standard applications (Resistivity >1 MΩ-cm and TOC ≤ 50 ppb)

6. REAGENTS, CHEMICALS, AND EQUIPMENT

- 6.1. 150 mL Beaker (Nalgene, Cat # Z130966 or equivalent)
- 6.2. Calibration Buffers pH 4.0, 7.0, 10.0 (VWR, Cat # 97007-790 kit or equivalent)
- 6.3. Electrode, pH, InLab Routine Pro-ISM (Mettler Toledo, Cat # 51344155)
- 6.4. Electrolyte Solution (VWR, Cat # 97007-764 or equivalent)
- 6.5. pH meter (Mettler Toledo, SevenCompact, Cat # 30019028 or equivalent)
- 6.6. Printer, RS-P25 (Mettler Toledo, Cat # 11124300 or equivalent)
- 6.7. Tubes, Conical, 50mL (Warehouse, Cat # 66401493 or equivalent)
- 6.8. Wash Bottle, Standard Spout (VWR, Cat # 414004-226 or equivalent)
- 6.9. Water, Type II
- 6.10. Wipe, No-Lint, Kimwipes (VWR, Cat # 21905-026 or equivalent)

7. HEALTH AND SAFETY CONSIDERATIONS

- 7.1. Proper safety precautions should be taken while working in a laboratory setting. This includes, but is not limited to, proper protective equipment such as lab coats, safety glasses, closed-toe shoes, and non-latex gloves.
- 7.2. Refer to the respective Safety Data Sheet (SDS) when working with any chemicals.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 4 of 16

Supersedes

2.0

Effective Date: 12Aug21

- 7.3. Refer to “HSL_GL_001: Waste Disposal at the Advanced Technology Research Facility,” “EHS-WM-1: Disposal and Minimization of Chemical Waste,” and “EHS-WM-2: Biological Waste Handling and Disposal” for waste disposal processes.


8. GENERAL USE

- 8.1. Press the On/Off key to turn pH Meter on. Press and hold the On/Off key for 2 seconds to turn pH Meter off after use.
- 8.2. pH Meter MUST be standardized/calibrated daily prior to use with three (3) calibration buffers (pH 4, 7, 10). Make fresh aliquots daily of each buffer into labeled 50 mL conical tubes prior to calibration. Calibration buffers can be used in any order.

Note: If the reagent is acidic, and the slope passes during calibration for the pH 4 and 7 calibrators (B1/B2, $\geq 93\%$), but not the pH 7 and 10 calibrators (B2/B3, $< 93\%$), then the analyst may proceed with the pH of the reagent. However, if the pH 4 and 7 calibrator slope failed (B1/B2, $< 93\%$), then the calibration must be repeated, or pH meter fixed before use. The same scenario holds true for a basic solution. The pH 7 and 10 calibrator slope must pass before use (B2/B3, $\geq 93\%$).

- 8.3. Always rinse the pH electrode with Type II Water before, in-between, and after measurements. May use Wash Bottle with a standard spout when rinsing electrode. Collect liquid run-off in beaker for proper disposal per “HSL_GL_001: Waste Disposal at the Advanced Technology Research Facility.”

Image 1: Rinse Electrode with Type II Water



- 8.4. When blotting electrode dry with Kimwipe, do not dry or touch the bottom/underside of the sensor tip. (See Attachment 1: InLab Routine Pro-ISM pH Sensor)
- 8.5. Place sensor tip as close to the middle of the sample as possible and use a large enough sample container so the sensor tip does not touch the bottom or sides of the container.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 5 of 16

Supersedes

2.0

Effective Date: 12Aug21

Image 2: Placement of Electrode





Image 3: Check Electrolyte Level If Needed For pH Electrode



- 8.7. Verify that the filling hole on the electrode is plugged during storage and use. (See Attachment 1)

Image 4: Ensure Filling Hole is Secured



- 8.8. Verify that the sensor tip is filled with electrolyte solution and does not have any bubbles. If bubbles are present, flick the sensor until bubbles disappears. (See Attachment 1)

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 6 of 16

Supersedes

2.0

Effective Date: 12Aug21

Image 5: Shake Electrode While Holding Top to Remove Bubbles



- 8.9. Do not allow external portion of pH electrode to dry out, which will damage the membrane, nor store in Type II Water, which will deplete the ion rich reference electrolyte. Wetting cap, filled with 3 M KCl solution, is to be placed over the end (tip) of the pH electrode during electrode storage and inactivity.

Image 6: Wetting Cap



- 8.10. The integrated temperature probe measures pH with temperature compensation. See Attachment 4: pH Electrode.

9. EQUIPMENT USE

9.1. Standardization/Calibration

9.1.1. Rinse pH electrode with Type II Water then gently blot dry with Kimwipe.

9.1.2. Place electrode in pH 4.0 calibration buffer/standard and press **CAL**.

⇒ **Cal 1** appears on the display. The meter endpoints according to preselected endpoint mode after signal has stabilized or after pressing **READ**. The relevant buffer value is shown on the display.

9.1.3. Rinse electrode with Type II Water then gently blot dry with Kimwipe.

9.1.4. Place electrode in pH 7.0 calibration buffer and press **CAL**.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 7 of 16

Supersedes

2.0

Effective Date: 12Aug21



Cal 2 appears on the display.

9.1.5. Rinse electrode with Type II Water then gently blot dry with Kimwipe.

9.1.6. Place electrode in pH 10.0 calibration buffer and press **CAL**.



Cal 3 appears on the display.

9.1.7. Rinse electrode with Type II Water then gently blot dry with Kimwipe.

9.1.8. Press **Save** to keep calibration or **EXIT** to reject calibration.

Note: Calibration data will automatically print when pressing **Save**. If it does not print automatically, follow data transfer procedure in section 10.

9.1.9. Record calibration information on "26011-01: pH Meter Daily Use Calibration Form" and attach printout.

9.1.10. Verify the B1/B2 slope is $\geq 93\%$ and B2/B3 slope is $\geq 93\%$; otherwise, repeat calibration with new pH calibration fluids. If after the second calibration event the slopes fail, regenerate electrode with Electrolyte solution (3 M KCl) and repeat calibration.

9.2. Sample Measurements

9.2.1. Rinse electrode with Type II Water then gently blot dry with Kimwipe.

9.2.2. Place sensor in the sample (see Image 2 for placement in container) and press **READ** to start a measurement.

9.2.3. The display shows the readings of the sample. The endpoint format blinks, indicating a measurement is in progress.

9.2.4. As soon as measurement is stable according to the selected stability criterion (Automatic Endpoint), the **Stability** icon appears.

Note: "Automatic endpoint" format is selected, measurement stops automatically once the **Stability** icon appears.

9.2.5. Data automatically prints after measurement is taken. If not, follow data transfer procedure in section 10.

10. DATA TRANSFER PROCEDURE

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.

This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 8 of 16

Supersedes

2.0

Effective Date: 12Aug21

10.1. Connect printer cable to pH Meter and the corresponding interface on backside of printer.

Note: Ensure printer's green light is not flashing. A steady green light indicates the printer is synced and ready.

10.2. Select interface "**printer**" in data transfer settings menu.

10.3. Start transfer in data menu. The data will automatically print.

11. MAINTENANCE

11.1. Monthly Maintenance

11.1.1. Verify pH electrode is always completely covered by appropriate filling solution (3 M KCl) using the wetting cap and do not allow electrode to dry out.

11.1.2. Verify pH electrode is always kept filled with appropriate filling solution (3 M KCl).

11.1.3. Verify that the filling hole on the electrode is plugged during storage and use. (See Attachment 1)

11.2. As Needed Maintenance

11.2.1. Meter Maintenance

11.2.1.1. Spray cavicide on a low-lint towel and wipe exterior surfaces.

11.2.2. External Electrode Maintenance

11.2.2.1. Ensure pH electrode is always completely covered by appropriate filling solution (3 M KCl) and do not allow electrode to dry out.

11.2.2.2. Rinse pH electrode using Type II Water as needed.

11.2.3. Spills

11.2.3.1. Turn pH Meter off.

11.2.3.2. Spray cavicide on a low-lint towel and wipe exterior surfaces.

11.2.3.3. Use low-lint absorbent wipe to dry components.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 9 of 16

Supersedes

2.0

Effective Date: 12Aug21

11.2.4. Document As Needed Maintenance in its respective section on form "26011-02 pH Maintenance Form."

11.3. Non-Routine Maintenance

11.3.1. In the case that the pH Meter is not operating correctly, transition processes being performed to another unit (when applicable), post a sign stating the equipment is out of service and initiate non-routine maintenance documentation per "10007: Non-Routine Equipment Maintenance."

11.3.2. Document the nature of any failures or malfunctions, how and when it was discovered, and the personnel involved on "10007-01: Non-Routine Equipment Maintenance Form."

11.3.3. Initiate a service request and complete the non-routine maintenance process following 10007.

12. SETTINGS

12.1. pH stability reading is set at "Automatic Endpoint."

13. ATTACHMENTS

13.1. Attachment 1: InLab Routine Pro-ISM pH Sensor

13.2. Attachment 2: 26011-01: pH Meter Daily Use Calibration Form

13.3. Attachment 3: How to Detect When a pH Sensor Needs to Be Replaced

13.4. Attachment 4: pH Electrode

13.5. Attachment 5: pH Maintenance Form

**Frederick National Laboratory
for Cancer Research**

sponsored by the National Cancer Institute

Vaccine, Immunity and Cancer Directorate
Standard Operating Procedure

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 10 of 16

Supersedes

2.0

Effective Date: 12Aug21

14. REVISION HISTORY

Revision Start Date	Version #	Changes	Reasons
15Mar17	New	Create new SOP for Use and Maintenance of a pH meter	Currently no SOP
25Mar19	2.0	Updated Maintenance Procedures	Reflect current GDP practices
04Aug21	3.0	Added Images and edited procedure	Clarity

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0


Page 11 of 16

Supersedes

2.0

Effective Date: 12Aug21

Attachment 1: InLab Routine Pro-ISM pH Sensor



Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

**Frederick National Laboratory
for Cancer Research**

sponsored by the National Cancer Institute

**Vaccine, Immunity and Cancer Directorate
Standard Operating Procedure**

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 12 of 16

Supersedes

2.0

Effective Date: 12Aug21

Attachment 2: 26011-01: pH Meter Daily Use Calibration Form

Frederick National Laboratory for Cancer Research <small>sponsored by the National Cancer Institute</small>		HPV Serology Laboratory Standard Operating Procedure Form	
Form Title: pH Meter Daily Use Calibration Form			
Document ID: 26011-01		Version:	3.0
Associated SOP: 26011		Effective Date:	12Aug21
Supersedes Version:	2.0	Page 1 of 1	

Equipment ID:	HSL_					
Calibration Buffer pH	Part Number	Lot Number	Expiration Date	B1/B2 Slope (%) (≥93%)	B2/B3 Slope (%) (≥93%)	Comments:
<input type="checkbox"/> N/A 4						<input type="checkbox"/> Pass <input type="checkbox"/> Fail <input type="checkbox"/> Other:
7						
<input type="checkbox"/> N/A 10						
Print out: <input type="checkbox"/> N/A						

Performed by/date: _____

Reviewed by/date: _____

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.

This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0













Page 13 of 16

Supersedes

2.0

Effective Date: 12Aug21

Attachment 3: How to Detect When a pH Sensor Needs to Be Replaced.

Slope (% / mV)	Offset (mV)		
	$\pm 0 - 20$	$\pm 20 - 35$	$> \pm 35$
95-105 (56.20-62.12)			
90-95 (53.24-56.20)			
85-90 (50.29-53.24)			
<85 or >105 (<50.29 or >62.12)			

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0


Page 14 of 16

Supersedes

2.0

Effective Date: 12Aug21

Attachment 4: pH Electrode



Typical combination pH electrode with inner pH sensor and outer reference element

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

Frederick National Laboratory for Cancer Research

sponsored by the National Cancer Institute

Vaccine, Immunity and Cancer Directorate Standard Operating Procedure

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 15 of 16

Supersedes

2.0

Effective Date: 12Aug21

Attachment 5: pH Maintenance Form

Frederick National Laboratory for Cancer Research <small>sponsored by the National Cancer Institute</small>		HPV Serology Laboratory Standard Operating Procedure Form	
Form Title: pH Maintenance Form			
Document ID: 26011-02		Version:	3.0
Associated SOP: 26011		Effective Date:	12Aug21
Supersedes Version:	2.0	Page 1 of 2	

Maintenance Year:		Equipment ID:	
--------------------------	--	----------------------	--

Monthly Maintenance						
Month	January	February	March	April	May	June
Verify Wetting Cap filled and placed on electrode	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl
Verify electrode filled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl
Verify filling hole is plugged during storage	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES
Performed By/Date:						
Reviewed By/Date:						
Month	July	August	September	October	November	December
Verify Wetting Cap filled and placed on electrode	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl
Verify electrode filled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl	<input type="checkbox"/> YES <input type="checkbox"/> Refilled with 3 M KCl
Verify filling hole is plugged during storage	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES	<input type="checkbox"/> YES
Performed By/Date:						
Reviewed By/Date:						

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.

This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

**Frederick National Laboratory
for Cancer Research**

sponsored by the National Cancer Institute

Vaccine, Immunity and Cancer Directorate
Standard Operating Procedure

SOP Title: Use and Maintenance of a pH Meter

Document ID: 26011

Version

3.0

Page 16 of 16

Supersedes

2.0

Effective Date: 12Aug21

**Frederick National Laboratory
for Cancer Research**
sponsored by the National Cancer Institute

HPV Serology Laboratory
Standard Operating Procedure
Form

Form Title: pH Maintenance Form

Document ID: 26011-02

Version:

3.0

Associated SOP: 26011

Effective Date:

12Aug21

Supersedes Version:

2.0

Page 2 of 2

As Needed Maintenance: ☐ N/A

Date	Activity Performed	Recorded by/date	Reviewed by/date
<input type="checkbox"/> N/A			
<input type="checkbox"/> N/A			

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.

This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.

Verify current version prior to use. Use of a superseded or obsolete document is prohibited.
This document contains confidential and proprietary information. Do not copy or distribute without prior, written permission.