National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 1 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
Author/Date:			
Approvals/Date:_			
Development Enginee BDP	r IV Director, Business Operations BDP	Director BQA	
SOP References: 11115, 11117		Supersedes: N	Α
<u>Purpose</u> : The objective of this procedure is to describe the step-by-step operations to program change values and setpoints for the operation of VHP® 1000ED-AB Generator.			rations to I000ED-AB
<u>Scope</u> : This procedure covers the VHP® 1000ED-AB Generator (BDP MEF 76720, Model number EDA-120, serial number 0126203-05) installed for the Biodecontamination of cGMP Virus Production in the BDP.			
<u>Contents</u> :			
1.0 Authority and Respo	nsibility		
2.0 Definitions			
3.0 Equipment			
4.0 Description of the VH	IP Generator Controls		
5.0 Controls Overview			
6.0 Controls Procedure			
7.0 Sequence of Operation for the VHP Generator			
8.0 Reference			

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 2 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-AB Generator			

1.0 Authority and Responsibility

- 1.1 The Director, Business Operations, Biopharmaceutical Development Program (BDP) has the authority to define this procedure.
- 1.2 The Development Engineer, BDP is responsible for training personnel in this procedure and for documenting this training to Biopharmaceutical Quality Assurance (BQA).
- 1.3 BDP personnel are responsible for implementation and performing this procedure.
- 1.4 BQA is responsible for quality oversight of this operation.

2.0 Definitions

- 2.1 Setpoints: The setpoints are used to control the general way the VHP® 1000ED-AB Generator operates. Unit Setpoints include time/date, out-of-cycle values, and setup values.
- 2.2 Cycle Parameters: The cycle parameters are used to control the decontamination cycle performing by a VHP® 1000ED-AB Generator.
- 2.3 VHP: Vapor Hydrogen Peroxide.

3.0 Equipment

3.1 Steris VHP® 1000ED-AB Generator.

4.0 Description of the Generator Controls

The VHP® 1000ED-AB Generator (henceforth will be referred to as VHP Generator) has two power switches, a main power disconnect switch, and a power ON/OFF switch. The control panel is operated as follows:

- 4.1 The main power disconnect switch is located on the back of the unit. This switch controls power supply to the generator and its control system.
- 4.2 The power ON/OFF switch is located on the back of the unit. This switch supplies power to the generator and its controls.
- 4.3 The operator can control cycle operation, program cycle, and unit operating parameters, as well as monitor cycle performance through the control panel.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 3 of 43	Revision 00

Title: Operating the VHP® 1000ED-AB Generator

- 4.4 Use of the PV300 display is normally self-explanatory. Many screens include the text "F4-Next/F8-Previous" to indicate that additional screens are available.
- 4.5 The following table identifies the special functions given to particular touch pads located on the PV300 display.

	Touch Pad Special Functions
Touchpad	Special Function
F1	Start the selected cycle (from Standby Mode)
F2	Reset from Cycle Complete phase back to Main Menu
F4	Displays the previous screen
F5	Abort cycle if a cycle is active
F6	Force a status print (demand print) from the Run Screen
F8	Displays the next screen when available
▲ ►	Move the cursor to the left/right or navigate through text
	selection screens
$\blacktriangle \blacksquare$	Select from a screen list or navigate through text selection
	screens
	Enter typed value at the cursor
	Clears the last character entered

- 4.6 The printer records the status of VHP 1000ED-AB Generator Power up, cycle parameters, and actual cycle data.
- **NOTE**: Refer Operator Manual (129383-120) page numbers 5-2 to 5-8 for sample printouts.

5.0 Controls Overview

The VHP Generator is provided with user-programmable cycles and sepoints, a complete operator interface, and a printout of generator cycles and critical data. Security for operation, system, and cycle setpoints access and service mode functions is provided with password requirements. The VHP Generator is operated by using one of the following modes:

- 5.1 Operator Mode enables selection of cycle and starting and aborting of cycles.
- 5.2 Supervisor Mode enables setting of cycle setpoints and system setpoints and options.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 4 of 43	Revision 00

Title: Operating the VHP® 1000ED-AB Generator

- 5.3 Service Mode enables I/O testing, device calibration functions, and other service-related functions.
- 5.4 Refer to the following flow charts for the basic structure of the VHP Generator screen sequence.







- 6.0 Controls Procedure
 - 6.1 Position the Power switch to ON. The display will advance to the Main Menu screen and the following printout will occur.



National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 8 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-AB Generator			

6.2 At the Main Menu Screen (shown below), the display prompts for selection of modes (Operator, Supervisor, or Service). Note that cycles may be initiated from the External Interface from this screen.

(date)	(time)
MAIN MENU F1-Operator F5-Service	PV01 F2-Sprvsr

- F1 = Enter Operator Mode
- F2 = Enter Supervisor Mode.
- F5 = Enter Service Mode.
- 6.3 Password Screen: Password screen appears when F1, F2, or F5 is pressed at the Main Menu.

Type in the Password and press <ENTER> xx

The Operator Mode requires entry of any of the programmable operator passwords, the supervisor password, or service password. The Supervisor Mode requires entry of the supervisor or service password. The Service Mode requires entry of the service password.

All Passwords are defaulted to the respective user number. Refer to the following table.

Password Defaults			
User	User	Accessible Medee	Default
Number	Туре	Accessible Modes	Password
1	Operator	Operator	1
2	Operator	Operator	2
3	Operator	Operator	3
4	Operator	Operator	4
5	Operator	Operator	5
6	Operator	Operator	6
7	Operator	Operator	7
8	Operator	Operator	8
9	Operator	Operator	9
10	Operator	Operator	10

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 9 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-AB Generator			

- 6.4 The Operator Control Sequence
 - 6.4.1 The Operator Mode is accessed from any valid password entry made at the Main Menu. The first screen, Operator Screen, displayed on the OP7 Display enables selection of the cycle to be started (by the operator), operation of the Reservoir Fill Cycle, the Run Screen, or initiation of a Regeneration Cycle.

Operator Scr	een
F1-Cycle Selection	
F2-H ₂ O ₂ Res. Fill	
F3-Run Screen	
F5-Regeneration Cycle	
F4-Previous	F8-Next

Press F8 and this screen is displayed:

Operator Screen
F1-Print Report
F2-Print Parameter
F4-Previous

NOTE: The alphanumeric label (i.e. F1) corresponds to the controller keypad for selecting the task as displayed on the screen. The printouts listed above may be generated even if the system printer option is disabled. The regeneration Cycle will be initiated immediately when F5 is pressed.

6.4.2 Press F1 to select a cycle. The Cycle Selection Screen will appear. Choose a cycle.

CYCLE SI	ELECTION nn
F1- #1-XXXXXX	
F2-#2-XXXXXX	
F3-#3-XXXXXX	
-4- Previous	F8-Next

Where: nn = current selected cycle XXXXXX = cycle name for the respective cycle Press F8 will display all the available cycle screens.

	Τ		
National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 10 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
6.4.3 Once the F	e a cycle has been selected, the F Run Screen	V300 Display will	advance to
	Cycle- #nn Time Rem.: 00:00:00 Phase: Standby F1-START F4- Previous F8-Next		
Where: r	nn = selected cycle, set per the C	ycle Selection Scr	reen
6.4.4 At an scree corre Typical s	by time while the Run Screen is a ens indicate VHP Generator actua sponding setpoints (S), if applical setpoint readings are shown in the	ctive, the scrolled al readings (A) and ole e displays below.	-down d their
Sc S: Inj S: F4	cale aaaa A: bbbb g ection Rate cccc A: dddd g/m -Previous F8-Next		
Where: a cycle bbbb cccc dddd <u>in gra</u> Press F8	aaaa = Reservior required weight , in grams. = Actual Scale reading, in grams = Setpoint Injection Rate, in gram = Actual Injection Rate (output ra ams/minute (updated every 15 sec 3 and this screen is displayed:	setpoint for the se n/minute. nte from the Reser conds)	elected rvoir),
Er S: Cy S: F4	ncl. Pressure (yyy) aaaa A: bbbb vcle Airflow (xxx) cccc A: dddd -Previous F8-Next		

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 11 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
Where: a bbbb cccc dddd Xxx = <u>yyy =</u> Press F8	aaaa = Setpoint Enclosure pressu = Actual Enclosure pressure. = Setpoint Cycle airflow. = Actual Cycle airflow. = cfm or cmh, depending upon air "WC or Pa depending upon press and this screen is displayed:	re. flow unit set sure unit set.	
Va S: Pre S: F4	por. Temp. aaaa A: bbbb xx e-Heater Temp. cccc A: dddd xx -Previous F8-Next		
Where: a bbbb cccc dddd <u>xx = (</u> Press F8	aaa = Setpoint Vaporizer temper = Actual Vaporizer temperature. = Setpoint Preheater temperature = Actual Preheater temperature C or F depending upon temperatu	ature. e ire unit set.	
Hu S: En A: F4	midity (%RH) aaa A: bbb closure Humidity ccc %RH -Previous F8-Next		
Where: a bbb = 	aa = Setpoint Air Inlet Relative H Actual Air Inlet Relative Humidity Actual Enclosure Relative Humic	umidity. y. lity.	
Press F8	and this last screen is displayed	:	

National Cancer Institute-Frede Frederick, MD	erick,	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Prog	gram		Page 12 of 43	Revision 00
Title: Operating the VHP [®] 100	0ED-A	B Generator		
6.4.5	$\begin{bmatrix} Dr, \\ S: \\ H_2 \\ S: \\ F4 \\ bbbb \\ cccc \\ dddd \\ xx = 0 \\ Press \\ Wher \\ and t \\ F1 \\ F5 \\ F4 \\ F5 \\ F4 \end{bmatrix}$	yer Temp. aaaa A: bbbb xx O_2 Return Temp. cccc A: dddd xx -Previous = Actual Dryer temperature. = Setpoint H ₂ O ₂ Return temperatur C or F depending upon temperatur C or F depending upon temperatur s F4 to return to previous screen. h back to the Operator Screen, se he following screen(s) is displayed H ₂ O ₂ Reservoir Fill -Start S: aaaa -Stop A: bbbb	ure. e. ire unit set. lect "F2-H22 Res d:	servoir Fill"
	Wher bbbb	e: aaaa = Setpoint in gram. = Actual in gram.	I	
6.4.6	On th scree Tir Ph F- F4 Press displa	e Operator Screen, select "F3-Ru en occur. rcle-#nn ne Rem.: hh:mm:ss ase: PHASE NAME Start -Previous F8-Next s F8 and all the series of screens ayed.	n Screen" and th	e following .4.4 are

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 13 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
6.4.7 On Op reger displa Tir Ph	perator Screen, select "F5-Regene neration cycle will be started imme ayed. rcle-RGN ne Rem.: hh:mm:ss pase: PHASE NAME	eration Cycle" and ediately with this s	d the screen
F-3	Start F5-Abort -Next		
6.5 Supervisor Op	perating Sequence		
6.5.1 Supe servic enab setpo F1	rvisor Mode: supervisor mode is a ce password entry made at the Ma les setting of the system setpoints ints. The display shows as follow -System Setup	accessed from su ain Menu. Superv s and options and s.	ipervisor or visor Mode I the cycle
F2 F3 F4	-Cycle Setup -Phase Advance -Previous		
F1 = Set F2 = Set F3 = If a	system setpoints, options, time a Cycle setpoints. cycle is active, advance to the ne	nd date, and pas ext phase.	swords.
6.5.2 Syste follow	em Setup (F1): For System Setup /s	, the display shov	vs as
F1 F2 F3	SYSTEM SETUP -Name F6-Print -Setpoints -Options		
F4	-Previous F8-Next		
F1 = Set F2 = Set F3 = Set F6 = Prir F8 = Scr	Generator system name (printed system setpoints. system configuration options. nt system setpoints and options. oll down screen.	at power-up).	

Press F8 and the this screen is displayed:

National Cancer Institute-Frederi Frederick, MD	ick,	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Progra	am		Page 14 of 43	Revision 00
Title: Operating the VHP [®] 1000	ED-A	B Generator		
F1 = F2 = F3 = 4.1.1 T 4 T G	F1 F2 F3 F4 Set Set Set fhe fo .1.1. he s Sene T etc a .1.1. he s ycles etpo cree alue	SYSTEM SETUP -Time/Date -Passwords -Utility Shutdown -Previous F8-Next time and date. operator and supervisor password oblowing points describe each sel 1 Set System Name (F1) system name is printed at the powrator. SYSTEM NAME (Sys Name 1) (Sys Name 2) the user may set two lines of 20 c ntered must be limited to those o o 7F hex (generally upper and low nd basic typographic symbols). 2 Set System Setpoints (F2) system setpoints apply to the VHF s. Pressing F8 scrolls the screen ints. Pressing F4 scrolls the screen st shown on the following screens	rds. vn and restart. ection. ver-up print of the characters. Chara f the ASCII set fr ver case letters, n P Generator, not (down) to additic en back (up) to p s are the typical of	e VHP

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 15 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
This s during 999.9 NOTI Deco	SYSTEM SETPOINTS Reservoir Fill 600 grams F4-Previous F8-N setpoint (above) is the weight to v g the H ₂ O ₂ Reservoir Fill Cycle ar grams. E: That this applies only to the stantamination Cycle, the reservoir i	ext which the reservo nd may be set fro and-alone fill cycle s filled to the nec	ir is filled m 0 to e. During a essary
This s be inj	SYSTEM SETPOINTS Dryer Capacity 900 grams F4-Previous F8-N setpoint (above) is the amount of ected before the dryer must be re 0.0 to 999.9 grams.	ext hydrogen peroxic egenerated and n	de that can nay be set
This s	SYSTEM SETPOINTS Vap. Temp. (Out Cycle) 25.0°C F4-Previous F8-N setpoint (above) is the temperaturer is maintained when the VHP G	ext re at which the va	aporizer -cvcle. and
may l appro heat l <u>for pr</u>	be set from 0-200.0°C/°F. This va oximately room temperature to co buildup within the VHP Generator oper vaporization of H ₂ O ₂ during	lue is normally so nserve energy ar This heater is r injection phases.	et to nd to reduce necessary



National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 17 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
This s the ni auton value NOTI	SYSTEM SETPOINTS Cycle Count 0 F4-Previous F8-N setpoint (above) is the VHP Gene umber of cycles run, and may be natically tracked by the VHP Gen of 1,000,000 before being reset of 1,000,000 before being reset E: That all started cycles are cons including aborted cycles.	erator running cyc set from 0-99999 erator and has a to "0". sidered in this cor	cle count, or 99. This is maximum unt,
	Flexible PID Params Prp. Gain: (Int. Time: 800 Der. Time: F4-Previous F8	0.0110 FLEX 70 8-Next	
These Deriv Flexit	e setpoints (above) are the Propo ative setpoints for the enclosure p ple Enclosure Type.	ortional, Integral, pressure PID loo	and p of the
These	SYSTEM SETPOIN Semi Rigid PID Params Prp. Gain: Int. Time: 800 Der. Time: F4-Previous F e setpoints (above) are the Propo ative setpoints for the enclosure of	NTS 0.0075 0 SEMI 70 F8-Next ortional, Integral, 5 pressure PID loop	and p of the
Semi	-Rigid Enclosure Type.		

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116				
Biopharmaceutical Development Program		Page 18 of 43	Revision 00				
Title: Operating the VHP [®] 1000ED-A	Title: Operating the VHP [®] 1000ED-AB Generator						
These Deriv Rigid 6.5.3. The S Gene setpo scree defau press the rig	SYSTEM SETPOINTS Rigid PID Params Prp. Gain: () Int. Time: 800 Der. Time: 800 F4-Previous F4 e setpoints (above) are the Properative setpoints for the enclosure percentive setpoints for the enclosure percentive setpoints for the enclosure percent set setpoints for the screen set setpoints. Pressing F4 scrolls the screen set settings shown on the screer set values. Pressing F1 enables the ing F5 disables the option. The certains is percent set set of the screen (i.e., "ENABLE SYSTEM OPTIONS System Printer F1-Enable ENABLE SYSTEM OPTIONS System Printer F1-Enable ENABLE K1-Disable F4-Previous System Printer option enables or option interface Module. Disabling term and power-up prints, or any othe natically by the VHP Generator. All the memory of VHP Generator of the Operator Mode after the cycle	20.0035 RIGD 20 8-Next ortional, Integral, a pressure PID loop options of the VH een (down) to add en back (up) to ph is below are the t e option shown, v urrent setting is s BLED BLED F8-Next disables the print his option disable r printouts genera any forced prints a cun cycle is alway controller and may e is completed.	and b of the P ditional revious typical while shown on LED"). er at the es only ated are always rs retained y be printed				

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 19 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
	SYSTEM OPTIONS Host PLC Interface F1-Enable DISAE K1-Disable F4-Previous	BLED F8-Next	
This Interf will s Host Interf inforr	option (above) enables or disable face to the Host PLC. If this option tart the VHP Generator cycles. If PLC is not able to start a cycle. F face procedure (STERIS Part Nur mation.	s the cycle start E n is enabled, the H this option is disa Refer to the Extern nber P129383-16	External Host PLC bled, the nal 1) for more
	SYSTEM OPTIONS Engineering Mode F1-Enable DIS/ F5-Disable F4-Previous	ABLE F8-Next	
This or op Gene setpo press enab	option (above) enables or disable tions intended for validation or ve erator. If this option is enabled, all bints may be changed during a cy sing ESC from the Run Screen). T led during normal cycle operation	s use of additiona rification use of th the parameters a cle run (permitted his option should	al features ne VHP and cycle l by l not be
	SYSTEM OPTIONS Airflow Units F1-cmh cmh F5-cfm F4-Previous F8-	Next	
This airflo and u chan	option (above) determines the dis w readings/calculations. Unit "cml unit "cfm" is cubic feet per minute. ged after airflow calibration is don	played/printed ur n" is cubic meters This option may e.	nits of per hour, be
This airflo and u chan	option (above) determines the dis w readings/calculations. Unit "cml unit "cfm" is cubic feet per minute. ged after airflow calibration is don	played/printed ur n" is cubic meters This option may le.	hits of per hour, be

National Cancer Institute-Frederi Frederick, MD	SK, STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Progra	n	Page 20 of 43	Revision 00
Title: Operating the VHP [®] 1000	D-AB Generator	•	
T pi "F (c te du te	SYSTEM OPTIONS Pressure Units F1-Pa Pa F5-"wc F4-Previous F4-Previous F8 nis option (above) determines the dis essure readings. Unit ""wc" is inches a" is Pascals. This option may be ch ptional enclosure pressure) calibrati SYSTEM OPTIONS Temperature Units F1-°C °C F5-°F F4-Previous °C soption (above) determines the dis egrees Fahrenheit. This option may mperature calibration is done.	3-Next splayed/printed ur s of water column hanged after press on is done 3-Next splayed/printed ur rees Celsius, and be changed after	hits of , and unit sure hits of I unit "°F" is
T E is 6. T ai pi	SYSTEM OPTIONS Volume Units F1-cfF1-cfcfF5-cmF4-PreviousF4-PreviousF8nclosure Volume readings. Unit "cf" iCubic Meters.5.3.4Print System Configurationnis printout documents system setpothe Operator Interface. This print isintout.	-Next splayed/printed ur s Cubic Feet, and (F6) ints and options t also provided in e	hits of d unit "cm" o the printer each cycle

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 21 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator	I	
6.5.3.	5 Set Time/Date (F1)		
The ti (hour monti date,	ime of day is set and operated in s:minutes:seconds) and the date h/day/year (m-d-y). After entering press F3 to set the values. Time hh:mm:ss (MIL) xx:xx:xx Date mm-dd-yyyy Fi M-D-Y xx-xx-xxxx S	military units is set and operat in the desired tin 3 ET	ed in format ne and/or
Passy set as corres only. Supe additi	words may be set for multiple ope s 4-digit numbers from 0-9999. Th spond to users 1-10 and allow en The supervisor password allows rvisor Mode. Press F4 / F8 to scr onal passwords.	erators and the sume operator passivitry to the Operati entry to Operator oll up/down for er	upervisor, words ng Mode Mode and htry of
	OPERATOR PASSWORD 1-xxxx 4-x 2-xxxx 5-xx 3-xxxx 6-xx F4-Previous OPERATOR PASSWORD 7-xxxx 9-xx	DS xxx xx xx F8-Next DS xxx	
	8-xxxx 10-> F4-Previous OPERATOR PASSW0	RDS	
	11-xxxx F4-Previous	F8-Next	
6.5.3	.7 Set Utility Shutdown (F3)		
Utility week turn c will st	Shutdown and Restart times ma . If the Utility Shutdown option is off from the shutdown time to the cart. Press F4/F8 to scroll up/dow	y be set for each Enabled, the vap restart time and r n for entry of add	day of the orizer will egeneration itional

setpoints.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116		
Biopharmaceutical Development Program		Page 22 of 43	Revision 00		
Title: Operating the VHP [®] 1000ED-A	AB Generator	_			
Scrol and F being	UTILITY SHUTDOWN Utility Shutdown F1-Enable DISA F5-Disable F4-Previous I Down. Use the numeric keypad Restart time. Time appears on the gentered. UTILITY SHUTDOWN Monday Shutdown hh:mm Restart hh:mm F4-Previous	ABLED F8-Next to enter the Utility e PV300 Display	/ Shutdown as it is		
Scrol each Durin scree	Scroll down to enter the Utility Shutdown and Restart times for each day of the week. During Utility Shutdown, the PV300 Display will show the following screens.				
	Time Rem.: hh:mm:ss Phase: Regen HeatUP F5-/	ABORT			
vvnie	Cycle-UTL Time Rem.: hh:mm:ss Phase: Regen CoolDown F5-/	ABORT			
While	e in Regeneration Cool Down. Utility Shutdown Restart Time Monday: hh:mm F5-ABORT n Regeneration is complete.				

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 23 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
6.5.3	.8 Cycle Setup Selection (F2)		
The C cycle availa	Cycle Setup Selection (F2) screer to be set up or edited. Thirteen s able.	n enables selection electable cycles a	n of the re
Γ	CYCLE SELECTION		
	F1-#1-xxxxx		
	F2-#2-xxxxx		
	F3-#3-xxxxx		
	r4-rievious	ro-inext	
	 F1 = Edit Cycle #1 setpoints. F2 = Edit Cycle #2 setpoints. F3 = Edit Cycle #3 setpoints. F8 = Scrool down screen. 	r cycle.	
Scro butto	oll-down screens allow for selections F1-F3:	on of additional cyc	cles, using
	CYCLE SELECTION nn F1-#4-xxxxxx F2-#5-xxxxxx F3-#6-xxxxxx	CYCLE SELEC F1-#13-xxxxxx F2-#14-xxxxxx F3-#15-xxxxxx	CTION nn
	F4-Previous F8-Next	F4-Previous	F8-Next
	CYCLE SELECTION nn F1-#7-xxxxxx F2-#8-xxxxxx F3-#9-xxxxxx	CYCLE SELEC F1-#16-xxxxx F2-#17-xxxxx F3-#18-xxxxx	CTION nn
	F4-Previous F8-Next	F4-Previous	F8-Next

B Generator CYCLE SELECTION nn	Page 24 of 43	Revision 00
B Generator CYCLE SELECTION nn		
CYCLE SELECTION nn		
F3-#10-XXXXX F2-#11-xxxxxx F3-#12-xxxxxx F4-Previous F8-Next F4-Previous F8-Next e: xxxxxx=cycle name for the response a cycle selection is made, the conscreen. 9 Cycle Setup setpoints are set from this screen a setpoints are set from this screen b Cycle Setup setpoints are set from this screen a cycle Name F2-Print Cycle F3-Leak Test F5-Dehumidify F4-Previous F8-Next Set cycle name (printed at the cycle Print cycle setpoints. Set cycle Leak Test setpoints. Set cycle Dehumidify phase setpo Scroll down screen. Scroll up screen Cycle Number Selected.	pective cycle. ntrol advances to n. Each of the ave	o the Cycle ailable er).
	2-#11-xxxxx 3-#12-xxxxx 4-Previous F8-Next : xxxxx=cycle name for the rest a cycle selection is made, the conscreen. Cycle Setup setpoints are set from this screen is set up in the same way. CYCLE SETUP ## F1-Cycle Name F2-Print Cycle F3-Leak Test F5-Dehumidify 4-Previous F8-Next et cycle name (printed at the cyclerint cycle setpoints. et cycle Leak Test setpoints. et cycle Dehumidify phase setports. et cycle Dehumidify phase setports. et cycle Name f2-Print Cycle SETUP ## F1-Cycle Name F2-Print Cycle F3-Leak Test F5-Dehumidify 4-Previous F8-Next et cycle name (printed at the cyclerint cycle setpoints. et cycle Dehumidify phase setports. et cycle Dehumidify phase setports. et cycle Number Selected.	2-#11-xxxxx 3-#12-xxxxx 4-Previous F8-Next : xxxxx=cycle name for the respective cycle. a cycle selection is made, the control advances to screen. Cycle Setup setpoints are set from this screen. Each of the av- is set up in the same way. CYCLE SETUP ## F1-Cycle Name F2-Print Cycle F3-Leak Test F5-Dehumidify 4-Previous F8-Next et cycle name (printed at the cycle printout header rint cycle setpoints. et cycle Leak Test setpoints. et cycle Dehumidify phase setpoints. croll down screen. croll up screen ycle Number Selected.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 25 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
The S	Scroll down screen shows the follo	owing.	
	CYCLE SETUP ## F1-Conditioning F2-Decontaminate F3-Aeration F5-Pressure Control		
	F4-Previous F8-Next		
F1 = F2 = F3 = F5 = F8 = F4 = ## =	Set cycle Condition phase setpoil Set cycle Decontaminate phase s Set cycle Aeration phase setpoin Set cycle Pressure Control phase Scroll down screen. Scroll up screen Cycle Number Selected.	nts. setpoints. ts. e setpoints.	
The S	Scroll down screen shows the follo	owing.	
	CYCLE SETUP ## F1-Enclosure Control F2-I/O Control F3-Vaporizer Setup F5-Preheater Setup		
	F4-Previous F8-Next		
F1 = F2 = F3 = F5 = F4 = ## =	Set cycle Enclosure Control phas Enable/disable phase inputs/outp Set cycle Vaporizer temperature. Set cycle Preheater temperature. Scroll up screen Cycle Number Selected.	e setpoints. uts.	
•	Set Cycle Name (F1) The cycle name is printed at the	e header of the cy	cle printout.
	CYCLE NAM	1E ## I	

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING	Effective Date	Procedure Number
	PROCEDURE	JAN 5 2005	11116
Biopharmaceutical Development Program		Page 26 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
The u	iser may set two lines of chara	icters.	
•	Print cycle (F2) This printout documents cycl in each cycle printout (refer to paragraph in Operator Mode Set Leak Test Setpoints (F3)	e setpoints and is als o Cycle Report Printo Section).	o provided out
	The cycle name is printed at LEAK TEST Leak Te F1-Enable F5-Disable F4-Previous	the header of the cyc F SETPOINTS est Enable DISABLED F8-Next	cle printout.
	This setpoint (above) enables portion of the cycle. The Lea the enclosure.	s or disables the Lea k Test is a pressure l	k Test eak test of
	LEAK TEST Enclosur 75	T SETPOINTS re Pressure .0 Pa	

0.00-622.00 Pa or "WC.



This setpoint (above) is the cycle airflow maintained during the pressurization portion of the Leak Test, and may be set from 0-999 cfm. The normal range is 8-20 scfm (14-34 cmh).

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 27 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
Title: Operating the VHP [®] 1000ED-A	B Generator LEAK TEST S Max. Press 5 Mini F4-Previous This setpoint (above) is the ma pressurization of the enclosure minutes. If the measured time p Leak Test Failure alarm occurs LEAK TEST S Maximum L 15.0 Pa F4-Previous This setpoint (above) is the max rate of the enclosure during the from 0.00-622.00 Pa/minute or " leak rate higher than this value f alarm and aborts the cycle. LEAK TEST S Pressure St 2 Mini F4-Previous This setpoint (above) is hold pe stopped before the timed Leak f may be set from 0-999 minutes. LEAK TEST S Print In 1.0 Min F4-Previous This setpoint (above) is the enclouse LEAK TEST S Print In 1.0 Min F4-Previous	ETPOINTS ure Time utes F8-Next ximum allowed tin and may be set f beriod exceeds thi and the cycle is a ETPOINTS eak Rate a/min F8-Next imum allowed pre Leak Test, and m WC/minute. A me forces a Leak Test ETPOINTS able Time utes F8-Next riod after pressuri F8-Next riod after pressuri Test phase is star ETPOINTS terval nutes F8-Next osure pressure st cle, and may be s	atus print
	 999.9 minutes. Set Dehumidify Phase Se The cycle name is printed printout. 	tpoints (F5) at the header of t	he cycle

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 28 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
	DEHUMIDIFY S TIM 00:0 HH:M F4-Previous This setpoint (above) is the Dehu	SETPOINTS E 0 IM F8-Next umidify phase tim	le, in hours
	and minutes, and may be set fro minutes:seconds.	m 00:00-99:59	
	DEHUMIDIFY Cycle Airflo 30 cm F4-Previous	SETPOINTS ow Rate nh F8-Next	
	This setpoint (above) is the contr during the Dehumidify phase, an or cmh. The normal operable rar 20 cfm (14 to 34 cmh).	rolled cycle airflov d may be set fror nge of cycle airflo	w rate m 0-999 cfm w is 8.0 to
	DEHUMIDIFY S Absolute H x.x (m F1 – 2.3 F2 – 4 F4-Previous	SETPOINTS lumidity g/l) .6 F3 – 6.9 F8-Next	
T H s a e	This setpoint (above) is the require lumidity level, and may be set to etpoint measures the amount of v specific temperature. Refer to th and of this manual for the Relative	ed Enclosure Abs 2.3, 4.6, 6.9 mg/L water in the enclo e Appendix locate Humidity equiva	olute This osure's air at ed at the lents.
	DEHUMIDIFY S Print Int 10.0 min F4-Previous	SETPOINTS erval nutes F8-Next	
Т С Г	This setpoint (above) is the status Dehumidify phase, and may be se Pertinent data is printed at each p	print interval duri t from 0.0-999.9 rintout.	ing the minutes.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 29 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
Title: Operating the VHP® 1000ED-A	AB Generator • Set Conditioning Phase Set The cycle name is printed printout. CONDITION S H2O2 Inject F1-Enable F5-Disable F4-Previous This setpoint (above) determade every minute of the fame of the family of the fame of the fame of the family of the family of the family of the fa	etpoints (F1) at the header of t ETPOINTS ion Print ENABLED F8-Next rmines whether a measured injection ETPOINTS E 0 IM F8-Next e Condition phase ay be set from 00 ETPOINTS ow Rate nh F8-Next e controlled cycle e, and may be set perable range of 0 34 cmh). ETPOINTS ion Rate min F8-Next	the cycle
	This setpoint (above) is the rate, and may be set from normal operable range of i grams/minute.	e controlled H ₂ O ₂ 0.0-99.9 grams/n njection rate is 1.	injection ninute. The 0 to 12.0

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 30 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
	CONDITION S Print Int 10.0 min F4-Previous This setpoint (above) is the the Condition phase, and r minutes. Pertinent data is Set Decontamination Phase name is printed at the heat DECONTAMINATION H2O2 Inject F1-Enable F5-Disable F4-Previous This setpoint (above) deter made every minute of the DECONTAMINATION TIM xx:x HH:M F4-Previous This setpoint (above) is the time, in hours and minutes 99:59 minutes:seconds. DECONTAMINATION Cycle Airfloo 30 cr F4-Previous This setpoint (above) is the during the decontamination from 0-999 cfm or cmh. Th cycle airflow is 8 to 20 cfm	ETPOINTS erval nutes F8-Next e status print internay be set from 0 printed at each pro- se Setpoints (F2) der of the cycle p ON SETPOINTS ion Print ENABLED F8-Next rmines whether a measured injection ON SETPOINTS E x IM F8-Next e decontamination , and may be set ON SETPOINTS by Rate nh F8-Next e controlled cycle n phase, and may be normal operabl (14 to 34 cmh).	rval during 0.0-999.9 rintout. The cycle rintout is printout is on rate. n phase from 00:00-

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 31 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator		
	DECONTAMINATIO H ₂ O ₂ Inject 12.0 g/ F4-Previous This setpoint (above) is the rate, and may be set from	ON SETPOINTS ion Rate /min F8-Next e controlled H2O2 0.0-99.9 grams/n	injection ninute. The
	normal operable range of i grams/minute.	njection rate is 1.	0 to 12.0
	DECONTAMINATIO Print Int 10.0 min F4-Previous	ON SETPOINTS erval nutes F8-Next	
	This setpoint (above) is the the decontamination phase 999.9 minutes. Pertinent d printout.	e status print inte e, and may be se ata is printed at e	rval during t from 0.0- each
	 Set Aeration Phase Setpoi The cycle name is printed printout. 	nts (F3) at the header of t	he cycle
	AERATION S TIM 00:1 HH:M F4-Previous	ETPOINTS E 5 IM F8-Next	
	This setpoint (above) is the hours and minutes, and ma minutes:seconds. It is reco period be set to no less tha cooling of the VHP Genera	e Aeration phase ay be set from 00 ommended that th an 15 minutes, to ator heaters.	time, in 0:00-99:59 nis time allow for

National Cancer Institute-Frederick, Frederick, MD	STAI OPE PRO	NDARD RATING CEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program			Page 32 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	B Generator			
	F	AERATION S Cycle Airflo 30 cr 4-Previous	ETPOINTS ow Rate nh F8-Next	
	This set during t cfm or c airflow i	point (above) is the he Aeration phase, mh. The normal o s 8 to 20 cfm (14 to	e controlled cycle and may be set perable range of o 34 cmh).	airflow rate from 0-999 cycle
		AERATION S Print Int 30.0 mir F4-Previous	ETPOINTS erval nutes F8-Next	
	This set the Aera minutes	point (above) is the ation phase, and m . Pertinent data is _l	e status print inter ay be set from 0. printed at each pr	rval during 0-999.9 intout.
		AERATION S Aux. Aerati 00:15 HH:M F4-Previous	ETPOINTS on TIME 5 1M F8-Next	
	This set time, in 99:59 m	point (above) is the hours and minutes ninutes:seconds.	e Auxiliary Aeratio , and may be set	on phase from 00:00-
		AERATION SI Blower O Aeration F1-Yes F4-Previous	ETPOINTS n After : YES F5-No F8-Next	
	Press F Aerate the cha (Auxilia	1 to select the char phase (Auxiliary Ae mber blower OFF a ry Aerate).	mber blower ON a erate), or press F& after the Aerate pl	after the 5 to have nase

PROCEDURE	JAN 5 2005	11116
nerator	Page 33 of 43	
nerator	1 age 55 01 +5	Revision 00
Set Pressure Control Setper The cycle name is printed a printout. ENCL. PRESS. CYCLE PRESS PRESS F1-PRES F2-VA F4-Previous Press F1 (to control at post control at negative pressure pressure control or press F control. ENCL. PRESS. Nominal Pre- 200.00 F4-Previous This setpoint (above) is the Pressure Differential, and r 999.99 Pa or "wc. Use the for the Nominal Pressure D range is 0.00 to 622.00 Pa NOTE: It is the Operator's selected Operating Pressure exceed the enclosure man recommendations. ENCL. PRESS. Low Pressure Setpoint: 37 Abort Tm: 3 F4-Previous This setpoint (above) is for the low limit that the enclos or vacuum before an alarm setpoint is closer to atmost	oints (F5) at the header of t SETPOINTS CONTROL URE AC F3-NO F8-Next itive pressure) or e) to select the cy 3 for NO cycle pl SETPOINTS ess. Diff: Pa F8-Next SETPOINTS ess. Diff: Pa F8-Next controlled cycle may be set from (touch pad to ente Differential. The n (0.00 to 2.50 "wo responsibility to e re (or Vacuum) d ufacturer's SETPOINTS re Alarm 70.00 Pa 5 sec. F8-Next the Low Pressur sure will maintain n occurs. The low	he cycle F2 (to ycle ressure Nominal 00.00- er a value ormal c). ensure the oes not e Alarm, a pressure alarm bh alarm
	Set Pressure Control Setpe The cycle name is printed a printout. ENCL. PRESS. CYCLE PRESS PRESS F1-PRES F2-VA F4-Previous Press F1 (to control at post control at negative pressure pressure control or press F control. ENCL. PRESS. Nominal Pre 200.00 F4-Previous This setpoint (above) is the Pressure Differential, and r 999.99 Pa or "wc. Use the for the Nominal Pressure E range is 0.00 to 622.00 Pa NOTE: It is the Operator's selected Operating Pressure exceed the enclosure man recommendations. ENCL. PRESS. Low Pressure setpoint: 37 Abort Tm: 3 F4-Previous This setpoint (above) is for the low limit that the enclose or vacuum before an alarm setpoint is closer to atmosp setpoint (see above). Once	Set Pressure Control Setpoints (F5) The cycle name is printed at the header of treprintout. ENCL. PRESS. SETPOINTS CYCLE PRESS. CONTROL PRESSURE F1-PRES F2-VAC F4-Previous F8-Next Press F1 (to control at positive pressure) or control at negative pressure) to select the cypressure control or press F3 for NO cycle pressure control or press F3 for NO cycle pressure control. ENCL. PRESS. SETPOINTS Nominal Press. Diff: 200.00 Pa F4-Previous F5-Next This setpoint (above) is the controlled cycle Pressure Differential, and may be set from 0 999.99 Pa or "wc. Use the touch pad to enter for the Nominal Pressure Differential. The n range is 0.00 to 622.00 Pa (0.00 to 2.50 "wc NOTE: It is the Operator's responsibility to e selected Operating Pressure (or Vacuum) d exceed the enclosure manufacturer's recommendations. ENCL. P

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program		Page 34 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-AB (Generator		
	aborted AFTER the LC has expired.	W Pressure Alarm Ab	ort timeout
	ENCL. PRE High Pre Setpoi Abort F4-Previous	SS. SETPOINTS essure Alarm int: 0.00 Pa Tm: 35 sec. F8-Next	

alarm setpoint (see below). Once the proper enclosure pressure has been reached, if the enclosure pressure exceeds the high alarm setpoint, the cycle will be aborted AFTER the HIGH Pressure Alarm Abort timeout has expired.

> ENCL. PRESS. SETPOINTS Pressurize Timeout xx:xx MM:SS F4-Previous F8-Next

Use the touch pad to enter a value for the Pressurize Timeout. This timeout (maximum time to reach pressure) is the maximum time for the VHP Generator to reach the enclosure pressure or vacuum. If this time is exceeded, the cycle ABORTS.

• Set Enclosure Control (F1) The cycle name is printed at the header of the cycle printout.

National Cancer Institute-Frederick, Frederick, MD	ST/ OP PR(ANDARD ERATING OCEDURE	Effective Date JAN 5 2005	Procedure Number 11116
Biopharmaceutical Development Program			Page 35 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-AB Generator				
	F1 = exam F2 = exam F3 = flexib F4 = This s used 0.0-9! • I/O C The c printc F1 = for thi F2 = for thi F3 = flexib F4 = F3 = flexib	ENCL. SET Enc. Type: F F1-Rigid F3-Semi Rigid F4-Previous If the enclosure is ma ple, glass). If the enclosure is ma le and rigid materials. Scroll down screen. ENCL. SET Enc. Volume F4-Previous setpoint (above) is the for the selected cycle 999 cm (0.0-9999 cf). ontrol Setpoints (F2) cycle name is printed a but. I/O CON Utilize I/O C YES F1-YES F4-Previous To enable the use of is cycle. To disable the use of is cycle. To disable the use of is cycle. If the enclosure is ma le and rigid materials. Scroll up screen. Scroll down screen	POINTS LEXIBLE F2-Flexible F8-Next de of a rigid mate de of a flexible m de of a flexible m de of a combinat POINTS e: xxx cm F8-Next e enclosure volum the header of t TROL Control? F2-NO F8-Next Phase Inputs and de of a combinat	erial (for aterial (for ion of he to be ne range is he cycle d Outputs d Outputs ion of

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116			
Biopharmaceutical Development Program		Page 36 of 43	Revision 00			
Title: Operating the VHP [®] 1000ED-A	Title: Operating the VHP [®] 1000ED-AB Generator					
	I/O CON Dehumidify Phase Output = No F1-YES F4-Previous Next	TROL O F5-NO F8-				
F1 = To enable the use of Phase Output for the Dehumidify Phase. F5 = To disable the use of Phase Output for the Dehumidify Phase. F4 = Scroll up screen. F8 = Scroll down screen						
	NOTE : If Phase Output is selected for a phase (YES), the output will only be ON if I/O Control is enabled (YES).					
	Use F4/F8 Keys to scroll between screens and enable or disable Phase Outputs for the following phases.					
 Dehumidify Condition Decontaminate Aeration Auxiliary Aeration Out of Cycle 						
	<u>NOTE</u> : The external outpu VDC, maximum 1.0 Amp.	t may be 24 VAC	or 24			
	Continue to scroll down to Inputs.	find settings for F	Phase			

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116			
Biopharmaceutical Development Program		Page 37 of 43	Revision 00			
Title: Operating the VHP [®] 1000ED-AB Generator						
			_			
	I/O CON Dehumidify Phase Input = NO F1-YES F4-Previous	F5-NO F8-Next				
	F1 = To have the external input detected for the Dehumidify Phase. F5 = To ignore the external input for the Dehumidify Phase. F4 = Scroll up screen. F8 = Scroll down screen.					
	NOTE : If Phase Input is selected for a phase (YES), the input will only be detected if I/O Control is enabled (YES) for the cycle. Also, if the Phase Input is enabled for a phase, and the external input is open during the selected phase, the cycle automatically ABORTS.					
	Use F4/F8 Keys to scroll between screens and enable or disable Phase Inputs for the following phases.					
 Dehumidify Condition Decontaminate Aeration 						
	 Vaporizer Setpoints (F3) The cycle name is printed at the header of the cycle printout. 					
	VAPORIZER SETPOINTS Vap. Temp. (In Cycle) 100.0 C					
	F4-Previous	F8-Next				
This setpoint (above) is the temperature at which the vaporizer heater is maintained when the VHP Generator is in-cycle, and may be set from 0.0- 300.0° C/ $^{\circ}$ F. This value is normally set to 100.0° C to ensure proper vaporization of the H ₂ O ₂ during injection phases.						

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116			
Biopharmaceutical Development Program		Page 38 of 43	Revision 00			
Title: Operating the VHP [®] 1000ED-A	Title: Operating the VHP® 1000ED-AB Generator					
 Preheater Setpoints (F5) The cycle name is printed at the header of the cycle printout. PREHEATER SETPOINTS Pre. Temp. (In Cycle) 90.0 C F4-Previous F8-Next 						
This setpoint (above) is the temperature at which the preheater is maintained when the VHP Generator is in- cycle, and may be set slightly below the "Vap. Temp. (In Cycle)" value to aid the vaporizer in maintaining its setpoint temperature during injection phases.						
6.5.3.10 Save Cycle Setup Changes Press F4 at one of the Cycle Setup screens and after made changes the following screen appears.						
	Save the Changes to Cycle nn? F2-E F1-Save F5-E	dit xit				
Where: nn=the selected cycle number (1-13) F2 = Return to Cycle Setup for the selected cycle. F1 = Save changes made to selected cycle and exit Cycle Setup. F5 = Ignore changes made to selected cycle and exit Cycle Setup.						
7.0 Sequence of Operation for the VHP Generator						
7.1 Installing Sterilant Cartridge The control displays a message when the sterilant cartridge needs replacement. Also, the windowed door on the sterilant cartridge compartment allows to easily check the cartridge to see when it is empty.						
7.1.1 Put prot	on chemical splash goggles, viny ective clothing.	/l or neoprene glo	oves and			

National Cancer Institute-Frederic Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number		
Rightermoontical Development Brown		JAN 5 2005	11116		
		Page 39 of 43	Revision 00		
Title: Operating the VHP [®] 1000E	Title: Operating the VHP [®] 1000ED-AB Generator				
7.1.2 T p	urn the cartridge control knob clock osition and open the cartridge comp	wise 180° to the partment door.	REPLACE		
7.1.3 C ti h	ently grasp the spent cartridge with ne cartridge about 1/16" and then re older.	both hands. Car move the cartride	efully lift ge from the		
7.1.4 V c tl	/ith water running in a sink, empty a ontents into the sink. When the cart noroughly rinse the cartridge with ta	all remaining cartı ridge is empty, ca p water before di	ridge arefully and sposal.		
7.1.5 F d e r	emove a new cartridge from the ca ate on the new cartridge. Do not us xpiration date listed on the label. Ac eleased VHP cartridge.	move a new cartridge from the carton. Check the expiration te on the new cartridge. Do not use the cartridge if beyond its piration date listed on the label. Additionally, only use QC eased VHP cartridge.			
7.1.6 F	emove the vented shipping cap from the new cartridge.				
7.1.7 C g ir h	ently grasp the new cartridge with l roove in the cartridge holder, and ca to the holder. Lower the cartridge in older.	ooth hands, align arefully slide it all nto the 1/16" rece	it with the the way ess in the		
7.1.8 T p	urn the cartridge control knob clock osition. Do not turn the knob past th	wise 180° to the lis position.	ENGAGE		
7.1.9 T ti C ti ti	he cartridge is now locked in place ne cartridge control knob is turned to once the cartridge control knob is tu ne previously used cartridge is rend ne amount of liquid remaining.	and cannot be re o the REPLACE p rned to the replac ered unusable rep	moved until position. ce position, gardless of		
7.1.10 C c a	lose the cartridge compartment do ommand from the PLC and the VH utomatically fill to the programmed	or. Select "Reser P Generator unit level.	voir Fill" will now		
NC cai Do ins Ru	TE: Always turn the cartridge contr tridge is not properly installed if it do not use the cartridge beyond its ex tallation. nning a Decontamination Cycle	ol knob clockwise oes not rest in the piration date or 4	e. The e recess. 5 days after		

Nation	ational Cancer Institute-Frederick, Frederick, MD		STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116		
Biopha	rmaceutical Development I	Program		Page 40 of 43	Revision 00		
Title:							
	7.2 Running a Decontamination Cycle						
	7.0	й 1 Тра	, avala proceede op abown in the	following table. C)porotor		
	1.2.	acti	on may be required during the Re	eservoir Fill phase	e and at		
		Cyc	le Complete. A cycle printout is g	enerated if the sy	/stem		
		prin	ter in enabled.				
			Decontamination Cycle Sequence	<u>a</u>			
	Phase Name		Process Description	Phase A	dvance		
			ľ	Require	ement		
	Standby or Main Menu	Idle stat	te, ready to run cycle	Start Cycle cor	nmand		
	Reservoir Fill	If the R	eservoir does not contain enough	Required reser	voir weight		
		H ₂ O ₂ to	complete the selected cycle, the	(per the select	ed cycle) is		
	luis et Drive e	reservo	ir is filled from the $H_2O_2Cartridge$.	met			
	Inject. Prime	The inje	ection line is primed with H_2O_2 liquid.	Setpoint time e	xpired		
		line in p	reparation for injection and also				
		purges	the fill line.				
	Dehumidity	Airflow	through VHP Generator is maintaine	d Humidity level	meets		
		at Dehu	imidify phase cycle airflow setpoint.	setpoint and se	etpoint time		
		Phase of humidity	operates by time and by relative	expired			
	PreHTR Warm-up	The pre	heater is heated to setpoint	Temperature n	neets setpoint		
	· · · · · · · · · · · · · · · · · · ·	tempera	ature.				
	Vapor Warm-up	The vap	porizer is heated to setpoint.	Temperature n	neets setpoint		
	Condition	Airflow	control is prointained per Condition	temperature.	waired		
	Condition	AIMOW phase s	control is maintained per Condition	Selpoint time e	xpired		
		priase s	e begins and is controlled according				
		to set ra	ate for the Condition phase.				
	Decontaminate	Airflow	control is maintained per	Setpoint time e	xpired		
		Deconta	amination phase setpoint. Injection c	of			
		vapor h	ydrogen peroxide begins and is				
			eu according to set rate for the				
	Aeration	Airflow	control is maintained perAerate	Setpoint time e	expired		
		phase s	setpoint and injection is stopped.		·		
	Auxiliary Aeration	Airflow	control is maintained per Aerate	Setpoint time e	xpired or		
		phase s	etpoint if the blower is enabled	manual phase	advance.		
	Cuelo Complete	during t	nis phase.	Operator calin	owlodamont		
		AILIOW	is stopped. Cycle is completed.		Jwieugineni		

National Cancer Institute-Fredo Frederick, MD	erick,	STANDARD OPERATING	Effective Date	Procedure Number
		PROCEDURE	JAN 5 2005	11116
Biopharmaceutical Development Pro	gram		Page 41 of 43	Revision 00
Title: Operating the VHP [®] 100	00ED-AB (Generator		
7.2.2	During Genera of H ₂ O than th	the cycle Reservoir Fill pha ator determines the VHP G 2 required for the cycle and e required amount, the res	ase of the cycle, the venerator determines , if the Reservoir con ervoir filled with H ₂ O ₂	/HP the amount tains less
7.2.3	The re amoun injectio	quired amount is 20 grams t per the Condition and Ste on rates.	more than the neces rilize phase setpoint	sary times and
7.2.4	Filling	occurs by means of the foll	owing procedure.	
	7.2.4.1	Before filling begins, the check of the filling system	VHP Generator perfo า.	rms a brief
	7.2.4.2	If the filling system check not successful, the cycle message is shown.	is successful, filling l is aborted and an ala	oegins. If Irm
	7.2.4.3	If the VHP Generator dete and the cycle is not yet co with the message: "H ₂ O ₂ stopped. Proceed as follo	ects the H ₂ O ₂ cartrido omplete, an alarm is i SUPPLY EMPTY" ar ows.	ge is empty indicated nd filling is
		 Set the Cartridge Control remove the cartridge. 	ol Knob to REPLACE	and
		 If the cartridge is not en filling system or with the install a new cartridge. 	npty, there is a proble e cartridge itself; othe	em in the rwise,
		 Set the knob to ENGAG and filling will resume. 	E, the alarm screen	will clear
	7.2.4.4	Filling is completed once the required amount.	the Reservoir weight	has met
7.2.5	Depen Genera signal	ding upon the External Inte ator will pause at certain cy is received from the interfac	rface options used, tl cle phases until a hai ce.	ne VHP ndshake
7.2.6	lf a cyc pause require	cle is started from the Host at the Vaporizer Warm-up ement is met, until a signal r	PLC, the VHP Gener phase after its norma received from the inte	rator will al phase erface.

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date	Procedure Number
		JAN 5 2005	11116
Biopharmaceutical Development Program		Page 42 of 43	Revision 00
Title: Operating the VHP [®] 1000ED-A	AB Generator		
7.2.7 This Ger equ this	s is done to allow external equipm nerator for the cycle at the same t lipment prepares the enclosure fo time, the display will show the fol	nent to prepare the ime that the exter r decontamination llowing:	e VHP nal n. During
Cy Tir Ph -A F4	rcle-#nn RD me Rem.: 00:00:00 hase: Vapor. Warm-up CTIVE- F5-ABORT -Previous F8-Next		
Where "RD" in the Ext	: nn = selected cycle, set per the idicates that the cycle is paused u ternal Interface.	Cycle Selection so Intil a signal is rec	creen. ceived from
7.2.8 At c ack PLC	completion of the Decontamination nowledgement is received from th C. The display shows the following	n cycle, the Unit h ne operator or fror g:	olds until m the Host
Cy Tir Ph F2 F4	vcle-#nn ne Rem.: 00:00:00 nase: Cycle Complete 2-RESET -Previous F8-Next		
Where: r	nn = selected cycle, set per the C	ycle Selection scr	een.
7.2.9 Ope One PLC	erator acknowledgement occurs b ce acknowledgement is received (C), the VHP Generator returns to t	y pressing the F2 (from the operato the Main Menu.	touchpad. r or Host
NO "AB	<u>TE</u> : If the cycle is an aborted cycl ORT" to the right of the cycle nur	e, the display indi nber indication.	cates
7.3 Running a Re	generation Cycle After Decontam	ination Cycle.	
7.3.1 The rea	e Run Screen will appear on PV30 d: "Reagent HEAT UP". Press F1	0 Display. Phase to START cycle.	should
	Cycle-#nn Time Rem.: 00:00:00 Phase: Regen HEAT UP F1-START F5-ABOR F4-Previous F8-Next	Т	

National Cancer Institute-Frederick, Frederick, MD	STANDARD OPERATING PROCEDURE	Effective Date JAN 5 2005	Procedure Number 11116	
Biopharmaceutical Development Program		Page 43 of 43	Revision 00	
Title: Operating the VHP [®] 1000ED-AB Generator				
7.3.2 On Sci	ce the cycle has started, use F8 to een.	o get to "Dryer Te	emp."	
	Dryer Temp. S:aaaa. A: bbbb xx H ₂ O ₂ Return Temp. S: cccc. A: dddd xx	(
Whe	re: aaaa=Setpoint Dryer temperat bbbb=Actual Dryer temperature. cccc=Setpoint H ₂ O ₂ Retyrn temper dddd=Actual H ₂ O ₂ Return temper Xx=C or F, depending upon temp	ure. erature. ature. erature unit set.		
7.3.3 The hea 212 Tei	e HEAT UP phase of cycle begins at up to at least 400°F (204°C). Th 2°F (100°C) is achieved by H ₂ O ₂ r mperature Detector).	. The Dryer temp ne phase will end eturn RTD (Resis	erature will when stive	
7.3.4 The coo act wh	e COOL DOWN phase will follow. of down. The phase will end when nieved by H ₂ O ₂ return RTD, or afte ichever occurs first.	The dryer tempe 130°F (54.4°C) c er two hours has c	rature will or less is elapsed,	
7.3.5 Wh Op Pre	en the cycle is complete, press F erator Screen for a complete dupl ess F2 (RESET) to advance to the	1 (PRINT REPOR icate print of the Main Operator S	RT) from the cycle. Screen.	
7.4 Manually Dra	in Reservoir			
7.4.1 The onl	e manually Drain Reservoir Cycle y. Only qualified Service Technicia	is available in Se an should perforn	ervice Mode n this cycle.	
8.0 <u>Reference</u>				
8.1 Operator Mar 120 (09/11/03	nual VHP® 1000ED-AB Biodecont 3).	tamination Syster	m P129383-	