



## Standard Operating Procedure

Biopharmaceutical Development Program

### Title: **Pierce Modified Lowry Protein Assay**

**SOP Number: 22146**

**Revision Number: 01**

**Supersedes: Revision 00**

**Effective Date: DEC 7 2010**

**Originator/Date:**

**Approval/Date:**

**Approval/Date:**

#### Table of Contents

- 1.0 Purpose
- 2.0 Scope
- 3.0 Authority and Responsibility
- 4.0 Materials and Equipment
- 5.0 Procedure
  - 5.1 Preparation of Solutions and Standards
  - 5.2 Sample Preparation
  - 5.3 Assay Protocol
- 6.0 Interpretation of Results
- 7.0 Documentation
- 8.0 References and Related Documents
- 9.0 Attachments

#### 1.0 Purpose

This procedure describes how to determine the amount of total protein in a sample using the Pierce Modified Lowry Protein Assay Kit. The Lowry method for determining protein concentration is essentially a biuret reaction, which incorporates the use of Folin-Ciocalteu reagent for enhanced color development. The procedure described below is the original Lowry protocol (see References) except that Reagent C (an alkaline copper solution that is a mixture of Reagent B and A) is replaced by Pierce's Modified Lowry Protein Assay Reagent.

**Title: Pierce Modified Lowry Protein Assay**

---

**2.0 Scope**

This procedure involves reaction of protein with cupric sulfate and tartrate in an alkaline solution, resulting in formation of tetradentate copper-protein complexes. When the Folin-Ciocalteu Reagent is added, it is effectively reduced in proportion to these chelated copper complexes, producing a water-soluble product whose blue color can be measured at 750 nm. The color response curves for the Modified Lowry Protein Assay and the original Lowry method have nearly 100% correlation. This procedure applies to Process Analytics staff who use this assay.

**3.0 Authority and Responsibility**

- 3.1 The Director, Process Analytics (PA) has the authority to define this procedure.
- 3.2 PA is responsible for training laboratory personnel and documenting this training to Biopharmaceutical Quality Assurance (BOA).
- 3.3 PA personnel are responsible for the performance of this procedure.
- 3.4 PA is responsible for reviewing the data and documentation of the results of this procedure.
- 3.5 BOA is responsible for quality oversight of this procedure.

**4.0 Materials and Equipment.**

- 4.1 Modified Lowry Protein Assay Reagent Kit (BDP PN 30283) containing:
  - Modified Lowry Protein Assay Reagent (Pierce No. 1856006);
  - 2 N Folin-Ciocalteu Reagent (Pierce No. 1856007); and
  - Bovine Serum Albumin (BSA) Standard (BOP PN 30060); BGG No. 30059 or equivalent reagents, which may be ordered separately.
- 4.2 Saline, phosphate buffered saline or other appropriate diluent (i.e., formulation buffer of test samples).
- 4.3 Beckman DU650i Spectrophotometer, or equivalent.
- 4.4 Disposable polystyrene cuvettes: 1.5 ml semi-micro (BDP PN 20070), 4.5 mL standard (Fisher PN 14-385-942), or BOP approved equivalent cuvettes.
- 4.5 Water, HPLC grade or better.
- 4.6 Assorted Micro-pipets, tips, culture tubes and laboratory glassware.

**5.0 Procedure****5.1 Preparation of Solutions and Standards**

- 5.1.1 Record the lot number of the Modified Lowry Protein Reagent, and its expiration date on the Lowry Protein Assay Data Sheet (Attachment 1).

5.1.2 Prepare 1X {1N} Folin-Ciocalteu reagent by diluting the supplied 2X reagent (2N) with an equal volume of water. Prepare fresh reagent each day as it is unstable in the diluted form. Record the lot number of the Folin-Ciocalteu reagent on the Lowry Protein Assay Data Sheet.

5.1.3 Record the concentration and lot number of the BSA Standard on the Lowry Protein Assay Data Sheet.

5.1.4 Prepare a set of protein standards of known concentration by diluting the BSA standard solution provided with this kit, or any other suitable standard, in the same diluent as the unknown samples. The set of standards should consist of eight dilutions of BSA over the range of 1-1500 µg/ml and should cover the range of protein concentrations being used. Record the dilution details on the Lowry Protein Assay Data Sheet. Triplicate standard samples require at least 600 µL of each concentration.

5.1.5 Also, analyze a sample of diluent buffer and use it as a background blank to zero the spectrophotometer. Record the description, including the lot number of the diluent buffer, on the Lowry Protein Assay Data Sheet (Attachment 1).

## 5.2 Sample Preparation

5.2.1 The samples and the protein standards must be diluted in the same diluent buffer.

5.2.2 Dilute the samples so that the approximate concentration is in the range of 20 - 1000 pg/ml or make a series of two-fold dilutions so that two or more of the dilutions will fall in this range. For samples of unknown protein concentration, the assay may require a second set of dilutions when all the readings in the first set are out-of-range from the standard curve.

5.2.3 Assay triplicates require at least 600 µL of each concentration.

## 5.3 Assay Protocol

5.3.1 Assay triplicate preparations of each concentration of standard and unknown sample.

5.3.2 To 200 µL of each standard or unknown samples, add 1.0 mL of the Modified Lowry Protein Assay Reagent.

5.3.3 Vortex the samples.

5.3.4 Incubate the samples at room temperature (18-25°C) for 10 minutes.

5.3.5 Add 10\_0 µL of the 1 N Folin-Ciocalteu Reagent while vortexing.

5.3.6 Incubate the samples at room temperature (18-25°C) for 30 minutes.

5.3.7 Measure the absorbance of the standards and samples at 750 nm. Refer to **SOP 22158, Operation of the Beckman DU Series 600 Spectrophotometer** for operation of the Beckman DU Series 600 spectrophotometer and to the Beckman DU Series 600 Spectrophotometer Protein Analysis Operating Instructions (Manual 517303B) for specific instructions for the protein assay.

FOR TITRATION & REFLUENCE ANALYSIS

5.3.8 Prepare a standard curve by plotting the net (blank corrected) absorbance at 750 nm versus the protein concentration. Using this curve, determine the protein concentration of the unknown samples. The Beckman DU650i Spectrophotometer will perform this operation automatically when the PROTEIN operation is performed.

5.3.9 When using the PROTEIN operation on the Beckman DU650i Spectrophotometer, obtain printouts of the following data:

Protein Analysis Standards Window with  $A_{750}$  values;  
Standard Curve of Concentration versus  $A_{700}$  values; and  
Protein Analysis Samples Window with  $A_{700}$  values.

5.3.10 Alternatively, the raw  $A_{750\text{nm}}$  data may be analyzed using approved spreadsheet software such as Origin, SAS or Microsoft Excel.

5.3.11 When using second-party software for data analysis, printouts containing the same graphic data plots as in 5.3.9 are required.

### Interpretation of Results

6.1 The protein concentration of the sample is obtained by comparison of the  $A_{15'}$  value of the sample with the standard curve. Solving for unknown's concentration "x."

$$fx = (A_{100} - \text{STD } Y \text{ INI} \therefore \text{RCEPT} / \text{STD.SLOPE}]$$

6.2 2 Since the Folin-Ciocalteu reagent reacts with tyrosine and tryptophan residues, the Lowry assay is dependent on the unknown and the standard protein having a similar content of tyrosine + tryptophan/ $\mu\text{g}$  protein. This method is **NOT** appropriate for proteins without tyrosine and tryptophan residues, or proteins known to contain significant proportions of tyrosine and tryptophan residues.

### 6.: Acceptance/Rejection Criteria

Standard Curve -.. If the PROTEIN operation on the Beckman DU650 is used, after all the designated standards are analyzed, the standard curve can be displayed. Curve fit coefficients are displayed above the standard curve. Clicking on <ViewStats> from the Standards window will display the Standard Stats window which can be used to determine which, if any, of the standards could be rerun or deleted to improve the accuracy of the standard curve.

No more than 4 individual standard readings may be deleted (out of 24 total).

Samples - If the value(s) obtained for an unknown sample is outside the range of the standard curve, repeat the assay with an appropriate set of dilutions to obtain values within the correct range.

6.4 The resulting standard curve  $R^2$  value must be  $\geq 0.90$  for the assay to be valid.

## 7.0 Documentation

7.1 Following the run, the data should be available:

Protein Analysis Standards  $A_{750}$  Values;  
Standard Curve; and  
Protein Analysis Samples  $A_{100}$  Values.

7.2 Prepare a final report that shall include for each test article the printout of the  $A_{750}$  values of the standards, the standard curve and associated graphs, the  $A_{100}$  values of the test samples and the calculated protein concentrations of the test samples.

7.3 BQA Documentation will maintain the raw data, final report, documentation and protocols.

## 8.0 References and Related Documents

8.1 **SOP 22158** *Operation of the Beckman DU Series 600 Spectrophotometer*

8.2 Lowry, O.H., et al, J. Biol. Chem. 193, 267 - 275(1951).ist all other non BOP references.

8.3 Pierce Modified Lowry Protein Assay Reagent Kit Instructions (Manual 0389).

8.4 Beckman DU Series 600 Spectrophotometer Operation Instructions (Manual 5173008).

8.5 Beckman DU Series 600 Spectrophotometer Protein Analysis Operating Instructions (Manual 5173038).

## 9.0 Attachments

9.1 **Attachment 1** Form 22146-01, Process Analytics Laboratory Lowry Protein Assay datasheet

**u**  
**S**  
**T**  
**R**  
**o**

**&**  
**C**  
**P**

Process AflaMics  
Laboratory Lowry Protein Assay Oatasheet  
For Sample\$ Analyzod with the Pierce Modlified Lowry Protein Msay

Descrição .....

M.Y.:Wia:1 Lowry Protein Reagent BOP Lot No. .... Exp. ....  
2 N Folin-Coultter Reagent ODP Lot No. .... Exl .....  
Standard Compound: 80 P. 1 Ne ..... Exl .....

0"11Cr iio:1.

<p>Name: <u>me_e</u> -1</p> <p>Strom: <u>0</u> / <u>1</u> -1</p> <p>Obwohl: <u>n</u> / <u>1</u> -1</p> <p>Conc: <u>1</u> / <u>1</u> -1</p>	<p>Strom: <u>0</u> / <u>1</u> -1</p> <p>Obwohl: <u>n</u> / <u>1</u> -1</p> <p>Conc: <u>1</u> / <u>1</u> -1</p>	<p>Strom: <u>0</u> / <u>1</u> -1</p> <p>Obwohl: <u>n</u> / <u>1</u> -1</p> <p>Conc: <u>1</u> / <u>1</u> -1</p>	<p>Strom: <u>0</u> / <u>1</u> -1</p> <p>Obwohl: <u>n</u> / <u>1</u> -1</p> <p>Conc: <u>1</u> / <u>1</u> -1</p>
--	--	--	--

He 'l e-ma By/Oate: