Frederick National Laboratory for Cancer Research, Frederick, MD

**BiopharmaceutIcal Development Program** 

#### Standard Operating Procedure

# Title: Quantitative PCR (qPCR) Method for Detection of Replication Competent Lentivirus

SOP Number: 22209

Revision Number: 00

Supersedes: Revision N/A

Effective Date: MAY 24 2019



Approval/Date:

Approval/Date:

## Table of Contents

- 1.0 Purpose
- 2.0 Scope
- 3.0 Background
- 4.0 Authority and Responsibility
- 5.0 Equipment, Materials and Reagents
- 6.0 Procedure
- 7.0 Validity of the Assay
- 8.0 Analysis of Results
- 9.0 Documentation
- 10.0 References and Related Documents
- 11.0 Attachments

## 1.0 Purpose

This SOP describes the materials and methods needed to detect replication-competent lentivirus (RCL) in a sample by quantitative PCR (qPCR).

#### 2.0 Scope

This method relies upon the ability of the primer/probe sequences to bind to VSV-G sequence that would be present if a replication-competentlentivirus were present in a sample. **If** the VSV-G plasmid used in the generation of the lentivirus vector possesses a VSV-G sequence that differs from the primer/probe sequences used in this method, this SOP would not be appropriate for detecting replication competent lentivirus in that sample.

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

### 3.0 Background

Lentiviral vectors are engineered to be replication defective, however, replication competent lentiviruses could be generated through recombination events among the plasmids used in the lentivirus vector production. One of these plasmids used in lentivirus vector production is the envelope plasmid encoding VSV-G which is used for virus pseudotyping. The replication competent lentivirus qPCR assay targets the VSV-G sequence to determine whether replication competent lentivirus was generated.

#### 4.0 Authority and Responsibility

- 4.1 The Director, Process Analytics/Quality Control (PA/QC), has the authority to define this procedure.
- 4.2 PA/QC is responsible for training laboratory personnel and documenting this training to Biopharmaceutical Quality Assurance (BQA)
- 4.3 PA/QC personnel are responsible for the performance of this procedure.
- 4.4 PA/QC is responsible for reviewing the data and documentation of the results of this procedure.
- 4.5 BQA is responsible for quality oversight of this operation.

## 5.0 Equipment, Materials and Reagents

- 5.1 Clinical Center VSVG Forward Primer (BDP PN 31162)Sequence: 5' TGCAAGGAAAGCATTGAACAA 3'
- 5.2 Clinical Center VSVG Reverse Primer (BDP PN 31163)
  Sequence: 5' GAGGAGTCACCTGGACAATCACT 3'
- 5.3 Clinical Center VSVG Probe (BDP PN 31164)
  Sequence: 5' FAM-AGGAACTTGGCTGAATCCAGGCTTCC-QSY 3'
- 5.4 TaqMan Universal PCR Master Mix (2X) (BDP PN 30268)
- 5.5 MicroAmp Optical 96 Well Reaction Plate (BDP PN 21141 or BDP-approved equivalent)
- 5.6 MicroAmp Optical Adhesive Film (BDP PN 21142 or BDP-approved equivalent)
- 5.7 Distilled Water, DNase Free RNase Free (BDP PN 10189 or BDP-approved equivalent)
- 5.8 1X TE (BDP PN 30267 or BDP-approved equivalent)
- 5.9 1X PBS (BDP PN 30007 or BDP-approved equivalent)
- 5.10 VSV-G DNA Control (BDP PN 31238)
- 5.11 qPCR Instrument QuantStudio 6, SDS 7900HT or BDP-approved equivalent

#### 6.0 Procedure

- 6.1 Extract genomic DNA from each sample as described in SOP 22212 Purification of DNA Using the DNeasy Blood and Tissue Kit or SOP 22972 Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions.
  - 6.1.1 Perform two or more replicate extractions for each transduced test sample.
  - 6.1.2 If available, a non-transduced sample should be extracted in parallel with the transduced samples.

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

- 6.1.3 An Extraction Negative Control should be performed in which 200 μL of 1X PBS are extracted in parallel with the test samples.
- 6.1.4 An Extraction Positive Control should be performed in which 20 μL of VSV-G DNA Control at 1.00E+03 copies/μL is added to 180 μL of buffer blank or 1X PBS.
  - **NOTE:** The Extraction Positive Control should not be centrifuged prior to the addition of Proteinase K and Buffer AL.
- 6.2 Quantify the amount of DNA in each test sample using the NanoDrop ND1000 spectrophotometer according to SOP 22946 The Operation of the NanoDrop ND1000 Spectrophotometer.

**NOTE:** It is not necessary to perform a spectrophotometry analysis of the Extraction Negative Control and Extraction Positive Control.

- 6.3 Dilute an aliquot of each sample to 8 ng/ $\mu$ L with 1X TE in a total volume of 100  $\mu$ L (total DNA is 800 ng).
- 6.4 Replication Competent Lentivirus qPCR Master Mix Preparation
  - 6.4.1 Record all reagent part numbers, lot numbers and expiration dates as well as any additional comments on **Form 22209-01**.
  - 6.4.2 Vortex the reagents.

Reagent	Volume for 1 Reaction	Final Concentration
TaqMan Universal PCR Master Mix (2X)	25.0	1X
Clinical Center VSVG Forward Primer (10 <b>µM)</b>	1.0	200 nM
Clinical Center VSVG Reverse Primer (10 µM)	1.0	200 nM
Clinical Center VSVG Probe (FAM/QSY) (10 µM)	0.5	100 nM
Water	10	N/A

6.4.3 Prepare a qPCR master mix as described below:

- 6.4.4 Prepare a sufficient volume of master mix to account for all test samples, standard curve points and controls with 10% excess to account for losses due to repeated pipetting.
- 6.4.5 Record the volumes used on **Form 22209-02**.

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

- 6.4.6 Mix the master mix by pipetting up and down.
- 6.4.7 Add 37.5 μL of master mix to each well of the 96-well PCR plate that will be used.
- 6.4.8 Add a 12.5 μL aliquot of 1X TE to each of three No Template Control (NTC) wells.
- 6.4.9 Seal the plate with adhesive film. The plate may be stored at 2-8°C for up to 1 hour while samples are prepared for addition to the plate.
- 6.5 Sample addition to the PCR plate
  - 6.5.1 Test Sample: Add 12.5 µL of each extracted sample (5 replicate wells), Extraction Negative Control and Extraction Positive Control to the appropriate wells of the PCR plate.
  - 6.5.2 Inhibition-control Reaction: Add 12.5 μL of non-transducedcontrol sample+ 2 μL of VSV-G DNA Control at 1.00E+03 copies/μL.
    - **NOTE:** If a non-transduced control sample is unavailable, add 2 µL of VSV-G DNA Control at 1.00E+03 copies/µL to 12.5 µL of transduced sample.
    - NOTE : The PCR inhibition spike control samples will contain a final volume greater than 50 μL due to the extra spike volume.
  - 6.5.3 Dilute an aliquot of VSV-G DNA Control in 1X TE over a range of concentrations from 1.00E+06 to 1.00E+01 copies/μL. Add 12.5 μL of each VSV-G DNA Control standard curve point to the appropriate wells of the plate.
- 6.6 Seal the plate with adhesive film. The plate may be kept at roomtemperature for up to 30 minutesprior to amplification on a qPCR instrument.
- 6.7 Place the sealed plate into the qPCR instrument.
- 6.8 Program the qPCR instrument to detect a FAM reporter and Non-fluorescent quencher with the following amplification thermal profile:

Temperature	Time (Min:Sec)	Repetition
50°C	02:00	N/A
95°c	10:00	N/A
95°c	00:15	10
60°C	01:00	40

- 6.9 Upon completion of the qPCR run, unload and discard the plate. **Do not remove the adhesive cover to prevent laboratory contamination with the target amplicon.** Retain and reuse the compression pad.
- 6.10 Save, analyze and print the run results using a manual Rn threshold of 0.2.

This procedue is made available through federal funds from the National Cancer Institute, NIH, under contract

**NOTE**: Use of other threshold values may be reported with PA Supervisor approval.

## 7.0 Validity of the Assay

- 7.1 Record the assay control results on **Form 22209-03**.
- 7.2 For the assay to be valid, the following criteria must be met:
  - 7.2.1 The slope of the standard curve must be  $\geq$  -3.7 and  $\leq$  -3.1.
  - 7.2.2 The R2 for the standard curve must be  $\geq$  0.990.
  - 7.2.3 For the Extraction Positive Control,  $\geq 10\%$  of the spiked DNA must be recovered.
  - 7.2.4 For the Inhibition Control,  $\geq 10\%$  of the inhibition spike must be detected.
- 7.3 If the run is valid, proceed to Step 7.1 to determine whether the sample is positive or negative for replication competent lentivirus.
- 7.4 If the run is invalid, the assay may be repeated with supervisor approval and indicate in the comments section of **Form 22209-03**.

## 8.0 Analysis of Results

- 8.1 Assay results are reported as either "Positive" or "Negative" for the presence of replication competent lentivirus on **Form 22209-04**.
- 8.2 Record the qPCR data on **Form 22209-05**.
- 8.3 A sample replicate is reported RCL-positive if at least three out of five wells are positive; if less than three wells are positive, the sample is reported RCL-negative.
- 8.4 If both sample replicates are positive, the sample is reported RCL-positive.
- 8.5 If both sample replicates are negative, the sample is reported RCL-negative.
- 8.6 If one replicate is positive and one replicate is negative, the assay must be repeated (indicate in the comments section of **Form 22209-04**).
  - 8.6.1 If both replicates in the repeat assay are negative, the sample is reported RCLnegative.
  - 8.6.2 If one or both replicates in the repeat assay are positive, the sample is reported RCL-positive.

## 9.0 Documentation

Record all reagent part numbers, lot numbers and expiration dates on **Form 22209-01** (Attachment 1). Record the 96-well plate layout and master mix preparation volumes on **Form 22209-02** (Attachment 2). Record control results on **Form 22209-03** (Attachment 3). Record sample analysis results on **Form 22209-04** (Attachment 4). Record RCL-qPCR data on **Form 22209-05** (Attachment 5).

## 10.0 References and Related Documents

- 10.1 SOP 22212 Purification of DNA Using the DNeasy Blood and Tissue Kit
- 10.2 **SOP 22972** Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions
- 10.3 **SOP 22946** The Operation of the NanoDrop ND1000 Spectrophotometer

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

11.1	Attachment 1	Form 22209-01, Quantitative PCR Reagents
11.2	Attachment 2	Form 22209-02, 96-Well qPCR Plate Key, Master Mix Preparation and Dilution of Standard
11.3	Attachment 3	Form 22209-03, RCL-qPCR Assay Validity Criteria Analysis
11.4	Attachment 4	Form 22209-04, RCL-qPCR Assay Sample Analysis
11.5	Attachment 5	Form 22209-05, RCL-qPCR Data

## Attachment 1 Form 22209-1, Quantitative PCR Reagents

FNLCR, BDP Form No.: 22209-01 SOP No.: 22209 Revision 00: MAY 24 2019

**Quantitative PCR Reagents** 

Reagent	Part Number	Lot #	Expiration
TaqMan Universal PCR Master Mix (2X)			
Clinical Center VSVG Forward Primer			
Clinical Center VSVG Reverse Primer			
CC VSVG Probe			
96-Well Optical Plates			
Optical Adhesive Covers			
Water/Buffer (type)			
Standard Control DNA			
1X TE			

Instrument	MEF Number	Calibration/PM Due Date
qPCR Instrument		

Comments:

Form Completed By:\_\_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_\_Date: \_\_\_\_

Reviewed and Approved By:\_\_\_\_\_

QC Test	Request #	
Page	of	

Date:

Initial / Date\_\_\_\_\_

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

#### Attachment 2 Form 22209-2, 96-Well qPCR Plate Key

Form No.: 22208-02 SOP No.: 22209 Revision 00: MAY 24 2019

96-Well qPCR Plate Key

28	12	2	3	4	5	6	7	8	9	10	15	12
A		o										
в												
с												
D												
E				2				0		o		с. 
F												2 8
a												
H												

Master Mix Preparation					
Reagent	Volume for 1 Reaction (µl)	(Number) Reactions (ul)			
Water	10				
TaqMan Universal PCR Master Mbr (2X)	25	4			
Clinical Center VSVG Forward Primer (10 µM)	1				
Clinical Center VSVG Reverse Primer (10 µM)	1				
Clinical Center VSVG Probe (FAM/QSY) (10 µM)	0.5				
Template	12.5				
Total Volume	50				

	Diluti	on Series of Sta	ndard	131
Dilution	Starting Concentration of Standard	Volume of Standard (µl)	Volume of Diluent (µl)	Final Concentration of Standard
1	5	8		8
2		0 20		15 E
3				1
4		1		1
5				1.00E+06
6				1.00E+05
7				1.00E+04
8		2		1.00E+03
9		02		1.00E+02
10		1		1.00E+01

Pipette(s) Used:	MEF#/Sorial#	Calibration Due Date
32523		
2		

Date:

enformed by:

viewed by:

QC Test Request #\_\_\_\_\_ Page\_\_\_\_\_ of \_\_\_\_\_ Initial / Date\_\_\_\_\_

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

#### Attachment 3 Form 22209-3, RCL-qPCR Assay Validity Criteria Analysis

rundig orderin villeli 213	
R^2 Value	
R*2 Value	
Extraction Negative Control Result (Pass/Fail):	
R*2 value result is Pass if the val	ue is ≥ 0.990.
Extraction Positive Control Results	
(Pre-Extraction Spike)	
Spike copy number used	
Copy number recovered	
Average recovery %:	
Extraction Positive Control Result (Pass/Fail):	
Extraction Positive Control result is pass if the amount recov	ered is ≥10%.
	R*2 Value      R*2 Value      Extraction Negative Control Result (Pass/Fail):      R*2 value result is Pass if the value <u>Pre-Extraction Positive Control Results</u> (Pre-Extraction Spike)      Spike copy number used      Copy number recovered      Average recovery %:      Extraction Positive Control Result      Extraction Positive Control Result      Extraction Positive Control result is pass if the amount recovered

Initial / Date\_\_\_\_

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

### Attachment 4 Form 22209-4, RCL-qPCR Assay Sample Analysis

		RCL-qPC	R Assay Sample A	narysis		
	Replicate A	nalysis		Sample	Analysis	
	Result	# Positive Wells per 5 Total Wells		Result	# Positive Replicates	
	Positive	≥3	1 Г	Positive	2	
	Negative	<3	t r	Negative	0	
	and the later		· [	Repeat*	1	
s this a repeated run (Yes	/No)?			"If 1 out of 2 replicates is the sample	positive in the repeated run, e is Positive.	
Sample ID	QC Test Request #	Replicate #	# of Positive Wells	Replica (Positive	te Result /Negative)	Sample Result (Positive/Negative)
				(r contro	in guare)	(i ostatoriograno)
			1			
ienta:						

QC Test Request #\_\_\_\_\_

Page\_\_\_\_\_ of \_\_\_\_\_ Initial / Date\_\_\_\_\_

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract

#### Attachment 5 Form 22209-5, RCL-qPCR Data

FNLCR, BDP Form No.: 22209-05 SOP No.: 22209 Revision 00: MAY 24 2019

RCL-qPCR Data

Sample ID: \_\_\_\_\_

Paste qPCR Data Below:

Comments:	
Form Completed By:	Date:
Reviewed and Approved By:	Date:
	QC Test Request #_
	Page of
	Initial / Date

Page 11 of 11

This procedure is made available through federal funds from the National Cancer Institute, NIH, under contract