

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	[REDACTED]	<input type="text"/>	<input type="text"/>	1 of 17




Biopharmaceutical Development Program (BDP)
 National Cancer Institute at Frederick
 SAIC-Frederick, Inc., P.O. Box B
 Frederick, MD 21702-1201

Master Production Record Approval	
Author Approval: _____	Date: _____
Purification Manager Approval: _____	Date: _____
Project Scientist Approval: _____	Date: _____
Biopharmaceutical Quality Assurance (BQA) Approval: _____	Date: _____
Comparison of Copy to Master Document	
This document is an accurate reproduction of MPR-P- [REDACTED] section [REDACTED] , as found in the Master Document File.	
Checked by: _____	Date: _____
Post-Manufacturing Document Review	
This completed production record has been reviewed and has been found to be complete, correct, and in conformance with relevant SOP's and other documents.	
Reviewed by: _____	Date: _____
Quality Assurance: _____	Date: _____

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>		<input type="text"/>	<input type="text"/>	2 of 17

Document the personnel involved in the production process in the table below.

Operator (Print Name)	Signature	Initials

Master Production Record (MPR)


Production of [redacted]

Section [redacted] : Bulk Aliquot, Sampling and Storage of [redacted]

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P-[redacted]	[redacted]	[redacted]	[redacted]	[redacted]	[redacted]	3 of 17

1.0 Production of [redacted]

Scheduled Tests



Critical Parameters: [redacted]

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	4 of 17

1.1 Distribution of Product

1.1.1 Distribution must take place in a class 100 area. Prior to distributing the “”, verify that EM was performed in the room/BSC where filling will take place. Assemble settling plates in the back two corners of the BSC for environmental monitoring during the fill process.

Room Number: _____ BSC BDP #: _____
EM Performed (Y/N): _____ QC Request # for EM: _____
Settling Plates Part #: _____ Lot #: _____ Exp. Date: _____

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

1.1.2 The “”, from section , will be distributed into tared containers with a fill volume of approximately mL per container for storage at °C until further processing. Weigh the bag of “” and calculate the number of L containers needed. Round the containers needed to the nearest whole number. Record the requested information in the table below. . Include the balance printouts as Attachment .

Bulk Aliquot, Sampling and Storage Step Start Date/Time: _____
Balance BDP #: _____
Net Wt. of <input type="text"/> = Gross Weight of <input type="text"/> : _____ g – Tare Weight of <input type="text"/> L Bag, sec <input type="text"/> ,: _____ g = _____ g
Number of <input type="text"/> Containers Needed = Net Wt. of <input type="text"/> g / <input type="text"/> mL = _____ containers

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	5 of 17

1.1.3 Using the tubing on the bag that was not used during sterile/final filtration, aseptically distribute, in a class 100 area, the product into the calculated number of tared ████████ containers. Weigh each container after transfer and calculate the net weight. Label each container with “████████”, Container Number, Fill Volume, Date, For Further Processing”. Include the balance printouts as Attachment ████████.

Balance BDP #: <input type="text"/> ████████ Container Part #: <input type="text"/> Lot #: <input type="text"/> Expiration Date: <input type="text"/>
Container # 1: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 2: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 3: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 4: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 5: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 6: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 7: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 8 Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 9: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g
Container # 10: Net Wt. = Gross Wt. : <input type="text"/> g –Tare Wt. of Container: <input type="text"/> g = <input type="text"/> g

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	[REDACTED]	<input type="text"/>	<input type="text"/>	6 of 17

Container # 11: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 12: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 13: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 14: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 15: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 16: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 17: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 18: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 19: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g
Container # 20: Net Wt. = Gross Wt. : _____ g –Tare Wt. of Container: _____ g = _____ g

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	[REDACTED]	<input type="text"/>	<input type="text"/>	7 of 17

1.1.4 Upon completion of the distribution of the product, individuals that performed the operation must perform personnel monitoring by submitting touch plates to QC. Submit the settling plates assembled in Section [REDACTED] to QC for environmental monitoring. Include the QC Request Verification for personnel and environmental monitoring as Attachment [REDACTED].

Persons Touch Plated (Print Names): _____ QC Request # for Personnel Monitoring: _____ QC Request # for EM: _____
--

Performed By:	Date:	Verified By:	Date:
---------------	-------	--------------	-------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	8 of 17

Comments:

Reviewed By Purification Manager:	Date:
-----------------------------------	-------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	9 of 17

1.2 Sampling and Sample Distribution

1.2.1 Aseptically remove samples from ONE of the containers filled in Section ██████████ using a new pipette with each sampling. Distribute the samples as described in the table below in a class 100 area. Label the samples as per Section ██████████ of this BPR.

NOTE: Sterility samples, if applicable, must be removed in equal portions from each of the containers filled in Section ██████████, i.e. if there are five containers, remove 2 X 2mL from each container, using a new pipette each time.

Pipette BDP Lot #: _____ BDP Part #: _____ Expiration Date: _____			
Container Used for Sampling of Non-Sterility Samples: _____ Sterility Samples Taken (Y/N): _____			
Volume of Sample Removed from Each Container for Sterility: _____ mL			
Number of Samples	Sample Vol., mL	Sample Volume Removed, mL	Purpose
██████████	██████████	██████████	Process Retains
██████████	1.0 ± 0.1mL Glass Vial	██████████	Corporate Retains
██████████	██████████	██████████	Assays in Section ██████████
██████████	██████████	██████████	Assays in Section ██████████

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	10 of 17

1.2.1.1 If the QC samples and/or retain samples, from section ██████████, are temporarily stored until submitted to QC/MMIC, record the requested information in the table below.

Refrigerator/Freezer BDP #: _____	Temperature: _____ °C
Retain Sample Stored (Y/N): _____	Date of Storage, if applicable: _____
QC Sample Stored (Y/N): _____	Date of Storage, if applicable: _____

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

1.2.2 Submit the samples (Part # ██████████) to QC for testing and enter the QC test request numbers in the table below. Store the samples at ██████████ °C (if applicable, sterility sample must be stored at 2-8°C). Include the QC Request Verification as described in the table below.

Test	SOP Number	QC Test Request #	Sample Volume, mL	Attachment #
██████████	██████████		██████████	██████████
██████████	██████████		██████████	██████████
██████████	██████████		██████████	██████████
██████████	██████████		██████████	██████████

Submitted By: _____	Date: _____
---------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	11 of 17

1.2.3 Submit the retain samples to MMIC for storage at °C and include a copy of the Sample Input Form as Attachment .

Refrigerator/Freezer ID number: _____ Cal. Exp. Date: _____	
Number of Samples	Sample Volume, mL
<input type="text"/>	<input type="text"/>

Submitted By: _____	Date: _____
---------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	12 of 17

1.3 Weighing/Labeling of "██████████"

1.3.1 Weigh each of the **sampled containers** of "██████████" and calculate the final volume. Include the balance printouts as Attachment ██████████.

Balance BDP #: _____	
Net Wt. ██████████ Container #1 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #2 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #3 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #4 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #5 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #6 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #7 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #8 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #9 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g
Net Wt. ██████████ Container #10 = Gross Wt.: _____ g – Tare Wt. ██████████ L PETG, Sec. ██████████ :	_____ g = _____ g

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	13 of 17

Net Wt. <input type="text"/>	Container #11 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #12 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #13 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #14 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #15 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #16 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #17 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #18 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #19 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g
Net Wt. <input type="text"/>	Container #20 = Gross Wt.: _____ g – Tare Wt. <input type="text"/>	L PETG, Sec. <input type="text"/>	: _____ g = _____ g

Performed By:	Date:	Verified By:	Date:
---------------	-------	--------------	-------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	[REDACTED]	<input type="text"/>	<input type="text"/>	14 of 17

1.3.2 Request for a label to be printed for each container of the "[REDACTED]". The label must contain, at a minimum, the information below. Include a copy of the actual label(s) as Attachment [REDACTED].

Note: If the product is being shipped to an outside vendor, the label must be requested from QA and a copy of the label galley must be included as Attachment [REDACTED].

Product Name: [REDACTED]
BDP Lot # _____
Store At: [REDACTED] °C Container Volume: _____ mL
Buffer: _____
Concentration: _____ mg/ml Fill Date: _____
Container # _____ of _____
** FOR FURTHER MANUFACTURING USE ONLY **
BDP NCI-Frederick

Performed By:	Date:	Verified By:	Date:
Purification Manager Approval:			Date:

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	15 of 17

1.3.3 Apply ONE label to each of the containers of "██████" and record the requested information in the table below.

Note: If the label, from section ██████, was requested from QA, include a copy of the label reconciliation to this BPR as Attachment ██████.

Room # for Labeling: _____
Container Cracks/Leaking Observed (Y/N): _____
Container Deterioration Observed (Y/N): _____
Container Capped and Sealed Properly (Y/N): _____
Label Correctly and Securely Applied (one per Container) (Y/N): _____
Label Information Legible and Correct (Y/N): _____

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

1.4 Storage of "██████"

1.4.1 Store the containers of "██████" at ██████ °C until transferred to MMIC for long term storage. Record the requested information in the table below.

Refrigerator/Freezer BDP #: _____ Temperature: _____ °C
Date/Time of Storage: _____

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>	██████████	<input type="text"/>	<input type="text"/>	16 of 17

1.4.2 Transfer the “” containers to Materials Management and Inventory Control for storage at °C. Include a copy of Form 20303-01, “MMIC CGMP Manufacturing Product Inventory”, as Attachment .

NOTE: Prior to transferring the containers to MMIC, it is ideal to individually shrink wrap each of the containers.

Shrink Wrap Part #: _____	Lot #: _____	Exp. Date: _____
Refrigerator/Freezer BDP #: _____	Calibration Exp. Date: _____	Temperature: _____ °C
Date/Time of Storage: _____		

Submitted By: _____	Date: _____	Received By: _____	Date: _____
---------------------	-------------	--------------------	-------------

1.4.3 Record the end date/time of the Bulk Aliquot, Sampling and Storage step. Calculate the duration of the entire step, using the step start date/time from section , in the table below.


Step Start Date/Time, sec. <input type="text"/> : _____	Step End Date/Time: _____	Duration: _____ hrs/min.
---	---------------------------	--------------------------

Performed By: _____	Date: _____	Verified By: _____	Date: _____
---------------------	-------------	--------------------	-------------

Master Production Record (MPR)

Production of

Section : Bulk Aliquot, Sampling and Storage of

Document No.	Project No.	Lot No.	Location	Revision	Effective	Page
MPR-P- <input type="text"/>	<input type="text"/>	<input type="text"/>		<input type="text"/>	<input type="text"/>	17 of 17

Comments:

Reviewed By Purification Manager:	Date:
-----------------------------------	-------