



BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

TABLE OF CONTENTS

1. PURPOSE	1
2. SCOPE	1
3. RESPONSIBILITIES	1
4. SAFETY	2
5. MATERIALS AND REAGENTS	2
6. EQUIPMENT	2
7. PROCEDURE	2
8. VALIDITY OF RESULTS	6
9. ANALYSIS OF RESULTS	7
10. DOCUMENTATION AND RECORDS	8
11. REFERENCES AND RELATED DOCUMENTS	8

1. PURPOSE

This procedure describes the materials and methods for quantitative PCR (qPCR) amplification and detection of various species of Mycoplasma that may be present in a sample using the MycoSEQ Mycoplasma Real-Time PCR Detection Kit.

2. SCOPE

The use of the MycoSEQ Mycoplasma Real-Time PCR Detection Kit as described in this SOP provides the ability to detect more than 90 species, including *M. genitalium*, *M. pirum*, *Acholeplasma* and *Spiroplasma*.

3. RESPONSIBILITIES

3.1 Director, Process Analytics/Quality Control (PA/QC)

- Defines the procedure.

3.2 Process Analytics/Quality Control (PA/QC)

- Trains Laboratory personnel.
- Performs the procedure.
- Reviews data.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

3.3 Biopharmaceutical Quality Assurance (BQA)

- Provides quality oversight.

4. SAFETY

Treat all samples suspected or known to contain Mycoplasma per **SOP 26101 - Labeling, Transport, Submission, Storage, and Handling of Biohazardous Materials Within the BDP.**

5. MATERIALS AND REAGENTS

Part Number	Description	BDP Approved Substitution Permitted?
31167	MycoSEQ Mycoplasma Real-Time PCR Detection Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
21141	MicroAmp Optical 96 Well Reaction Plates	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
21142	MicroAmp Optical Adhesive Film	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
30267	1X TE	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
10189	Distilled Water, DNase Free RNase Free	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
30007	1X PBS	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO

6. EQUIPMENT

- QuantStudio 6 qPCR System
- Plate Centrifuge

7. PROCEDURE

7.1 Determine the sample composition to determine whether a nucleic acid extraction is required.

7.1.1 Previously purified DNA in low salt (< 100mM) typically does not require extraction prior to amplification.

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

- 7.1.2 Samples containing concentrated protein, intact cells and/or high salt or other PCR inhibitors must be extracted prior to amplification.
- 7.2 An extraction negative control should be performed in which 200 μ L of 1X PBS are extracted in parallel with the test samples, as described in **SOP 22212 - Purification of DNA Using the DNeasy Blood and Tissue Kit** or **SOP 22972 – Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions**.
- 7.3 An extraction spike control may be performed to control for the potential loss of target mycoplasma DNA during the sample extraction step. The MycoSEQ Discriminatory Positive/Extraction Control (or BDP-approved equivalent) may be used to spike a buffer blank sample or 1X PBS as described below.
 - 7.3.1 Completely thaw the MycoSEQ Mycoplasma Real-Time PCR Detection Kit Box 2 MycoSEQ Discriminatory Positive/Extraction Control, 1,000 copies/ μ L.
 - 7.3.2 Add a 5 μ L aliquot of MycoSEQ Discriminatory Positive/Extraction Control to 195 μ L of buffer blank sample or 1X PBS.
 - 7.3.3 Extract the extraction control sample in parallel with the test samples, as described in **SOP 22212 - Purification of DNA Using the DNeasy Blood and Tissue Kit** or **SOP 22972 – Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions**.
- 7.4 Extract all test samples and controls according to **SOP 22212 - Purification of DNA Using the DNeasy Blood and Tissue Kit** or **SOP 22972 – Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions**.
- 7.5 qPCR Master Mix Preparation – Perform in the “DNA-free” PCR Set-up Laboratory.
 - 7.5.1 Completely thaw the MycoSEQ Mycoplasma Real-Time PCR Detection Kit reagents.
 - 7.5.2 Record the reagent and material part numbers, lot numbers and expiration dates as well as any relevant additional comments on **Form 22208-01**.
 - 7.5.3 Vortex the reagents.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

7.5.4 Prepare a Premix Solution as follows.

Component for Premix Solution	Volume for One 30- μ L Reaction
Power SYBR Green PCR Master Mix, 2X	15.0 μ L
Mycoplasma Real-Time PCR Primer Mix, 10X	3.0 μ L
Water	2 μ L
Total Premix Solution Volume	20.0 μ L

7.5.5 Prepare sufficient volume of Premix Solution to account for all test samples (3 replicate wells) and controls (2 or 3 wells) with 10% excess to account for losses due to repeated pipetting.

7.5.6 Record the volumes used on **Form 22208-02**.

7.5.7 Mix the Premix Solution by pipetting up and down.

7.5.8 Add 20 μ L of Premix Solution to each well of the 96-well PCR plate that will be used.

7.5.9 Kit Negative Control (water): Add 10 μ L of Negative Control (water) to each of 2 wells.

7.6 Sample addition to PCR plate – Perform in the Positive Control Set-up Laboratory

7.6.1 To the Premix Solution dispensed into each well, add the following:

7.6.1.1 Test Sample: Add 10 μ L of unknown sample to each of 3 wells.

7.6.1.2 Inhibition-control Reaction: Add 10 μ L of unknown sample +2 μ L of Discriminatory Positive Control to each of 2 wells.

NOTE: The PCR inhibition spike control samples will contain a final volume greater than 30 μ L due to the extra spike volume.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

- 7.6.1.3 Positive Control Reaction: Add 2 μ L of Discriminatory Positive Control + 8 μ L of Negative Control (water) to each of 2 wells.
- 7.6.1.4 Extraction Positive Control Reaction: Add 10 μ L of extraction positive control sample to each of 3 wells.
- 7.6.1.5 Extraction Negative Control Reaction: Add 10 μ L of extraction negative control sample to each of 2 wells.
- 7.7 Seal the plate with adhesive film. The plate may be kept at room temperature for up to 30 minutes prior to qPCR amplification on the QuantStudio 6 (**SOP 22973 – Operation and Maintenance of the QuantStudio 6 qPCR System**) or a BDP-approved equivalent instrument.
- NOTE:** Perform in the Positive Control Set-up Laboratory.
- 7.8 The sealed plate may be centrifuged at 2,000 xg for 1 minute to remove or displace bubbles from the bottom of the well.
- 7.9 Place the sealed plate into the qPCR instrument ensuring that the plate is appropriately placed in the carriage by lining up the notched corner of the plate to the notched corner of the plate holder.
- 7.10 Select SYBR Green Reagents as the reagent type.
- 7.11 Program the QuantStudio 6 per **SOP 22973 – Operation and Maintenance of the QuantStudio 6 qPCR System**, using the following amplification thermal profile:

Temperature	Time (Min:Sec)	Repetition
95°C	10:00	N/A
95°C	00:15	40
60°C	01:00	
95°C	00:15	Melt Curve
60°C	01:00	
95°C	00:15	

- 7.12 Upon completion of the qPCR run, unload and discard the plate. Do not remove the adhesive cover to prevent laboratory contamination with the target amplicon!

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

7.13 Save, analyze, and print the run results according to **SOP 22973 – Operation and Maintenance of the QuantStudio 6 qPCR System**, using a manual Rn threshold of 0.2.

NOTE: You must click the green “Analyze” button after setting the manual threshold to apply the change to the data.

NOTE: Use of other threshold values may be reported with PA Supervisor approval.

8. VALIDITY OF RESULTS

8.1 Record the assay control results on **Form 22208-03**.

8.2 The assay validity criteria are as follows:

Control	Ct	Tm
Extraction Negative (Blank) Control	≥36.00	<82°C*
No Template Control (NTC)	≥36.00	<82°C*
Extraction Spike Control	<36.00	82°C-86°C
PCR Positive Control	<36.00	82°C-86°C
Inhibition Control	ΔCt < 3	82°C-86°C

***NOTE:** Any detected amplification products with a Tm below 82°C may be ignored. These amplification products represent non-specific amplification or primer-dimers and do not interfere with the final results.

8.3 An “Undetermined” or “No Ct” result is a valid result for an Extraction Negative (Blank) Control and a No Template Control (NTC). An “Undetermined” or “No Ct” result occurs when the amplification plot never crosses the cycle threshold.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

- 8.4 The derivative value (DV) from the melt curve analysis should be greater than background for the Extraction Spike Control, PCR Positive Control and Inhibition Control. Include the melt curve for these controls on **Form 22208-04**.
- 8.5 The Inhibition Control ΔCt is defined as the mean Ct value for the Inhibition control sample(s) minus the mean Ct value for the qPCR Positive Control.
- 8.6 At least one of the extraction positive control qPCR replicates must be lower than Ct 36.
- 8.7 If the run is valid, proceed to Step 9.1 to determine whether the sample is positive or negative for mycoplasma.
- 8.8 If the run is invalid, the assay may be repeated with supervisor approval.

9. ANALYSIS OF RESULTS

- 9.1 Assay results are reported as either "Positive" or "Negative" for the presence of Mycoplasma genomic DNA on **Form 22208-05**.
- 9.2 The criteria for determining whether a sample is Mycoplasma Positive or Negative are indicated in the table below:

Result	Ct	Tm
Positive	<36.00	75°C-86°C
Negative	≥36.00	<75°C

- 9.3 For a sample to be called positive, the Ct value must be <36 and the Tm for the amplification product must be in the range of 75°C to 86°C.
 - 9.4 For a sample to be called negative, the Ct value must be ≥36 or the Tm must be less than 75°C or both.
- NOTE:** An Undetermined or No Ct result is also considered negative.
- 9.5 The derivative value (DV) from the melt curve analysis should be greater than background for a sample to be called positive. Include the melt curve for each sample on **Form 22208-04**.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Method for Detection of Mycoplasma Species Using the MycoSEQ Detection Kit
SOP Number: 22208
Revision: 04

9.6 If the sample is negative, but the inhibition control exhibits a ΔC_t that is ≥ 3 , the qPCR reaction was likely inhibited. Indicate this in the comments section of **Form 22208-05** and re-purify and re-test the sample. If the sample is mycoplasma-positive but the inhibition control exhibits a ΔC_t that is ≥ 3 , the sample is reported as positive for mycoplasma. Indicate this result in the comments section of **22208-05**.

10. DOCUMENTATION AND RECORDS

- 10.1 Record all reagent part numbers, lot numbers and expiration dates on **Form 22208-01**.
- 10.2 Record the 96-well plate layout and master mix preparation volumes on Form 22208-02. Record control results on **Form 22208-03**.
- 10.3 Record Melt Curves on **Form 22208-04**
- 10.4 Record sample analysis results on **Form 22208-05**

11. REFERENCES AND RELATED DOCUMENTS

Document Number	Title
22212	Purification of DNA Using the DNeasy Blood and Tissue Kit
22972	Operation and Maintenance of the MagNA Pure 24 System for Nucleic Acid Extractions
22973	Operation and Maintenance of the QuantStudio 6 qPCR System
26101	Labeling, Transport, Submission, Storage, and Handling of Biohazardous Materials Within the BDP
22208-01	Quantitative PCR Reagents
22208-02	96-Well qPCR Plate Map and Master Mix Preparation
22208-03	MycoSEQ qPCR Assay Validity Criteria Analysis
22208-04	Melt Curve Analysis
22208-05	MycoSEQ qPCR Assay Sample Analysis