



BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

TABLE OF CONTENTS

1. PURPOSE	1
2. SCOPE	1
3. RESPONSIBILITIES	1
4. MATERIALS AND REAGENTS.....	2
5. EQUIPMENT	4
6. SAFETY.....	4
7. PROCEDURE.....	4
8. VALIDITY OF RESULTS	9
9. ANALYSIS OF RESULTS.....	11
10. DOCUMENTATION AND RECORDS.....	12
11. REFERENCES AND RELATED DOCUMENTS.....	13
12. ATTACHMENTS	14

1. PURPOSE

This procedure describes the general materials and methods for quantitative PCR (qPCR) amplification and detection of various nucleic acid sequences that may be present in a test sample

2. SCOPE

This Standard Operating Procedure (SOP) is intended for general qPCR assays that require the characterization and/or quantitation of specific nucleic acid target sequences in a test sample. Multiple qPCR reactions can be performed to elucidate more than one species or target sequences that may be present in a given sample. Specific reagents, primers, and probes necessary for each qPCR reaction are individually listed in the attachments associated with this SOP. Specific SOPs are available for some qPCR methods including Mycoplasma detection, viral contamination, and detection of RCA; the requirements of specific assay SOPs (if applicable) supersede those of this SOP.

3. RESPONSIBILITIES

3.1 Director, Process Analytics/Quality Control (PA/QC)

- Defines this procedure.
- Controls samples submitted for qPCR analysis.
- Reviews and approves result of test analyses.



BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 3.2 PA/QC Laboratory Personnel / Analyst
- Implements this procedure.
 - Follows this procedure.
- 3.3 PA/QC Supervisor
- Approves the use of alternate threshold values.
- 3.4 Quality Assurance
- Provides quality oversight.

4. MATERIALS AND REAGENTS

Part Number	Description	BDP Approved Substitution Permitted?
General qPCR Reagents, as required		
21141	Applied Biosystems, Inc. (ABI) 96-Well Plate	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
21142	ABI Adhesive Optical Plate Cover	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30268	TaqMan Universal PCR Master Mix	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30265	Super Script III RT/Platinum 1-step kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
31244	RNaseOUT	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
31129	Microcentrifuge Tube, Low Binding (1.7 mL)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
31130	Microcentrifuge Tube, Low Binding (0.65 mL)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
31486	Eppendorf Protein Low Binding tubes (5.0 mL)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
Buffer Reagents, as required		
10189	Nuclease-Free Water	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
30007	PBS (pH 7.4)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO



BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Part Number	Description	BDP Approved Substitution Permitted?
30267	TE (1X)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
30295	WFI, USP	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
30266	DEPC Treated Water	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
10042	Tris HCl (pH 7.4)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
DNase Reagents as required		
31243	DNaseI (New England Biolabs)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
10098	EDTA (0.5M, pH 8.0)	<input checked="" type="checkbox"/> YES <input type="checkbox"/> NO
DNA/RNA Extraction Reagents, as required		
30442	QIAGEN DNA Mini Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30444	QIAGEN Blood Mini Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30450	QIAGEN Viral RNA Mini Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30443	QIAGEN DNeasy Blood & Tissue Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30446	QIAGEN ATL lysis Buffer	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30447	QIAGEN AVL lysis Buffer	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
30445	QIAGEN AL lysis Buffer	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
31236	Roche MagNA Pure 24 Total NA Isolation Kit	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
qPCR Reaction Specific Reagent (see Attachments for PNR)		
Assay Specific	Nucleic Acid Standards	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Part Number	Description	BDP Approved Substitution Permitted?
Assay Specific	Forward Primer	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
Assay Specific	TaqMan Probe	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
Assay Specific	Reverse Primer	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
Assay Specific	Positive Reference Control Reagents	<input type="checkbox"/> YES <input checked="" type="checkbox"/> NO

5. EQUIPMENT

- QuantStudio™ 6 qPCR System
- Roche MagNA Pure 24

6. SAFETY

- 6.1 Treat all viral and biohazardous samples per **SOP 26101 - Labeling, Transport, Submission, Storage, and Handling of Biohazardous Materials within the BDP.**
- 6.2 All viral samples must be inactivated in QIAGEN buffer ATL or AL or by MagNA Pure extraction prior to use.

7. PROCEDURE

- 7.1 Determine the sample composition in order to determine whether a nucleic acid extraction is required. Refer to the assay-specific attachments for a list of the approved extraction kits to be used with each amplicon.
- 7.1.1 Previously purified DNA or RNA (such as a plasmid preparation) in low salt (< 100 mM) and cleaning swabs typically do not require extraction prior to amplification.
- 7.1.2 Samples containing concentrated protein, intact cells, and/or high salt or other PCR inhibitors (including some culture medias) will typically require extraction prior to amplification.
- 7.2 To control for potential loss of target DNA/RNA during the sample extraction step, perform an extraction spike control with either the control target DNA/RNA or a heterologous control target using a second qPCR amplicon.

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 7.2.1 To determine the extraction and DNA recovery efficiency, separately extract the appropriate requestor or PA/QC-provided buffer blank sample after adding a spike within the assay's standard control range to the buffer sample aliquot.
- 7.2.2 Record the amount (in the final reaction) and type of spiked DNA used on PCR Inhibition Spike Results or/and Sample Recovery Spike results on the Form 22195-03.
- 7.3 Determine the appropriate qPCR or RT-qPCR amplification reagent kit to use. Refer to the assay-specific attachments for the approved amplification reagent kits that are to be used with each amplicon.
- 7.4 Determine the appropriate amount of qPCR Master Mix or RT-qPCR Master Mix to prepare as follows.
 - 7.4.1 Standard Curve reactions are performed in duplicate (R&D, cleaning, and in-process testing) or triplicate (product release testing only).
 - 7.4.2 No Test Control (NTC) – performed in single, duplicate, or triplicate reactions with nuclease-free (NFW), Diethyl pyrocarbonate (DEPC) treated, or WFI water as the control sample. Perform triplicate NTC reactions for release testing.
 - 7.4.3 Sentinel controls – optional water or buffer controls to monitor intra-assay and QC test contamination levels identical in setup to that of the NTC control and located near the corners of a 96-well reaction plate. Sentinel controls are most useful with samples known or expected to contain high target copy numbers.
 - 7.4.4 Extraction efficiency (spike) control – performed in duplicate or triplicate reactions. Spike reactions can utilize either the standard DNA/RNA target or a heterologous amplicon using a secondary reporter dye / detector channel and a non-fluorescent quenching dye.

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 7.4.5 PCR Inhibition (post-extraction spike) control – performed in single or duplicate reactions with representative test sample(s) that may be expected to have the greatest level of PCR inhibition. For example, a column load can be expected to have a greater probability of exhibiting inhibition than a column eluate sample and samples containing high salt (≥ 0.25 M), Phosphate (≥ 50 mM) and/or Ethylenediaminetetraacetic acid (EDTA) (≥ 25 mM) concentrations may experience PCR inhibition even after extraction. Where feasible, test each submitted sample for the presence of PCR inhibition. Spike reactions can utilize either the standard DNA/RNA target or a heterologous amplicon using a secondary reporter dye / detector channel and a non-fluorescent quenching dye.
- 7.4.6 Test sample reactions – duplicate or triplicate reactions for each test sample and/or sample dilution to be assayed. Product release tests must utilize triplicate reactions, while cleaning, in-process, and intermediate release samples may utilize duplicate sample replicates.
- 7.4.7 It is advantageous to test the highest possible concentration of sample not exhibiting inhibition in order to obtain the optimal assay sensitivity. The expected target concentration must be no more than 3.33 Ct lower than the standard curve's upper limit in order to maintain quantitative accuracy. Ct scores < 12 are progressively less accurate and should be avoided. DNA concentrations > 100 ng/rxn may inhibit qPCR. Nucleic acid concentrations > 500 ng/rxn will inhibit qPCR.
- 7.4.8 Total the number of individual reactions necessary for the assay, multiply by the amount of Master Mix and other reagents required (see 7.5 below) per reaction and prepare at least 110% this amount for aliquoting. Record the number of reactions required on Form 22195-02.
- 7.5 Prepare the calculated amount of qPCR or RT-qPCR Master Mix in the “DNA free” PCR set-up area () by adding sufficient reagents to generate the concentrations and volumes per individual reaction listed per the amplicon specific attachment. A general 50 μ L reaction design for a qPCR is provided in the table below.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

REAGENT	VOLUME PER REACTION IN μL
Universal 2X Master Mix	25 μL (1x final conc.)
DEPC-treated, Nuclease-free, or WFI Water	Assay Specific (Bring to a 30 μL volume after addition of primers / probe.)
Primers (Fwd and Rev)	Assay Specific (typically 5-20 pmol)
Probe	Assay Specific (typically 5-20 pmol)
Reagent Mix Volume	30 μL
Volume of Sample to be added	20 μL
Total Reaction Volume	50 μL

NOTE: The PCR inhibition spike control samples will contain a final volume greater than 50 μL due to the extra spike volume unless the spike volume is substituted for water in the master mix.

- 7.6 Ensure that the proper final pmol concentration of primers and probe is added per reaction based on the assay-specific attachment. The amplicon-specific primer and probe pmol concentrations and other indicated reagent volumes take precedent over the volume figures above as necessary.
- 7.7 Record the reagent and material part numbers, lots, and expiration dates as well as any additional comments on **Form 22195-01**. Record the volumes used per reaction on **Form 22195-02**.
- 7.8 Load the 96-well plate with prepared Master Mix in a dedicated “DNA free” PCR set-up area (master mix setup room) according to the desired 96-well qPCR plate key (Form 22195-02). Complete the NTC controls by adding NFW or WFI to the NTC wells. Temporarily seal or cover the 96-well plate for transport to the “DNA allowed” final set-up area. Plate transport between laboratories must be between +4 °C (on ice) and room temperature.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 7.9 In a dedicated PCR set-up area (PCR plate and DNA sample setup room), load the standard, test, sentinel (if needed), and spike control samples to their appropriate wells as designated on the 96-well qPCR plate key (**Form 22195-02**).
- 7.10 Seal the plate with an optical adhesive cover and ensure that all bubbles are removed from the samples in each well, by gently tapping the plate several times if necessary. The plate may be kept at room temperature for up to 30 minutes prior to amplification on the Applied Biosystems QuantStudio™ 6.
- 7.11 Place the sealed plate into the qPCR instrument, ensuring that the plate is appropriately placed in the carriage by lining up the notched corner of the plate to the notched corner of the plate holder.
- 7.12 Perform qPCR amplification according to **SOP 22973 – Operation and Maintenance of the QuantStudio™ 6 qPCR System** using the thermocycler parameters listed in the appendix for the specific qPCR assay to be used.
- 7.13 Refer to the amplicon-specific requirements (see assay specific attachment) for cycling profiles. A typical qPCR reaction profile is shown below.

STAGE	TEMP	TIME (MIN:SEC)	REPETITION
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	0:15 1:00	40 cycles

- 7.14 A typical RT-qPCR reaction profile is shown below.

STAGE	TEMP	TIME (MIN:SEC)	REPETITION
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C 60°C	0:15 1:00	40 cycles

- 7.15 Since most qPCR reactions utilize the same reaction profile, more than one amplicon can be on a plate. This is especially useful when a test sample requires multiple qPCR assays for different nucleic acid targets.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 7.16 Prior to starting the qPCR run, ensure that the assay file contains the appropriate Reporter dyes and cycle conditions for each assay performed on the plate.
- 7.17 Upon completion of the qPCR run, unload and discard the plate. For viral samples, treat the sealed plate as potentially infectious, biohazardous waste. Do not remove the adhesive cover to prevent laboratory contamination with the target amplicon!
- 7.18 Save, analyze, and print the run results according to **SOP 22973 – Operation and Maintenance of the QuantStudio™ 6 qPCR System**, using the Rn Threshold value suggested or indicated for the specific amplicon used (0.2 Rn is typical). Use of other threshold values, including software calculated threshold values, may be reported with PA/QC Supervisor approval. Software calculated Rn Threshold values are typically used with R&D samples and especially with samples containing large amounts of the target sequence.

8. VALIDITY OF RESULTS

- 8.1 Assay results are considered valid if:
 - 8.1.1 The absolute standard curve r^2 value should be ≥ 0.945 . The low copy (or low mass) point of the standard curve may be omitted from the best fit calculation at the PA/QC Supervisor's discretion. Document low copy (mass) point omission in the comments section of **Form 22195-03**.
 - 8.1.2 For release tests using 3 reactions per standard point, no more than 2 other outliers on the standard curve may be dropped using the appropriate Dixon-Q outlier test (refer to Form 22195-04) with a one-sided alpha value of 5%. Only replicates that meet or exceed the appropriate Dixon-Q P-score (at alpha 0.05) may be considered outliers and be dropped. Putative outliers that do not meet or exceed the P-score value cannot be dropped from the standard curve. Only a single replicate Ct may be dropped from any given point on the standard curve. If more than 2 replicates must be dropped from the standard curve prior to the best fit calculation, excluding the low copy point replicates, the assay is considered invalid. Dixon-Q scoring is not required for R&D, FIO, cleaning, EM, and In-process testing, due to the increased sample replicate variability inherent in these samples and the use of two standard replicates (rather than 3) per point.

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

- 8.1.3 The Dixon-Q one-sided outlier test is performed as follows. For the affected standard point containing the putative outlier replicate, rank the point's replicate Ct values from lowest to highest. If the putative outlier Ct falls below the other replicate Ct values apply the Low Value Tau calculation (i.e., the outlier is the lowest ranked replicate), if the putative outlier falls above the other replicate Ct values apply the High Value Tau calculation (i.e., the outlier is the highest ranked replicate).

Low Value Tau calculation: $(\text{second lowest replicate Ct} - \text{outlier Ct}) / (\text{highest replicate Ct} - \text{outlier Ct})$.

High Value Tau calculation: $(\text{outlier Ct} - \text{second highest replicate Ct}) / (\text{outlier Ct} - \text{lowest replicate Ct})$.

The tau value obtained is compared to the P-score for the number of replicates tested at $\alpha = 0.05$ (5%, one-sided test). Only tau values \geq the P-score are considered genuine outliers and may be dropped from the standard curve.

For example, at 3 replicates per point, the one-sided P-score for $\alpha 0.05$ (5%) is 0.941 (also see current USP). Only outliers confirmed to have tau-values ≥ 0.941 using the appropriate Dixon-Q calculation may be dropped from the standard curve best fit analysis.

Record the Dixon-Q analysis on **Form 22195-04**.

- 8.1.4 The average Ct score of the No Test Control(s) (NTC) must be ≥ 33 Ct (> 38 is typical). Values below 33 Ct indicate excessive cross-contamination on the plate or in the laboratory. Individual NTC Ct scores should not be dropped prior to averaging. Most amplicons will have NTC Ct scores < 33 Ct when a large amount (typically $> 1 \times 10^9$ copies or > 100 ng per reaction) of target is present on the test plate, such low NTC Ct scores only invalidate a release assay. R&D, FIO, and in-process tests with NTC Ct scores < 33 Ct can be invalidated, or allowed with the QC Supervisor's discretion, when explained in the comments section of Form 22195-03 or elsewhere in the report.
- 8.1.5 The Extraction spike control (if required) indicates that $\geq 10\%$ of the spiked DNA was recovered. Test samples that unexpectedly have large amounts of target DNA/RNA that substantially exceeds ($\geq 10X$) the amount of spiked control DNA will not be considered to invalidate the assay.

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

8.1.6 The PCR Inhibition control (if required) indicates that $\geq 10\%$ of the inhibition spike DNA is recovered when compared to the otherwise identical but unspiked sample. For samples containing large amounts of target sequence (i.e., containing $\geq 10X$ the amount of DNA spiked), PCR inhibition will be considered to have occurred only if the spiked control Ct score is > 3 Ct greater than the unspiked sample Ct score.

8.2 All assay reactions will be reviewed for indications of false amplification, PCR inhibition, cross-plate contamination, or other unusual effects via the Amplification Plot report page.

9. ANALYSIS OF RESULTS

9.1 The final assay results in copies or mass per reaction (rxn) are typically converted to reflect the results in copies/mL or mass/mL.

Calculation Method: $\left(\left(\left[\text{Average Sample copy/rxn based on Assay Standard Curve and averaged sample Ct} \right] - \left[\text{Average NTC copy number} \right] \right) / \left[\text{Sample Volume } \mu\text{L per rxn} \right] \times 1000 \mu\text{L/mL} \right)$.

This calculation is performed by Quant Studio Software and the result can be exported as excel sheet.

NOTE: This formula does not capture any relative difference in extraction to elution volume nor does it capture dilutions made to the sample by PA/QC or the test requestor prior to qPCR analysis. Any such dilutions to the sample performed by PA (such as on **SOP 22972 MagNAPure DNA extraction** or **SOP 23113 QIAamp DNA/RNA extraction**) must be incorporated in the final result calculation as well. Use the **Form 22195-05** for quantitative PCR data.

9.1.1 The Average NTC copy number is subtracted from the average test sample copy number values.

9.1.2 The Extraction recovery efficiency (if determined) is set as 100% if the recovery is $\geq 10\%$ with the corrected sample copy or mass/rxn value.

9.1.3 Any effective sample dilution or concentration made during the extraction (and subsequent elution) step are calculated and adjusted (i.e., if 400 μL of sample is extracted and a volume of 200 μL is eluted, the effective result is a 2X concentration of sample).

9.1.4 The volume per reaction is normalized to 1 mL by multiplying by 50 assuming 20 μL of test sample are used per reaction (20 μL sample per reaction $\times 50 = 1$ mL).

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

9.1.5 Any sample dilutions made by PA/QC prior to amplification are corrected for errors.

9.1.6 The use of viral lysis buffers necessitates the appropriate correction factor to the final copy number/mL based on sample volume used and elution volume.

9.2 For product release tests, the 95% Confidence Interval (CI) for the test sample replicates is calculated using the formula:

$$95\% \text{ CI} = \text{Average replicate value} \pm 1.96 (\sigma \text{ of replicates} / \sqrt{n})$$

Where “n” is the number of replicates per sample and the standard deviation (σ) of sample replicates is provided from QuantStudio Real time qPCR software.

9.3 Additional qPCR amplicons as appropriate may be used to clarify and further characterize the test sample as required.

9.4 Any additional analytical methods, such as melting curve analysis using SYBR Green dyes, TD-qPCR, etc., are performed according to this SOP or as directed by the Project Scientist, or QC requestor.

10. DOCUMENTATION AND RECORDS

10.1 Document the preparation of control/standard nucleic acids in the **Preparation of Control or Standard DNA/RNA Logbook (Form 22195-06)**.

10.2 Document the resuspension of primers and probes in the **Primer Resuspension Logbook (Form 22195-07)**.

10.3 Document the use of the PCR Workstations for Master Mix preparation and sample set-up in the **PCR Workstation Logbook (Form 21531-05)**.

10.4 All documentation must be performed in accordance with **SOP 21409 – Good Documentation Practices**.

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

10.5 Record all reagent and material suppliers (indicate BDP, otherwise indicate the vendor), part numbers, lot numbers, and expiration dates on **Form 22195-01**. Record the 96-well plate layout, Master Mix preparation volumes and the volumes/concentrations of Standard used to construct the standard curve on **Form 22195-02**. Record assay parameters and control sample data on **Form 22195-03**. Record test sample data and post-instrument analysis results on **Form 22195-05**. Record Dixon-Q standard curve analysis results (if required) on Form 22195-04. Record DNA extraction procedure on **Form 22972-01, 22972-02**. Record viral RNA extraction procedure on **Form 23113-01, 23113-03, 23113-04**. Refer to any amplicon, product, procedural, and instrument SOPs as necessary.

11. REFERENCES AND RELATED DOCUMENTS

Document Number	Title
21409	Good Documentation Practices
22195-01	Quantitative PCR Reagents
22195-02	96-Well qPCR Plate Key
22195-03	Quantitative PCR Analysis Report
22195-04	Two-Sided Dixon-Q Standard Outlier Test for qPCR Standard Curve Data
22195-05	Quantitative PCR Data
22195-06	Preparation of Control or Standard DNA/RNA
22195-07	Primer Working Stock Logbook
22972	MagNA Pure DNA Extraction
22972-01	Global Run Settings
22972-02	Magna Pure 24 Reagent and Sample Prep
22973	Operation and Maintenance of the QuantStudio™ 6 qPCR System
23113	QIAamp DNA RNA Extraction

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Document Number	Title
23113-01	Sample Preparation Equipment, Reagents, Materials
23113-03	Virus RNA Extraction
23113-04	DNase treatment of Viral RNA
26101	Labeling, Transport, Submission, Storage, and Handling of Biohazardous Materials within the BDP

12. ATTACHMENTS

Attachment 1-64 Assay-Specific Appendices

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 1

AMPLICON: *E. coli* genomic DNA

AMPLICON SEQUENCE:

CGGGCTTTCTTTTGCTATCGACCAGGCGTAGCGGCATATATTCCGGTTGCTCATGATTATCTTGATGCGC
 CCGATC

AMPLICON REGION: *E. coli* Polymerase I (polA) gene

AMPLICON SIZE: 77 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: EC1107F 5' CGGGCTTTCTTTTGCTATCG 3'

REVERSE PRIMER: EC1183R 5' GATCGGGCGCATCAAGATAA 3'

PROBE: EC1134P 5' [6-FAM] CGTAGCGGCATATATTCCGGTTGCTCAT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

EC1107F Forward Primer / EC1183R Reverse Primer BDP PN 30765

EC1134P Probe BDP PN 30761

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: *E. coli* total genomic DNA (ATCC 10798D)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (5 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	≥40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 2

AMPLICON: Human genomic DNA

AMPLICON SEQUENCE:

TGGTAGAGACAAAGGAGACACATTTTATCCATGGACCCAAAACCTCGGCGCCGGTCACGGACTGGGAAG

AMPLICON REGION: Human Endogenous Retrovirus HERV-H2 LTR
AMPLICON SIZE: 68 bp
COPY NUMBER/GENOME: ~ 862
FORWARD PRIMER: HERV 712F 5' TGGTAGAGACAAAGGAGACACATTTT 3'
REVERSE PRIMER: HERV 712R 5' CTTCCCAGTCCGTGACCG 3'
PROBE: HERV712P 5' [6~FAM] CCATGGACCCAAAACCTCGGCG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

HERV712F Forward Primer / HERV712R Reverse Primer BDP PN 30764

HERV712P Probe BDP PN 30763

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Total human genomic DNA (BDP PN 31222)

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 3

AMPLICON: Murine genomic DNA

AMPLICON SEQUENCE:

TACAGAAGTGGATGCTCACAGTCANCTATTGGATGGGTCACACAGGGCCCCCAATGGAGGAGCTAGAGA
 AAGTGCCCAAGGA

AMPLICON REGION: SINE-R Murine retro-transposonic repeat sequence
AMPLICON SIZE: 81 bp
COPY NUMBER/GENOME: ~ 0.6% of murine genome
FORWARD PRIMER: SHF 5' TACAGAAGTGGATGCTCACAGTCA 3'
REVERSE PRIMER: SHR 5' TCCTTGGGCACTTTCTCTAGCT 3'
PROBE: NOVAP 5' [6~FAM] CACAGGGCCCCCAATGGAGGA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

SHF Forward Primer / SHR Reverse Primer BDP PN 30766
NOVAP Probe BDP PN 30759
RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Total murine genomic DNA (Clontech 6650-1)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 4

AMPLICON: Tamarin (*saguinus oedipus*) genomic DNA

AMPLICON SEQUENCE:

TCCAGCAGGAAAAACCTTGGTCCCAAGTACCCCTAGGCTCATATGGCTACCTCATTCAAGAAAAGGGCC
TTATG

AMPLICON REGION: mspE intronic sequence specific for Tamarin
AMPLICON SIZE: 74 bp
COPY NUMBER/GENOME: 3
FORWARD PRIMER: TAMMSP1FWD 5'TCCAGCAGGAAAAACCTTGGT 3'
REVERSE PRIMER: TAMMSP1REV 5' CATAAGGCCCTTTTCTTGAATGA 3'
PROBE: TAMMSP1PROBE 5' [6~FAM] CCAAGTACCCCTAGGCTCATATGGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

Tamarin Forward and Reverse PCR primers BDP PN 30906
Tamarin PCR Probe BDP PN 30907

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Total Tamarin genomic DNA (liver) (BDPQC)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 5

AMPLICON: Human genomic DNA (β -Actin)

AMPLICON SEQUENCE:

tcaccacactgtgccatctacgaggggtatgccctccccatgccatcctgctgctggacctggctggccgggacctgactgactacctcatgaagatc
ctcaccgagcgggctacagcttcaccaccacggccgagcgggaaatcgctgctgacattaaggagaagctgtgctacgtcgccctggactcgagca
agagatggccacggctgctccagctcctcctggagaagagctacgagctgctgacggccaggctatcaccattggcaatgagcgggtccgctg

AMPLICON REGION: Human Beta-Actin gene (ABI commercial kit)
AMPLICON SIZE: 294 bp
COPY NUMBER/GENOME: 6
FORWARD PRIMER: fwd 5' TCACCCACACTGTGCCATCTACGA 3'
REVERSE PRIMER: rev 5' CAGCGGAACCGCTCATTGCCAATTG 3'
PROBE: probe 5' [6~FAM] ATGCC(T)CCCCATGCCATCCTGCG [TAMRA]
3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Total human genomic DNA (Invitrogen 14410-012)

REAGENT	VOLUME/REACTION (μ L)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	N/A
Primers	5 each (15 pmol each)
Probe	5 (10 pmol)
Reaction Mix Volume	40
Sample Volume / rxn	10
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	\geq 40

RESULTS ANALYSIS:

Threshold: 0.2 Δ Rn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 6

AMPLICON: MuLV RNA (RT-qPCR)

AMPLICON SEQUENCE:

AAAGGTGTCCTAACGCAAAAACCTGGGACCTTGGCGTCGGCCGGTGGCCTACCTGTCCAAAAAGCTAGAC
 CCAGTAGCA

AMPLICON REGION: Pol gene
AMPLICON SIZE: 77 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: MU3465F 5' AAAGGTGTCCTAACGCAAAAACCTG 3'
REVERSE PRIMER: MU3542R 5' TGCTACTGGGTCTAGCTTTTTTGG 3'
PROBE: MU3498P 5' [6~FAM] CGTCGGCCGGTGGCCTACCT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP 30450)

CONTROL GENOMIC RNA: X-MuLV pNFS Th-1 viral standard (BDP 30378)

Reagent	Volume/reaction (µl)
One step rtqpcr kit (30265) 2x mix	25 (1x)
Water	15.4
20x primer/probe set (or each primers and probe)	2.5 (1, 1, 0.5)
ROX reference dye	0.1
RNAse out (31244)	1
Superscript/taq mix	1
Reaction mix volume	45
Sample volume / rxn	5
Total reaction volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 7

AMPLICON: Adenovirus E1A region “Mark II” amplicon

AMPLICON SEQUENCE:

GAACCACCTACCCTTCACGAACTGTATGATTTAGACGTGACGGCCCCCGAAGATCCCAACGAGGAGGC

AMPLICON REGION: Ad5 E1A locus (present in HEK293FT and 293Vec RD114)

AMPLICON SIZE: 68 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: E1A158F 5' GAACCACCTACCCTTCACGAACT 3'

REVERSE PRIMER: E1A225R 5' GCCTCCTCGTTGGGATCTTC 3'

PROBE: E1A182P 5' [6~FAM] TATGATTTAGACGTGACGGCCCC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

E1A158F Forward Primer / E1A225R Reverse Primer BDP PN 30769

E1A182P Probe BDP PN 30762

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: E1A control DNA (BDP 31456)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	11
Primers	1 (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	37.5
Sample Volume / rxn	12.5
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn or Automatic

Positive DNA control (E1A control, 107 bp, BDP 31456)

5'- CCTCCTAGCCATTTTGAACCACCTACCCTTCACGAACTGTATGATTTAGACGTGACGGCCCC
 CGAAGATCCCAACGAGGAGGCGGTTTCGCAGATTTTTCCCGACTC -3'

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 8

AMPLICON: Adenovirus E1B region

AMPLICON SEQUENCE:

AGGTTTACTGGCCCCAATTTTAGCGGTACGGTTTTCTGGCCAATACCAACCTTATCCTACACGGTGTAAGCTTCTATGG

AMPLICON REGION: Ad5 E1B locus
AMPLICON SIZE: 80 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: E1B2736F 5' AGGTTTACTGGCCCCAATTTTAG 3'
REVERSE PRIMER: E1B2791R 5'CCATAGAAGCTTACACCGTGTAGGA3'
PROBE: E1B2764P 5' [6~FAM] CGGTTTTCTGGCCAAYACCAACCTT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Adenoviral reference material DNA

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 9

AMPLICON: Adenovirus Hexon region

AMPLICON SEQUENCE:

CCTACTCTGGCACTGCCTACAACGCCCTGGCTCCCAAGGGTGCCCCAAATCCTTGCGAATGGGATG

AMPLICON REGION: Ad5 Hexon gene
AMPLICON SIZE: 65 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: HEX568F 5' CCTACTCTGGCACTGCCTACAA 3'
REVERSE PRIMER: HEX633R 5' CATCCCATTCGCAAGGATTT 3'
PROBE: HEX593P 5' [6~FAM] CCTGGCTCCCAAGGGTGCCC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

Hex568F Forward Primer / Hex633R Reverse Primer BDP PN 30767

Hex593P Probe BDP PN 30757

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Ad5 Hexon DNA control BDP PN 31420

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

Positive control (Ad5 Hexon DNA, 110bp, BDP 31420)

5'- GGGGCCCTACTTTTAAGCCCTACTCTGGCACTGCCTACAACGCCCTGGCTCCCAAGGGTGC
 CCCAAATCCTTGCGAATGGGATGAAGCTGCTACTGCTCTTGAATAAAC -3'

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 10

AMPLICON: HEK293 Adenovirus E1a Left Hand Junction

AMPLICON SEQUENCE:

CACCATAACAACAGCTACTCAAAGTGTAAACCAGGATAACAAGTTGATGACTTGCCATCATCAATAATATACCTTATTTTGGATTGAAGCCAA

AMPLICON REGION: Left Hand Junction of HEK293-Ad5 E1 genomic DNA
AMPLICON SIZE: 109 bp
13. COPY NUMBER/GENOME: 1
FORWARD PRIMER: LHF 5' CACCATAACAACAGCTACTCAAAGTGT
REVERSE PRIMER: LHR 5' TTGGCTTCAATCCAAAATAAGGTAT 3'
PROBE: LHP 5' [6~FAM] AACCAGGATAACAAGTTGATGACTTGCCATCA 3' [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: HEK293 genomic DNA (ATCC CRL-1573)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 11

AMPLICON: HEK293 Adenovirus E1a Right Hand Junction

AMPLICON SEQUENCE:

TGGTGTGTTGTAGATGATCCAGTCGTAGCAGGAGCGCTGGGCGTGGTGCCTAAAAATGTCTTTCAGTAGCAAGCTGATTGCCAAAGATGTCAGAACAAGACTCCCCATCATGATAAGGCTCCCACCCCTCTTAACTGTCCTGCTCATGCCTG

AMPLICON REGION: Right Hand Junction of HEK293-Ad5 E1 genomic DNA
AMPLICON SIZE: 151 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: RHF 5' TGGTGTGTTG TAGATGATCCAGTCGTA 3'
REVERSE PRIMER: RHR 5' CAGGCATGAGCAAGGACAGTT 3'
PROBE: RHP 5' [6~FAM] TGATTGCCAAAGATGTCAGAACAAGACTCC 3' [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: HEK293 genomic DNA (ATCC CRL-1573)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	:15	40
	60°C	1:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 12

AMPLICON: Adenovirus E1a Δ 24 region (for detection of wild-type Ad5 E1a)

AMPLICON SEQUENCE:

GTCCGGTTTCTATGCCAAACCTTGTACCGGAGGTGATCGATCTTACCTGCCAGGAGGCTGGCTTTCCAC
 CCAGTGACGACGA

AMPLICON REGION: Ad5 E1a Δ 24 region for detection of wild-type E1a sequence Refer to SOP 22172

AMPLICON SIZE: 79 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: RGDF 5' CGGTTTCTATGCCAAACCTTGT 3'

REVERSE PRIMER: RGDR 5' TCGTCGTCACTGGGTGGAA 3'

PROBE: RDGP 5' [6~FAM] ATCGATCTTACCTGCCAGGAGGCTG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Adenoviral reference material DNA

REAGENT	VOLUME/REACTION (μ L)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	0.5
Primers	2 each (20 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 Δ Rn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 13

AMPLICON: Lambda phage DNA

AMPLICON SEQUENCE:

GGTAAACATGGCGCTGTACGTTTCGCCGATTGTTTCCGGTGAGGTTATCCGTTCCCGTGGCG

AMPLICON REGION: λ capsid E gene

AMPLICON SIZE: 60 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: LAMF 5' GGTAAACATGGCGCTGTACGT 3'

REVERSE PRIMER: LAMR 5' CGCCACGGGAACGGATA 3'

PROBE: LAMP 5' [6~FAM] TCGCCGATTGTTTCCGGTGAGG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: lambda phage DNA [clind1ts857 Sam7]

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 14

AMPLICON: Polio wild type (Sabin-1 strain) virus

AMPLICON SEQUENCE:

CGTGGCGGAACCGACTACTTTGGGGTGCCGTGTTTCCTTTATTTTATTGTGGCTGCTTATGGTGACA

AMPLICON REGION: 5' – UTR/IRES region Sabin-1 strain polio
AMPLICON SIZE: 68 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: PVS2F 5' CGTGGCGGAACCGACTAC 3'
REVERSE PRIMER: PVS2R 5' TGTCACCATAAGCAGCCACAA 3'
PROBE: PVS1P 5' [6~FAM] TGTCCTGTTTCCTTTATTTT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP 30450)
CONTROL GENOMIC RNA: Sabin strain 1 wild-type polio virus

Reagent	Volume/reaction (µl)
One step RTqPCR kit (30265) 2x mix	25 (1x)
Water	15.4
20x primer/probe set (or each primers and probe)	2.5 (1, 1, 0.5)
ROX reference dye	0.1
RNase out (31244)	1
Superscript/taq mix	1
Reaction mix volume	45
Sample volume / rxn	5
Total reaction volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 15

AMPLICON: PVS-RIPO Plasmid and virus (HRV-2 IRES)

AMPLICON SEQUENCE:

AACCCAATGTGTATCTAGTCGTAATGAGCAATTGCGGGATGGGACCAACTACTTTGGGTGTCCGTGTTCA

AMPLICON REGION: PVS-RIPO HRV-2 IRES sequence
AMPLICON SIZE: 71 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: PVS1F 5' AACCCAATGTGTATCTAGTCGTAATGA 3'
REVERSE PRIMER: PVS1R 5' TGAAACACGGACACCCAAAG 3'
PROBE: PVS2P 5' [6~FAM] CAATTGCGGGATGGGACCAACT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP 30450)

CONTROL GENOMIC DNA: PVS-RIPO plasmid (BDP L0305006)

Reagent	Volume/reaction (µl)
One step RTqPCR kit (30265) 2x mix	25 (1x)
Water	15.4
20x primer/probe set (or each primers and probe)	2.5 (1, 1, 0.5)
ROX reference dye	0.1
RNase out (31244)	1
Superscript/taq mix	1
Reaction mix volume	45
Sample volume / rxn	5
Total reaction volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 ΔR_n

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 16

AMPLICON: Wolinella Asparaginase AsnA

AMPLICON SEQUENCE:

CACCCCGATGATACTGATGTTTTAGTCAATGCAGCCCTTCAGGCAGGAGCCAAAGGAATCATCCATGCA
G

AMPLICON REGION: Wolinella Asparaginase gene [AsnA] in plasmid rWS-Asp pET28a

AMPLICON SIZE: 70 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: WOLF 5' CAC CCC GAT GAT ACT GAT GTT TT 3'

REVERSE PRIMER: WOLR 5' CTG CAT GGA TGA TTC CTT TGG 3'

PROBE: WOLPR 5' [6~FAM] TCA ATG CAG CCC TTC AGG CAG GA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: rWS-Asp/pET28a plasmid L0401018

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	:15	40
	60°C	1:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 17

AMPLICON: *E. coli* Asparaginase gene AsnB

AMPLICON SEQUENCE:

TACCCCGGC ACGTAAGCAT ACCAGCGACA CGCCATTCGA TGTCTCTAAG CTGAATGAAC TGCCGA

AMPLICON REGION: *E. coli* asparaginase AsnB gene (*E. coli* gene most homologous to Wolinella AsnA)

AMPLICON SIZE: 65 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: ECBF 5' TAC CCC GGC ACG TAA GCA T 3'

REVERSE PRIMER: ECBR 5' - TCG GCA GTT CAT TCA GCT TAG A 3'

PROBE: ECBPR 5' [6~FAM] CCA GCG ACA CGC CAT TCG ATG T [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: *E. coli* total genomic DNA (ATCC 10798D)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	:15	40
	60°C	1:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 18

AMPLICON: *E. coli* asparaginase gene AsnA

AMPLICON SEQUENCE:

TGACGTCGTGCGCAATTTTCTGCGCCAACCGGTGAAAGCATTGATTCTGCGCTCCTATGGC

AMPLICON REGION: *E. coli* asparaginase gene AsnA (low homology to Wol AsnA)
AMPLICON SIZE: 61 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: ECAF 5' TGA CGT CGT GCG CAA TTT T 3'
REVERSE PRIMER: ECAR 5' GCC ATA GGA GCG CAG AAT CA 3'
PROBE: ECAPR 5' [6~FAM] TGC GCC AAC CGG TGA AAG CA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: *E. coli* total genomic DNA (ATCC 10798D)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 19

AMPLICON: *E. coli* putative asparaginase [ybiK] gene

AMPLICON SEQUENCE:

TGGATTACGGCGGATTAAGTCTCGCGGAAGCCTGCGAGCGGGTAGTAATGGAAAACTCCCTGCG

AMPLICON REGION:

E. coli putative asparaginase [ybiK] gene
(low homology to Wol AsnA)

AMPLICON SIZE:

65 bp

COPY NUMBER/GENOME:

1

FORWARD PRIMER: YBIF

5' TGG ATT ACG GCG GAT TAA GTC T 3'

REVERSE PRIMER: YBIR

5' CGC AGG GAG TTT TTC CAT TAC T 3'

PROBE: YBIPR

5' [6~FAM] CGG AAG CCT GCG AGC GGG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT:

QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA:

E. coli total genomic DNA (ATCC 10798D)

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	:15	40
	60°C	1:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 20

AMPLICON: Herpes Simplex 1 (HSV-1) Major Capsid gene

AMPLICON SEQUENCE:

AACAGC ACGTTCTCAG TCACAAAGCG GTCCTGTCGG ACGACGGTGA ACCCAAACCC GGGATGGAGG
 CCCGTCTTGA GCTGATGATG CAAGG

AMPLICON REGION: Universal HSV-1 major capsid gene UL36
AMPLICON SIZE: 91 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: HSVCAPF 5' AACAGCACGTTCTCAGTCACAAA 3'
REVERSE PRIMER: HSVCAPR 5' CCTTGC ATCATCAGCTCAAGAC 3'
PROBE: HSVCAPPR 5' [6~FAM] ACCCAAACCCGGGATGGAGGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

HSV CAP forward and reverse primers BDP PN 30944
HSV CAP probe BDP PN 30945
RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: rHSV-1 DNA reference material ARR (BDP 02012004C)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 21

AMPLICON: HSV-1 Ribonucleotide reductase gene (UL39 ICP6)

AMPLICON SEQUENCE:

GCCTGTGTCTGGACGTTCCCTCCGGTCCCGCCGAACGCATACATGCCCTATTATCTCAGGGAGTATGTGACG

AMPLICON REGION: HSV-1 RR large subunit gene UL39
 Amplicon is not present in rRp450-HSV-1

AMPLICON SIZE: 71 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: HSVRRF 5' GCCTGTGTCTGGACGTTCCCT 3'

REVERSE PRIMER: HSVRRR 5' CGTCACATACTCCCTGAGATAATAGG 3'

PROBE: HSVRRP 5' [6~FAM] CGGTCCCGCCGAACGCAT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: wild-type HSV-1 DNA extract (ABi 08-705-000)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 22

AMPLICON: Acholeplasma 16S Ribosomal DNA

AMPLICON SEQUENCE:

AACAAA GGGCACACAG TGGATGCCTT GGCACCAAGA GGCGATGAAG GACGGAACTA ACACCGAAAT
 GCTCGGG

AMPLICON REGION: Mollicute *Acholeplasma laidlawii* 16S rDNA
AMPLICON SIZE: 73 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: ARF 5' AAC AAA GGG CAC ACA GTG GAT 3'
REVERSE PRIMER: ARR 5' CCC GAG CAT TTC GGT GTT AG 3'
PROBE: ARPR 5' [6~FAM] CAC CAA GAG GCG ATG AAG GAC GGA
 [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: *A. laidlawii* gDNA (ATCC 23206D)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: Variable ΔR_n

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 23

AMPLICON: Tn903 Kanamycin resistance gene [AphA]

AMPLICON SEQUENCE:

GATGTTGGACGAGTCGGAATCGCAGACCGATACCAGGATCTTGCCATCCTATGGAAGTGCCTCGGTG

AMPLICON REGION: Tn903 derived AphA kanamycin resistance gene Present in most KanR plasmids

AMPLICON SIZE: 67 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: KAN689F 5' GATGTTGGACGAGTCGGAATC 3'

REVERSE PRIMER: KAN755R 5' CACCGAGGCAGTTCATAGG 3'

PROBE: KAN711P 5' [6~FAM or VIC] CAGACCGATACCAGGATCTTGCCA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

Kanamycin TN903 dPCR Primers/probes (FAM) set BDP PN 31227

Kanamycin TN903 dPCR Primer/probe set (VIC/HEX) BDP PN 31232

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: plasmid pNGVL4a-CRT/E7(detox) L0506020

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 24

AMPLICON: Chloramphenicol resistance gene [CAT]

AMPLICON SEQUENCE:

AACGGCATGATGAACCTGAATCGCCAGCGGCATCAGCACCTTGTGCGCCTTGCGTATAATATTTGCCCAT
GGTGAA

AMPLICON REGION: [CAT] chloramphenicol acetyltransferase gene (CamR)
AMPLICON SIZE: 75 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: CAMP147F 5' AACGGCATGATGAACCTGAATC 3'
REVERSE PRIMER: CAMP221R 5' TTCACCATGGGCAAATATTATACG 3'
PROBE: CAMP174P 5' [6~FAM] CGGCATCAGCACCTTGTGCGCC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: HA22 (plasmid pRB698 VL) lot L0212018

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 25

AMPLICON: Tn5 Neomycin/Kanamycin resistance gene

AMPLICON SEQUENCE:

TGGATTGCACGCAGGTTCTCCGCGCCGCTTGGGTGGAGAGGCTATTCGGCTATGACTGGGCACA

AMPLICON REGION: Neomycin / Kanamycin resistance gene from Tn5
AMPLICON SIZE: 62 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: NEOF 5' TGGATTGCACGCAGGTTCT 3'
REVERSE PRIMER: NEOR 5' GTGCCAGTCATAGCCGAAT 3'
PROBE: NEOPR 5' [6~FAM] CGGCCGCTTGGGTGGAGAGG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

NEOF Primer BDP PN 31154
NEOR Primer BDP PN 31155
NEOPR Probe BDP PN 31156

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: plasmid pHGM-CSF L0312004

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 26

AMPLICON: VERO (and CV-1) cell line genomic DNA

AMPLICON SEQUENCE:

cctctgcc cagcgtgaag atcaccagg tcacgtggca gaagatcacc caagccacca atggctccaa gcagaacgtg gccatctaca
 accatccat gggcgtgtct gtg

AMPLICON REGION: *Cercopithecus aethiops* Nectin-1 gene
AMPLICON SIZE: 111 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: NECF-1 5' CCTCTGCCCAGCGTGAAG 3'
REVERSE PRIMER: NECR-1 5' CACAGACACGCCCATGGAT 3'
PROBE: NECPR-1 5' [6~FAM] CACCCAAGCCACCAATGGCTCCAA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

NECF Forward Primer / NECR Reverse Primer BDP PN 30770
NECPr Probe BDP PN 30760

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: CV-1 genomic DNA (WCB BDP 10205, EOP 10377)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	:15	40
	60°C	1:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 27

AMPLICON: Universal Polio virus – polyprotein 3C gene

AMPLICON SEQUENCE:

TTGGTGGGAACGGTTCACACGGGTTTGCAGCGGCCCTGAAGCGATCATACTTCACTCAGAGTCAAGGTG
A

AMPLICON REGION: Polio polyprotein 3C region in Sabin-1 strain and PVR-RIPO virus

AMPLICON SIZE: 70 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: PO1F 5' TTG GTG GGA ACG GTT CAC A 3'

REVERSE PRIMER: PO1R 5' TCA CCT TGA CTC TGA GTG AAG TAT GA 3'

PROBE: PO1PR 5' [6~FAM] TTG CAG CGG CCC TGA AGC G [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

Polio PPP primers PO1F and PO1R BDP PN 30947

Polio PPP PO1 Probe BDP PN 30948

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP 30450)

CONTROL GENOMIC RNA: Sabin strain 1 wild-type polio virus WHO std.

Reagent	Volume/reaction (µl)
One step RTqPCR kit (30265) 2x mix	25 (1x)
Water	15.4
20x primer/probe set (or each primers and probe)	2.5 (1, 1, 0.5)
ROX reference dye	0.1
RNAse out (31244)	1
Superscript/taq mix	1
Reaction mix volume	45
Sample volume / rxn	5
Total reaction volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 28

AMPLICON: Polio Type 1 and Type 2 specific IRES (Sabin and Mahoney)

AMPLICON SEQUENCE:

TTGGCGG CCTACCTATG GCTAACGCCA TGGGACGCTA GTTGTGAACA AGGTGTGAAGAGCCTATTGA
 GCTACATAAG AATCCTCCGG CCCCTGAATG CGGCTAATCC CA

AMPLICON REGION: Polio type 1 and type 2 IRES, will not amplify PVS-RIPO

AMPLICON SIZE: 109 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: POSAF 5' TTG GCG GCC TAC CTA TGG 3'

REVERSE PRIMER: POSAR 5' TGG GAT TAG CCG CAT TCA G 3'

PROBE: POSAPR 5' [6~FAM] AGC CTA TTG AGC TAC ATA AGA ATC CTC CGG C [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP 30450)

CONTROL GENOMIC RNA: Sabin strain 1 wild-type polio virus WHO std.

REAGENT	VOLUME/REACTION (ML)
ONE STEP RTQPCR KIT (30265) 2X MIX	25 (1X)
WATER	15.4
20X PRIMER/PROBE SET (OR EACH PRIMERS AND PROBE)	2.5 (1, 1, 0.5)
ROX REFERENCE DYE	0.1
RNASE OUT (31244)	1
SUPERSCRIPT/TAQ MIX	1
REACTION MIX VOLUME	45
SAMPLE VOLUME / RXN	5
TOTAL REACTION VOLUME	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 29

AMPLICON: CHO cell line genomic DNA - DHFR gene

AMPLICON SEQUENCE:

CCGGTGAA GTGAGTGGAG AAAGGGGATA CGAAGACAGC ATCCCACATG ACTGCTCCCA
 GTAAAGGCAA GGTCTTCATC CAT

AMPLICON REGION: Chinese Hamster Ovary cell line K1 DHFR (MetR)

AMPLICON SIZE: 81 bp

COPY NUMBER/GENOME: 1 (for wild-type Chinese Hamster and non-K1 cell lines),
 > 1000 copies (for K1 derived cell lines such as C400)

FORWARD PRIMER: CHO5243F 5' CCGGTGAAGT GAGTGGAGAAA 3'

REVERSE PRIMER: CHO5300R 5' ATGGATGAAGACCTTGCCTTTACT 3'

PROBE: CHO5264P 5' [6~FAM] TACGAAGACAGCATCCCACATGACTGCT [Quencher] 3

REAGENTS (and associated BDP and/or manufacturer part numbers):

CHO5243F Forward Primer BDP PN 31276

CHO5300R Reverse Primer BDP PN 31277

CHO5264P Probe BDP PN 31278

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: CHO-CC49(dhfr-) – one copy/haploid genome (BDP 012203)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 30

AMPLICON: CHO cell line genomic DNA–Endogenous type-C retrovirus sequence

AMPLICON SEQUENCE:

CCCG GACCCCTGAG TCACCGGACT GCATGGGTAC AGGGAGCTAC AGGCGGAAAG CAGTACCATT
 GGACTACAAA TCGGCAGCTC CAGCTCGCGA CTGGT

AMPLICON REGION: CHO cell type-C proviral (endogenous retroviral) pol-region

AMPLICON SIZE: 98 bp

COPY NUMBER/GENOME: ~ 100 – 300 copies

FORWARD PRIMER: CHF 5' CCC GGA CCC CTG AGT CAC 3'

REVERSE PRIMER: CHR 5' ACC AGT CGC GAG CTG GAG 3'

PROBE: CHPR 5' [6~FAM] TGG GTA CAG GGA GCT ACA GGC GGA AA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

CHF Forward Primer / CHR Reverse Primer BDP PN 30768

CHPr Probe BDP PN 30758

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236) For proviral RNA use Viral Mini Kit (30450)

CONTROL GENOMIC DNA: CHO-CC49(dhfr-) – one copy/haploid genome (BDP 012203)

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 31

AMPLICON: PER.C6 cell line genomic DNA HDEP Junction (Ad5 E1b-βlac)

AMPLICON SEQUENCE:

CTGGCCAATACCAACCTTATCCTACACGGTGTAAGCTTCTAAGGCCGTATCGTAGTTATCTACACGACGG
GGAGTCAGG

AMPLICON REGION: Ad5 E1b – plasmid β-lac junction sequence from
AMPLICON SIZE: 79 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: HDEPF 5' CTGGCCAATACCAACCTTATCC 3'
REVERSE PRIMER: HDEPR 5' CCTGACTCCCCGTCTGTAG 3'
PROBE: HDEPPR 5' [6~FAM] ACACGGTGTAAGCTTCTAAGGCCGTATCGTAGT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: plasmid pTOPO211-213 (CruCell E273-087, QC026038)

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 32

AMPLICON: Adenovirus 5 E1a (Mark I) amplicon for KD3 E1a RCA testing

AMPLICON SEQUENCE:

GACACCG GGACTGAAAA TGAGACATAT TATCTGCCAC GGAGGTGTTA TTACCGAAGA AATGGCCGC

AMPLICON REGION: Ad5 E1a region deleted in virus Ad-KD3, present in wtAd5

AMPLICON SIZE: 66 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: E1A20F 5' GAC ACC GGG ACT GAA AAT GAG 3'

REVERSE PRIMER: E1A93R 5' GCG GCC ATT TCT TCG GTA A 3'

PROBE: E1A50P 5' [6~FAM] ACA TAT TAT CTG CCA CGG AGG TGT T [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Adenoviral reference material DNA

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 33

AMPLICON: Adenovirus 5 wild-type fiber tropism (non-RGD)

AMPLICON SEQUENCE:

TTACTACTAAACGGGTACACAGGAAACAGGAGACACAAC TCCAAGTGCATACTCTATGTCATTTTCATGGGA CTGGTCTGGC

AMPLICON REGION: Wild-type Ad5 wild-type fiber knob sequence (for non-RGD virus)
AMPLICON SIZE: 80 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: WTFIBF 5' TTACTACTAAACGGGTACACAGGAAACAG 3'
REVERSE PRIMER: WTFIBR 5' GCC AGA CCA GTC CCA TGA AA 3'
PROBE: WTFIBPR 5' [6~FAM] CAC AAC TCC AAG TGC ATA CTC TAT GTC A [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Adenoviral reference material DNA

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 34

AMPLICON: Adenovirus 5 RGD modified fiber knob tropism

AMPLICON SEQUENCE:

CACACTAAACGGTACACAGGAAACAGGAGACACAACCTTGTGACTGCCGCGGAGACTGTTTCTGCCCA7C
TGCATACTCTATGTCATTTTCATGGGACTGGTCTGGC

AMPLICON REGION: RGD tropism modified Ad5 fiber knob sequence
 (does not amplify wild-type fiber sequence)

AMPLICON SIZE: 105 bp

COPY NUMBER/GENOME: 1

FORWARD PRIMER: LAMF 5' ACA CTA AAC GGT ACA CAG GAA ACA 3'

REVERSE PRIMER: LAMR 5' GCC AGA CCA GTC CCA TGA AA 3'

PROBE: LAMP 5' [6~FAM] CCG CGG AGA CTG TTT CTG CCC A [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: Ad5-D24-RGD virus L0405025

REAGENT	VOLUME/REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.5
Primers	1 each (10 pmol each)
Probe	0.5 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C 60°C	:15 1:00	40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 35

AMPLICON: AAV2

AMPLICON SEQUENCE:

CACAAGTTTGGGATCTAATACAATGGCTTCAGGCGGTGGCGCACCAATGGCAGACAATAACGAGGGTGC
CG

AMPLICON REGION: Rep 78 protein, protein VP1 (cap) genes
AMPLICON SIZE: 71 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: AAV F 5' ACAAGTTTGGGATCTAATACAATGG 3'
REVERSE PRIMER: AAV R 5' CGGCACCCTCGTTATTGTCT 3'
PROBE: AAV P 5' [6-FAM] TTCAGGCGGTGGCGCACCA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid Rep 78 protein / protein VP1 (capsid) genes, accession #AY695376.1

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 36

AMPLICON: B19

AMPLICON SEQUENCE:

GGCCTGGCAATGAGCTACAAGCTGGGCCCCCGCAAAGTGCTGTTGACAGTGCTGCAAGGATTCATGA

AMPLICON REGION: VP1 gene
AMPLICON SIZE: 67 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: B19 F 5' GGCCTGGCAATGAGCTACAA 3'
REVERSE PRIMER: B19 R 5' TCATGAATCCTTGCAGCACTGT 3'
PROBE: B19 P 5' [6-FAM] CTGGGCCCCCGCAAAGTGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid Capsid protein (VP), accession #Z68146

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 37

AMPLICON: BK/JC

AMPLICON SEQUENCE:

TCCTCAATGGATGTTGCCTTTACTTTTAGGGTTGTACGGGACTGTAACACCTGCTCTTGAAGCATATGAA
 GATGGCCC

AMPLICON REGION: VP2 gene

AMPLICON SIZE: 78 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: BKJC F 5' TCCTCAATGGATGTTGCCTTT 3'

REVERSE PRIMER: BKJC R 5' GGGCCATCTTCATATGCTTCAA 3'

PROBE: BKJC P 5' [6-FAM] TTTTAGGGTTGTACGGGACTGTAACACCTGCT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid Capsid protein VP2 region, spans the sequence of JCV Isolate 330A, #AY382188, and BKV isolate TW-2 #AB213487 VP2 region

REAGENT	VOLUME / REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 38

AMPLICON: CMV

AMPLICON SEQUENCE:

TCCGCGGTTGTCTCTGTGTATAACGTTTTTATTTCCGGTTCCGCGTTTGGTCGCCTGCCT

AMPLICON REGION: UL38 gene
AMPLICON SIZE: 60 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: CMV F 5' TCCGCGGTTGTCTCTGTGTA 3'
REVERSE PRIMER: CMV R 5' AGGCAGGCGACCAAACG 3'
PROBE: CMV P 5' [6-FAM] CGTTTTTATTTCCGGTTCCG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid UL38 gene (pp65 lower matrix protein) from ~120,656 to 122,341 in #AY446894, Merlin strain

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 39

AMPLICON: EBV

AMPLICON SEQUENCE:

CGGCCGTGATGGAGGCTATGACGGCGGCGAGTGACTACGCGCGTGGCCTGGGCGTGAAGCTGACCTT
 TGGCTCGGCCTCCTGCCCCGAGACCGGCTCGTCCGCCTCCAACCTTCATGACCGTGGTGGCCTCTGTCT

AMPLICON REGION: BNR1 gene
AMPLICON SIZE: 134 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: EBV F 5' CGGCCGTGATGGAGGCTATG 3'
REVERSE PRIMER: EBV R 5' AGACAGAGGCCACCACGG 3'

PROBE: EBV P 5' [6-FAM] TGACCTTTGGCTCGGCCTCCTGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid BNR1 gene for p140, accession #X67777.1

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 40
AMPLICON: HHV 6 A/B

AMPLICON SEQUENCE:

GGCGCGACCTCTGACAGTGCTCACGGAGAGGACAGAGGACGTCTCTAGAAAACATCTCATCCACCTGCT
 TAATATAGTCCGTATCTTGCCTAAAC

AMPLICON REGION: Immediate early gene
AMPLICON SIZE: 95 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: HHV6 F 5' GGCGCGACCTCTGACAGT 3'
REVERSE PRIMER: HHV6 R 5' GTTTACGCAAGATACGGACTATATTAAGC 3'
PROBE: HHV6 P 5' [6-FAM] ACAGAGGACGTCTCTAGAAAACATCTCATCCACC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid U25 / EPLF3 gene, accession #L25528 and NC001664, and AF157706 B-strain Z29, and AB021506 B-strain HST

REAGENT	VOLUME / REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 41

AMPLICON: HHV 7

AMPLICON SEQUENCE:

CCGTAAAGCGTGTGGGATTGTCATCCAATGCTTTAACAAACAAATGGAAGACTCTTCATTCTGCACTTCC
 CATGAAAAAA

AMPLICON REGION: Glycoprotein B gene

AMPLICON SIZE: 79 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: HHV7 F 5' CCGTAAAGCGTGTGGGATTG 3'

REVERSE PRIMER: HHV7 R 5' TTTTTCATGGGAAGTGCAGAA 3'

PROBE: HHV7 P 5' [6-FAM] CATCCAATGCTTTAACAAACAAATGGAAGACTCTTC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid UL39 envelope glycoprotein (g)B, accession #AY192554 and NC_001716

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 42

AMPLICON: HHV 8

AMPLICON SEQUENCE:

CTCGAATCCAACGGATTTGACCCCGTGTTCCCCATGGTCGTGCCGCAGCAACTGGGGCACGCT

AMPLICON REGION: Glycoprotein B gene

AMPLICON SIZE: 63 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: HHV8 F 5' CTCGAATCCAACGGATTTGAC 3'

REVERSE PRIMER: HHV8 R 5' AGCGTCCCCAGTTGCT 3'

PROBE: HHV8 P 5' [6-FAM] CGTGTTCCCCATGGTCGTGCC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid ORF 26, minor capsid protein gene, accession #AY219444

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 43

AMPLICON: HIV 1

AMPLICON SEQUENCE:

AGCAGCCATGCAAATGTAAAAGAGACCATCAATGAGGAAGCTGCAGAATGGGATAGATTGCATCCAGTGCATGC

AMPLICON REGION: gag gene
AMPLICON SIZE: 74 bp
COPY NUMBER / GENOME 1
FORWARD PRIMER: HIV1 F 5' AGCAGCCATGCAAATGTAAA 3'
REVERSE PRIMER: HIV1 R 5' GCATGCACTGGATGCACTCT 3'
PROBE: HIV1 P 5' [6-FAM] CATCAATGAGGAAGCTGCAGAATGGGA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid gag gene, accession #308760

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 44

AMPLICON: HIV 2

AMPLICON SEQUENCE:

ATGGGCGCGAGAACTCCGTCTTGAGAGGGAAAAAGCAGATGAATTAGAAAAATTAGGTTACGGCCC
 GGCGGAAAGAAAAAGTACAAGTTAAACATATTGTGTGGGCAGCGAATGAATTGGACAGATTCGGATTAG
 CAAAGAGCCTGTTGGA

AMPLICON REGION: 3' LTR gene
AMPLICON SIZE: 154 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: H2DG F 5' ATGGGCGCGAGAARCTCCG 3'
REVERSE PRIMER: H2DG R 5' TCCAACAGGCTYTCTGCYAATCC 3'
PROBE: H2DG P 5' [6-FAM] YAGGTTACGGCCRRRCGGRAA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid gag gene, accession #HIV2CAM2

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 45

AMPLICON: MMV

AMPLICON SEQUENCE:

GCCCATGATTTGTGCTTGGTTGGTAAAGAATGGTTACCAATCTACCATGGCAAGCTACTGTGCTAAATGG
 GGCAAAGTTCCT

AMPLICON REGION: NS1gene

AMPLICON SIZE: 82 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: MVM F 5' GCCCATGATTTGTGCTTGGTT 3'

REVERSE PRIMER: MVM R 5' AGGAACTTTGCCCATTTAGC 3'

PROBE: MVM P 5' [6-FAM] AAAGAATGGTTACCAATCTACCATGGCAAGCTACT [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

MVM forward and reverse primers BDP PN 30902

MVM Probe BDP PN 30903

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid Major non-structural protein gene, accession #FJ445512.1 and equivalent to most homologous mouse parvovirus strains' sequences.

REAGENT	VOLUME / REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 46

AMPLICON: SMRV

AMPLICON SEQUENCE:

GGCTACTACTTCCCTTGCAAAGATTAGACAAGGCCCGATGAGTCATACAGTGATTTTGTAAAGCCGCCTC
CAGG

AMPLICON REGION: gag protein gene

AMPLICON SIZE: 74 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: SMRV F 5' GGCTACTACTTCCCTTGCAAAGAT 3'

REVERSE PRIMER: SMRV R 5' CCTGGAGGCGGCTTACAAA 3'

PROBE: SMRV P 5' [6-FAM] ACAAGGCCCGATGAGTCATACAGTGA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid gag protein gene, accession #U23805

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 47

AMPLICON: SRV 1/2

AMPLICON SEQUENCE:

AATGAATGCAAATCCAAAAGTATGATAGTCAAGGAAACCCACTACCACCCCATCAGGGAAACGGACTGAGG GGC

AMPLICON REGION: gag gene
AMPLICON SIZE: 71 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: SRV F 5' AATGAATGCAAATCCAAAAGTATG 3'
REVERSE PRIMER: SRV R 5' GCCCCTCAGTCCGTTTCC 3'
PROBE: SRV P 5' [6-FAM] CAAGGAAACCCACTACCACCCCATCA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid gag gene region, SRV-1 accession #M11841 and SRV 2 accession #M16605.

REAGENT	VOLUME / REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 48

AMPLICON: SV 40

AMPLICON SEQUENCE:

CACAGGCCTATGCTGTGATATCTGGGGCTCCTGCTGCTATAGCTGGATTTGCAGCTTTACTGCAAACCTGT
 GACTGGTGTGAGCGCTGTTGCTCAAGTGGGGTATAGATTTTT

AMPLICON REGION: VP2 gene

AMPLICON SIZE: 112 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: SV40 F 5' CACAGGCCTATGCTGTGATATCTG 3'

REVERSE PRIMER: SV40 R 5' AAAAATCTATACCCCACTTGAGCAA 3'

PROBE: SV40 P 5' [6-FAM] AGCTTTACTGCAAACCTGTGACTGGTGTGAGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmid VP2 gene, accession #AF316141

REAGENT	VOLUME / REACTION (µL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 49

AMPLICON: HTLV 1

AMPLICON SEQUENCE:

GGACAAGCATATATTGTCACAAAGATCATTCCCCCTTCCGCCACCGCACAAAGTCGGCCCAA
 CGGGCCGAAGTTCTCGGACTTTTGCATGG

AMPLICON REGION: gag pol gene
AMPLICON SIZE: 89 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: HTLV1 F 5' CCATGCAAAAGTCCGAGAAGTT 3'
REVERSE PRIMER: HTLV1 R 5' GGACAAGCATATATTGTCACAAAGATC 3'
PROBE: HTLV1 P 5' [6-FAM] CGTTGGGCCGACTTGTGCGG [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmids cloned gag/pol genes for each strain, various accession #'s.

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 50

AMPLICON: HTLV 2

AMPLICON SEQUENCE:

TTGCGACTTCCTGAGGAACTCCCCTCATCCAAGTGTCGGCATCCTCATTACCCACATGGGTTCGATTCCAT
AACCTTGGCAGCCA

AMPLICON REGION: rex protein gene

AMPLICON SIZE: 83 bp

COPY NUMBER / GENOME: 1

FORWARD PRIMER: HTLV2 F 5' TGGCTGCCAAGGTTATGGA 3'

REVERSE PRIMER: HTLV2 R 5' TTGCGACTTCCTGAGGAACTC 3'

PROBE: HTLV2 P 5' [6-FAM] CGACCCATGTGGTGAATGAGGATGC [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

(Note – Viral detection qPCR assays for cell line characterization require target input quantitation, e.g. total DNA mass, genomic copy number, cell count, or similar measure)

CONTROL GENOMIC DNA: Plasmids cloned gag/pol genes for each strain, various accession #'s.

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	25 (1x)
NFW or WFI Water (various)	2.0
Primers	1 each (5 pmol each)
Probe	1 (5 pmol)
Reaction Mix Volume	30
Sample Volume / rxn	20
Total Reaction Volume	50

REACTION CONDITIONS: This method detects proviral DNA when used without an RT step.

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 51

AMPLICON: Human AADC Exon Junction Amplicon from rAAV2-CMV-hAADC

AMPLICON SEQUENCE:

CTGGCCTGATTCTTTCTTTATGGTTGCCACCCTGGGGACCACAACATGCTGCTCCTTTGACAATCTCTT
 AGAAGTCGGTCTATCTGCAACAAGGAAGACATATG

AMPLICON REGION: Human AADC from rAAV2-CMV-hAADC
AMPLICON SIZE: 106 bp
COPY NUMBER/GENOME: 1
FORWARD PRIMER: AADC Exon Jxn F 5' CTGGCCTGATTCTTTCTTTATGG 3'
REVERSE PRIMER: AADC Exon Jxn R 5' CATATGCTTCCTTGTTGCAGATAGG 3'
PROBE: AADC Exon Jxn P 5' GGGACCACAACATGCTGCTCCTTTGAC 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

AADC Exon Junction Forward Primer BDP PN 31120
AADC Exon Junction Reverse Primer BDP PN 31121
AADC Exon Junction Probe BDP PN 31122
RECOMMENDED DNaseI ENZYME: DNaseI, RNase Free, Thermo Scientific (BDP 30695)
RECOMMENDED DNaseI METHOD: Combine 150 µL aliquot of sample with 20 µL 1xPBS, 20 µL 10x Reaction Buffer and 10 µL of DNaseI, RNase-free. Incubate at 37°C for 30 minutes, then 65°C for 10 minutes (to inactivate the DNaseI). Perform prior to QIAGEN extraction.

RECOMMENDED EXTRACTION KIT: QIAGEN Blood Mini Kit (BDP 30444)

CONTROL GENOMIC DNA: pAAV2-CMV-hAADC-KanR

REAGENT	VOLUME/REACTION (µL)
Universal Mastermix (BDP 30268)	12.5 (1x)
NFW or WFI Water (various)	6.25
Primers	0.5 each (10 pmol each)
Probe	0.25 (5 pmol)
Reaction Mix Volume	20
Sample Volume / rxn	5
Total Reaction Volume	25

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.05 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 52

AMPLICON: Human CD33 CAR Amplicon from CD33.28CAR Vector

AMPLICON SEQUENCE:

GCTCCGACATTCAGATGACCCAGAGCCCTAGCAGCCTGAGCGCTTCCGTGGGAGACAGGGTGACCATC
 ACATGCAGGGCCTCCGAGAGCGTGGACAATTACGGCATCAGCTTCATGAACT

AMPLICON SIZE: 120 bases
COPY NUMBER/VIRAL VECTOR GENOME: 1
FORWARD PRIMER: 5' CGA CAT TCA GAT GAC CCA GAG 3'
REVERSE PRIMER: 5' CAT GAA GCT GAT GCC GTA ATT G 3'

PROBE: 5'/5HEX/AGACAGGGT/ZEN/GACCATCACATGCAG/Quencher 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

CD33Fc dPCR Primers/probe set (VIC/HEX) 20X BDP PN 31224

RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP PN 30450)

RECOMMENDED DNaseI ENZYME: NEB (BDP PN 31243)

RECOMMENDED DNaseI METHOD: Elute RNA from Viral RNA Mini Kit in 60.0 µL of Buffer AVE. Combine 39.0 µL 1X PBS, 5.0 µL eluted RNA, 5.0 µL 10x DNase I Buffer and 1.0 µL of DNaseI. Incubate at 37°C for 15 minutes. Add 4.0 µL of RNase-free water and 6.0 µL of 50 mM EDTA and incubate at 75°C for 10 minutes (to inactivate the DNaseI).

CONTROL GENOMIC DNA: CD33Fc RNA Control (BDP PN 31234)

REAGENT	VOLUME/REACTION (µL)
DEPC-Treated Water or RNase-Free Water	16.0
2X Reaction Mix (BDP 30265)	25.0 (1X Final Concentration)
CD33Fc-VIC	1.0 (200 nM Final Concentration)
ROX Reference Dye	1.0
RNaseOUT	1.0
SuperScript III RT/Platinum Taq Mix	1.0
Reaction Mix Volume	45.0
Sample Volume / rxn	5.0
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	01:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 53

AMPLICON: GFP

AMPLICON SEQUENCE:

AGGACGACGGCAACTACAAGACCCGCGCCGAGGTGAAGTTCGAGGGCGACACCCTGGTGAACCGCATC
 GAGCTGAAGGGCATCGACTTCAAGGAGGACGGCAACATCCTGGGGCACAAGCTGGAGTACAACACTACAA
 CAGCCACAACGTCTATATCATGGCCGACAAGCAGAA

AMPLICON REGION: GFP gene
AMPLICON SIZE: 172 bp
COPY NUMBER / GENOME: 1
FORWARD PRIMER: BDP GFP Forward 5' AGGACGACGGCAACTACAAG 3'
REVERSE PRIMER: BDP GFP Reverse 5' TTCTGCTTGTGGCCATGAT 3'
PROBE: BDP GFP Probe 5' [6-FAM] AGGACGGCAACATCCTGGGGCACAA [Quencher] 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

BDP GFP Forward BDP PN 31270
BDP GFP Reverse BDP PN 31271
BDP GFP Probe BDP PN 31272

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP 30442) or QIAGEN Blood Mini Kit (BDP 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP 31236)

CONTROL GENOMIC DNA: GFP Gene Block DNA Control (BDP PN 31273)

REAGENT	VOLUME / REACTION (μL)
Universal Mastermix (BDP 30268)	15 (1x)
NFW or WFI Water (various)	3.5
Primers	0.6 each (5 pmol each)
Probe	0.3 (5 pmol)
Reaction Mix Volume	20
Sample Volume / rxn	10
Total Reaction Volume	30

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50 C	2:00	N/A
2	95 C	10:00	N/A
3	95 C 60 C	:15 1:00	≥ 40

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 54 **AMPLICON:** Human GD2 CAR Amplicon

AMPLICON SEQUENCE:

CTCTGTGATGATCTCCTGCAAGGCCAGCGGCAGCTCCTTCACCGGCTACAACATGAACTGGGTGCGCCA
 GAACATCGGCAAGAGCCTGGAATGGATCG

AMPLICON SIZE: 98 bases
COPY NUMBER/VIRAL VECTOR GENOME: 1
FORWARD PRIMER: GD2 Forward Primer 5' CTCTGTGATGATCTCCTGCAA 3'
REVERSE PRIMER: GD2 Reverse Primer 5' CGATCCATTCCAGGCTCTT 3'
PROBE: GD2 Probe 5'/5HEX/TCATGTTGT/ZEN/AGCCGGTGAAGGAGC/3IABkFQ 3'

REAGENTS (and associated BDP and/or manufacturer part numbers):

GD2 CAR Primer/Probe set (VIC/HEX) 20X BDP PN 31260
RECOMMENDED EXTRACTION KIT: QIAGEN Viral RNA Mini Kit (BDP PN 30450)
RECOMMENDED DNaseI ENZYME: NEB (BDP PN 31243)

RECOMMENDED DNaseI METHOD: Elute RNA from Viral RNA Mini Kit in 60.0 µL of Buffer AVE. Combine 39.0 µL 1X PBS, 5.0 µL eluted RNA, 5.0 µL 10x DNase I Buffer and 1.0 µL of DNaseI. Incubate at 37°C for 15 minutes. Add 4.0 µL of RNase-free water and 6.0 µL of 50 mM EDTA and incubate at 75°C for 10 minutes (to inactivate the DNaseI).

CONTROL GENOMIC DNA: GD2 RNA Control

REAGENT	VOLUME/REACTION (µL)
DEPC-Treated Water or RNase-Free Water	16.0
2X Reaction Mix (BDP 30265)	25.0 (1X Final Concentration)
GD2CAR-VIC	1.0 (200 nM Final Concentration)
ROX Reference Dye	1.0
RNaseOUT	1.0
SuperScript III RT/Platinum Taq Mix	1.0
Reaction Mix Volume	45.0
Sample Volume / rxn	5.0
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	01:00	

RESULTS ANALYSIS:

Threshold: 0.2 ΔRn

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 55 **AMPLICON:** SV40 Large T Antigen

AMPLICON SEQUENCE:

ATAGTGGCTGGGCTGTTCTTTTTAATACATTTTAAACACATTTCAAACACTGTAAGTAAATTCCAAGTACAT
 CCCAAGCAATAACAACACATCATCACATTTTGTTCATTGCATACTCTGTTACAAGCTTCCAGGACACTT
 GTTTAGTTTCTCTGCTTCTTCTGGATTAAAATCATGCTCCTTTAACCCACCTG

AMPLICON REGION: SV40 Large T Antigen
AMPLICON SIZE: 197 bases
SV40 TAg Primer/Probe (FAM) set BDP PN 31313
FORWARD PRIMER: 5' ATAGTGGCTGGGCTGTTCTTT 3'
REVERSE PRIMER: 5' CAGGTGGGTTAAAGGAGCATGA 3'
PROBE: 5' [FAM] TGA TGA TGT GTT GTT ATT GCT TGG GAT GTA CTT GG [Quencher] 3'
RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP PN 30442) or QIAGEN Blood Mini Kit (BDP PN 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP PN 31236)

CONTROL DNA: SV40 TAg Positive Control Gene Fragment (BDP PN 31314)

REAGENT	VOLUME/REACTION (µL)
Universal Master Mix (BDP 30268)	25.0 (1X Final Concentration)
Water	10
SV40 TAg Primer/Probe Set	2.5
Reaction Mix Volume	37.5
Sample Volume / rxn	12.5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	00:15	40
	60°C	01:00	

RESULTS ANALYSIS:

Threshold: 0.05 ΔR or Automatic

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 56 **Residual plasmid DNA assay**

AMPLICON: b-Lactamase (Ampicillin resistance gene)

AMPLICON SEQUENCE:

5'-TTGTTGCCGGGAAGCTAGAGTAAGTAGTTCGCCAGTTAATAGTTTGCGCAACGTTGTTGCC
 ATTGCTACAGGCATCGTGGTGTACGCTCGTCGTTTGGTATGGCTTCATTCAGCTCCGGTTCC
 CAACGATCAAGGCGAGTTACATGATCCCCCATGTTGTGCAA-3'

AMPLICON REGION: Ampicillin resistance gene
AMPLICON SIZE: 165 bases
COPY NUMBER/GENOME: 1

b-Lactamase Primer/Probe (FAM) set BDP PN 31357
FORWARD PRIMER: 5' - TTGTTGCCGGGAAGCTAGAG -3'

REVERSE PRIMER: 5' - TTGCACAACATGGGGGATCA - 3'

PROBE: 5'- /56- FAM/TCCGGTTCC /ZEN/ CAACGATCAAGGCGA/3IIABkFQ/ -3'

RECOMMENDED EXTRACTION KIT: MagNA Pure 24 Total NA Isolation Kit (BDP PN 31236)
CONTROL DNA: bLactam_control (BDP PN 31368)

REAGENT	VOLUME/REACTION (µL)
Universal Master Mix (BDP 30268)	25 (1X)
Water	10
20x primer/probe set	2.5
Reaction Mix Volume	37.5
Sample Volume / rxn	12.5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	00:15	40
	60°C	01:00	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (bLactam_control) 184bp (BDP PN 31368)
CAGTCTATTAATTGTTGCCGGGAAGCTAGAGTAAGTAGTTCGCCAGTTAATAGTTTGCGCAACGTTGTTG
CCATTGCTACAGGCATCGTGGTGTACGCTCGTCGTTTGGTATGGCTTCATTCAGCTCCGGTTCCCAAC
GATCAAGGCGAGTTACATGATCCCCCATGTTGTGCAAAAAGCGG

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 57

Virus genome copy assay (RTqPCR)

AMPLICON: hYP218_VH (hYP218 CAR lentivirus VH_CDR3)

AMPLICON SEQUENCE: (114 bp)

5'-GCCGCGATAACTCCAAGAACACGCTCTACTTGCAAATGAATAGCCTTAGGGCAGAGGACAC
 AGCAGTATACTACTGCGCGCGGAGTACGGCGAACACCAGGTCCACATACTATC-3'

AMPLICON REGION: hYP218_VH (CAR lentivirus)

AMPLICON SIZE: 114 bases

COPY NUMBER/GENOME: 1

YP218_VH Primer/Probe (HEX(VIC)) set BDP PN 31342

FORWARD PRIMER: 5' – GCCGCGATAACTCCAAGAA-3'

REVERSE PRIMER: 5'- GATAGTATGTGGACCTGGTGTTC -3'

PROBE: 5'- /5HEX/TTAGGGCAG/Zen/AGGACACAGCAGTAT/3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: YP218_RNA control (BDP 31344)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (YP218_RNA control, BDP PN 31344) 120 bp

TTAGCCGCGATAACTCCAAGAACACGCTCTACTTGCAAATGAATAGCCTTAGGGCAGAGGACACAGCAG
 TATACTACTGCGCGCGGAGTACGGCGAACACCAGGTCCACATACTATCTCA

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 58

Virus genome copy assay (RTqPCR)

AMPLICON: STEAP1_VH (STEAP1 CAR lentivirus VH_CDR3)

AMPLICON SEQUENCE: (92bp)

5'- CAGTCTGAGAGCGGAAGATACAGCAGTGTATTATTGTGCGCGCGAAAGGAACTATGATTA TGACGAC TATTATTATGCCATGGATTACTGGG-3'

AMPLICON REGION: STEAP1_VH (CAR lentivirus)

AMPLICON SIZE: 92 bases

COPY NUMBER/GENOME: 1

STEAP1_VH Primer/Probe (HEX(VIC)) set BDP PN 31350

FORWARD PRIMER: 5' – CAGTCTGAGAGCGGAAGATAC-3'

REVERSE PRIMER: 5'- CCCAGTAATCCATGGCATAATAATAG -3'

PROBE: 5'- /5HEX/TTGCGCGCGC/Zen/ACAATAATACTGC/3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: STEAP1_RNA control(BDP PN 31366)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (STEAP1_RNA control, BDP PN 31366) 120 bp

GTACCTCCAGATGAACAGTCTGAGAGCGGAAGATACAGCAGTGTATTATTGTGCGCGCGAAAGGAACTA TGATTATGACGACTATTATTATGCCATGGATTACTGGGGTCAGGGCACCCCT

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 59

Virus genome copy assay (RTqPCR)

AMPLICON: CD123 VL (CD123xCD3 BiTE retrovirus VL)

AMPLICON SEQUENCE: (91bp)

5'-CGAGGACTTTGCCATGTACTATTGCCAGCAGCACACAACAAGTACCCTTACACCTTCGGCGGAGGCA
 CCAAGCTGGAAATCAAGAGCGGAGGG -3'

AMPLICON REGION: CD123 VL (CD123xCD3 BiTE retrovirus)

AMPLICON SIZE: 91 bases

COPY NUMBER/GENOME: 1

CD123_VL Primer/Probe (HEX(VIC)) set BDP PN 31369

FORWARD PRIMER: 5' – CGAGGACTTTGCCATGTACTATT-3'

REVERSE PRIMER: 5'- CCCTCCGCTCTTGATTCC -3'

PROBE: 5'- /5HEX/AACAAGTAC/ZEN/CCTTACACCTTCGGCG/3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: CD123_VL_RNA control (BDP PN 31370)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (CD123_VL_RNA control, BDP PN 31370) 120 bp
 TCAGCTCCCTGGAACCCGAGGACTTTGCCATGTACTATTGCCAGCAGCACACAACAAGTACCCTTACACCTT
 CGGCGGAGGCACCAAGCTGGAAATCAAGAGCGGAGGGGGCGGATCCGACA

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 60

Virus genome copy assay (RTqPCR)

AMPLICON: GPC2Re_VL (GPC2 retrovirus CDR VL)

AMPLICON SEQUENCE: (88bp)

5'- CCTCCAGCCTGAAGATGTTGCAACATATTACTGTCAACAGTATGATAATCTCCCGATCACCT
 TCGGCCAAGGGACCAAGCTGGAAATC -3'

AMPLICON REGION: GPC2 Retrovirus CDR3 VL
AMPLICON SIZE: 88 bases
COPY NUMBER/GENOME: 1
GPC2Re_VL Primer/Probe (HEX(VIC)) set BDP PN 31355

FORWARD PRIMER: 5' - CCTCCAGCCTGAAGATGTTG -3'
REVERSE PRIMER: 5' - GATTTCCAGCTTGGTCCCTT - 3'
PROBE: 5' - /5HEX/ATAATCTCC/ZEN/ CGATCACCTTCGGCC/ 3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: GPC2Re_VL_RNA control (BDP PN 31356)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (GPC2Re_VL_RNA control, BDP PN 31356) 120 bp
 CTTTCACCATCAGCAGCCTCCAGCCTGAAGATGTTGCAACATATTACTGTCAACAGTATGATAATCTCCC
 GATCACCTTCGGCCAAGGGACCAAGCTGGAAATCAAACGTGGAGGTGGCG

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 61

Virus genome copy assay (RTqPCR)

AMPLICON: GPC2Lenti_VL (GPC2Lenti CAR lentivirus VL_CDR3)

AMPLICON SEQUENCE: (102bp)

5'- TTCTCTCACAATCAGCAGCATGGAGGGTGAAGATGCTGCCACTTATTACTGCCAGCAGTTTT
 CTAGTTCCCCATCCACGTTCCGGTACTGGGACCAAGCTGGA-3'

AMPLICON REGION: GPC2Lenti_VL (CAR lentivirus)
AMPLICON SIZE: 102 bases
COPY NUMBER/GENOME: 1
GPC2Lenti_VL Primer/Probe (HEX(VIC)) set BDP PN 31378

FORWARD PRIMER: 5' – TTCTCTCACAATCAGCAGCAT-3'
REVERSE PRIMER: 5'- TCCAGCTTGGTCCCAGTA -3'
PROBE: 5'- /5HEX/TAAGTGGCA/ZEN/GCATCTTCACCCTCC/3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: GPC2Lenti_VL_RNA control

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (GPC2Lenti_RNA control, BDP P/N 31379) 120 bp
 ACTCTTATTCTCTCACAATCAGCAGCATGGAGGGTGAAGATGCTGCCACTTATTACTGCCAGCAGTTTTC
 TAGTCCCCATCCACGTTCCGGTACTGGGACCAAGCTGGAGCTGAAAAC TA

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 62

Virus genome copy assay (RTqPCR)

AMPLICON: mTRAC (CD22TCR retrovirus)

AMPLICON SEQUENCE: (83 bp)

5'-
 ACAAGTGCCTGCTGGATATGAAGGCTATGGACAGCAAGTCCAACGGCGCCATCGCCTGGTCTAATCAGA
 CCAGCTTCACATGC -3'

AMPLICON REGION: CD22TCR mTRAC (retrovirus)
AMPLICON SIZE: 83 bases
COPY NUMBER/GENOME: 1
mTRAC Primer/Probe (HEX(VIC)) set BDP PN 31459

FORWARD PRIMER: 5' – ACAAGTGCCTGCTGGATATG -3'
REVERSE PRIMER: 5'- GCATGTGAAGCTGGTCTGATTA -3'
PROBE: 5' -/5HEX/ AGGCTATGG/ZEN/ACAGCAAGTCCAACG /3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: mTRAC_RNA control (BDP 31460)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (mTRAC_RNA control, BDP P/N 31460) 120 bp

CTGGCACCTTTATCACAGACAAGTGCGTGGATATGAAGGCTATGGACAGCAAGTCCAACGGCGC
 CATCGCCTGGTCTAATCAGACCAGCTTCACATGCCAGGATATCTTTAAGGAGA

BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 63

Virus genome copy assay (RTqPCR)

AMPLICON: BCMA_VL (PRODH2 CAR lentivirus, targeting BCMA VL CDR3)

AMPLICON SEQUENCE: (81 bp)

5'- AGGACTTCGCCACCTACTACTGCCAGCAGTACAGGAAGCTCCCCTGGACTTTCGGCCAGG
 GCACCAAACCTGGAGATCAAGC -3'

AMPLICON REGION: BCMA_VL
AMPLICON SIZE: 81 bases
COPY NUMBER/GENOME: 1
BCMA_VL Primer/Probe (HEX(VIC)) set BDP PN 31462

FORWARD PRIMER: 5' – AGGACTTCGCCACCTACTA -3'
REVERSE PRIMER: 5'- GCTTGATCTCCAGTTTGGTG -3'
PROBE: 5'- /5HEX/ TGCCAGCAG/ZEN/TACAGGAAGCTCC /3IABkFQ/-3'

RECOMMENDED EXTRACTION KIT: QIAmp Viral RNA kit (BDP 30450)

CONTROL RNA: BCMA VL RNA control (BDP 31463)

REAGENT	VOLUME/REACTION (µL)
One Step RTqPCR kit (30265) 2x mix	25 (1X)
Water	15.4
20x primer/probe set	2.5
ROX reference dye	0.1
RNAse out (31244)	1
SuperScript/Taq mix	1
Reaction Mix Volume	45
Sample Volume / rxn	5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	15:00	N/A
2	95°C	2:00	N/A
3	95°C	00:15	40
	60°C	00:30	

RESULTS ANALYSIS:

Threshold: 0.2 or automatic

Positive fragment (BCMA_VL_RNA control, BDP 31463) 119 bp
CCAGCCTGCAGCCCGAGGACTTCGCCACCTACTACTGCCAGCAGTACAGGAAGCTCCCCTGGACTTTC
GGCCAGGGCACCAAACCTGGAGATCAAGCGTGGTGGAGGAGGTAGCGGAGGA



BIOPHARMACEUTICAL DEVELOPMENT PROGRAM

SOP Title: Quantitative PCR (qPCR) Methods for Detection and Quantification of Nucleic Acids
SOP Number: 22195
Revision: 09

Attachment 64

AMPLICON: SV40 gp6 (Large T Antigen)

AMPLICON SEQUENCE:

CGTTCAGGCAATGCTTTAAATAATCTTTGGGCCTAAAATCTATTTGTTTTACAAATCTGGCCTGCAGTGT
 TTTAGGCACACTGTACTCATTTCATGGTGACTATTCCAGGG

AMPLICON REGION: SV40 gp6 gene (Large T Antigen)

AMPLICON SIZE: 111 bases

SV40 gp6 Primer/Probe (FAM) set BDP PN 31489

FORWARD PRIMER: 5' CGTTCAGGCAATGCTTTAAATA 3'

REVERSE PRIMER: 5' CCCTGGAATAGTCACCATGAAT 3'

PROBE: 5' [FAM] TTACAAATC/ZEN/TGGCCTGCAGTGTTT
 [Quencher] 3'

RECOMMENDED EXTRACTION KIT: QIAGEN DNA Mini Kit (BDP PN 30442) or QIAGEN Blood Mini Kit (BDP PN 30444) or MagNA Pure 24 Total NA Isolation Kit (BDP PN 31236)

CONTROL DNA: SV40gp6 DNA Control (BDP PN 31490)

REAGENT	VOLUME/REACTION (µL)
Universal Master Mix (BDP 30268)	25.0 (1X Final Concentration)
Water	10
SV40 gp6 Primer/Probe Set	2.5
Reaction Mix Volume	37.5
Sample Volume / rxn	12.5
Total Reaction Volume	50.0

REACTION CONDITIONS:

STAGE	TEMP	TIME	REPEATS
1	50°C	2:00	N/A
2	95°C	10:00	N/A
3	95°C	00:15	40
	60°C	01:00	

RESULTS ANALYSIS:

Threshold: 0.05 ΔR or Automatic

Positive fragment (SV40gp6 DNA control, BDP PN 31490)

CTCACTGCGTTCAGGCAATGCTTTAAATAATCTTTGGGCCTAAAATCTATTTGTTTTACAAATCTGGCCT
 GCAGTGTTTTAGGCACACTGTACTCATTTCATGGTGACTATTCCAGGGGAAATATTTGAGTTC