

# Protein Science - Generation of Recombinant Bacmid DNAs

## Purpose

To generate bacmid DNAs containing the gene you wish to overexpress in insect cells.

## Scope

This protocol details the construction of bacmid DNAs where a baculovirus-expression-system-style expression clone is transformed into an *E. coli* strain containing a baculovirus shuttle vector and transposases on a helper plasmid. Once the transposition is confirmed, bacmid DNA is extracted to be used for propagating recombinant baculovirus.

## Definitions

**PCR:** Polymerase Chain Reaction

## Materials and Equipment

- Oligonucleotides, 5  $\mu$ M (details below)
- Phusion High-Fidelity Master Mix (New England Biolabs, cat. #M0531S)
- Baculoviral expression clone (generated in SOP Combinatorial Assembly of Clone Libraries Using Site-Specific Recombination) or similar pFastBac-style plasmid containing your gene of interest
- Chemically competent cells of an *E. coli* bacmid strain, such as DH10Bac, DE95, etc.
- Lysogeny broth (LB) media (Fisher Scientific, cat. #BP1427-500)
- LB agar plates with 50  $\mu$ g/mL kanamycin, 10  $\mu$ g/mL tetracycline, 7  $\mu$ g/mL gentamycin, 40  $\mu$ g/mL IPTG, and 100  $\mu$ g/mL Bluo-gal (Teknova, SKU #L1919)
- 1.5 mL tubes (USA Scientific, item #1615-5500)
- 2.0 mL tubes (USA Scientific, item #1620-2700)
- Deionized water
- 10 $\times$  TE (Fisher BioReagents, cat. #BP2477-1) (diluted to 1 $\times$  with water, filter-sterilized)
- Ethanol
- Plasmid Prep Buffer A (25 mM Tris-HCl [pH 8.0], 10 mM EDTA, 0.9% (w/v) D-glucose)
  - Store at room temperature.
- Plasmid Prep Buffer B (0.20 M NaOH, 1.33% (w/v) sodium dodecyl sulfate [SDS])
  - Store at room temperature. SDS may precipitate at low temperatures.
  - If precipitation has occurred, redissolve before use by heating at 37°C.
  - Make a fresh batch after 6 months and be sure to cap immediately after use.
- Plasmid Prep Buffer C (7.5 M ammonium acetate)
  - Store at 4°C.
- Bacmid Resuspension Buffer (10 mM Tris-HCl [pH 8.0], 0.1 mM EDTA, 0.1  $\mu$ g/mL ribonuclease A [Sigma, cat. #12091021])
  - Store at room temperature for up to 6 months.
- 0.2 mL PCR tubes (USA Scientific, item #1402-8100)
- Falcon 2059 culture tubes (Fisher Scientific, item #14-959-11B)
- 10 $\times$  gel loading dye (100 mM Tris-HCl [pH 7.5], 100 mM EDTA [pH 8.0], 50 mM NaCl, 0.3% xylene cyanol [Sigma, SKU #X4126], 3% tartrazine [Sigma, SKU #T0388], 50% glycerol) (diluted to 2 $\times$  with deionized water)
- 1 kb Plus DNA Ladder (New England Biolabs, cat. #N3200L) (100  $\mu$ L DNA ladder, 20  $\mu$ L 10 $\times$  gel loading dye, 350  $\mu$ L 1 $\times$  TE, 5  $\mu$ L 0.5M EDTA [pH 8.0])
- Agarose gel 0.8% with ethidium bromide (Embi Tec, cat. #GE-3701)
- 10 $\times$  TAE (Embi Tec, cat. #EC-1016) (diluted to 1 $\times$  with water)
- Glycerol (Sigma, SKU #G5516) (diluted to 60% with water, filter-sterilized)
- 37°C incubator with shaking platform
- Benchtop centrifuge
- Thermocycler

- Electrophoresis equipment
- Gel imager with ultraviolet light and camera
- Spectrophotometer

## Safety Precautions

Ethidium bromide is a mutagen and should be handled with great care. Liquid waste as well as contaminated tips and other disposables must be disposed of properly in accordance with your organization's guidelines.

## Procedures

### A. Oligonucleotide Design

To confirm the transposition event, Thermo Fisher recommends using pUC/M13 forward and reverse primers to amplify across the entire insert from a pFastBac plasmid. We prefer to use four primers (Baccheck-R, Baccheck-F, BaccheckBtB-L, BaccheckBtB-R) to multiplex the two transposition junctions. Oligonucleotides can be synthesized at whatever concentration you prefer and then diluted to 5  $\mu$ M with 1 $\times$  TE. For the multiplex junction PCR, the oligonucleotides can be pooled at equal volumes.

Primer Name	Primer Sequence 5'-3'
pUC/M13 Forward	CCCAGTCACGACGTTGTAAAACG
pUC/M13 Reverse	AGCGGATAACAATTTACACACAGG
Baccheck-R	GTGCTGCAAGGCGATTAAGT
Baccheck-F	TGTGGAATTGTGAGCGGATA
Baccheck BtB-L	ATCAGCCGGACTCCGATTA
Baccheck BtB-R	CCCACACCTCCCCCTGAACCTG

### B. Transformation of Baculoviral Expression Clone into Bacmid Strain

1. Thaw on ice 50  $\mu$ L of DE95, DH10Bac, or a similar bacmid strain.
2. Add 1  $\mu$ L (100–200 ng) of the baculoviral expression clone to the thawed cells.
3. Incubate the mixture on ice for 15 minutes.
4. Heat-shock the cells at 42°C for 45 seconds and immediately place them back on ice for 2 minutes.
5. Add 450  $\mu$ L of LB and place all transformants into a Falcon 2059 tube.
6. Shake the tube for 4 hours at 37°C.
7. Add 48  $\mu$ L LB and 2  $\mu$ L of the transformants to an LB agar plate with 50  $\mu$ g/mL kanamycin, 10  $\mu$ g/mL tetracycline, 7  $\mu$ g/mL gentamycin,

40  $\mu$ g/mL IPTG, 100  $\mu$ g/mL Bluo-gal.

8. Grow the cells overnight at 37°C.
9. Pick two or three white colonies to Falcon 2059 tubes containing 3 mL LB media with 50  $\mu$ g/mL kanamycin and 7  $\mu$ g/mL gentamycin.
10. Grow the colonies overnight at 37°C, 250 rpm.

### C. Preparation of Bacmid DNAs

1. Pellet 2 mL of overnight culture in a benchtop centrifuge at max speed for 1 minute.
2. Aspirate the supernatant.
3. Add 250  $\mu$ L of Plasmid Prep Buffer A to the tube and resuspend the pellet by vortexing.
4. Add 250  $\mu$ L of Plasmid Prep Buffer B and resuspend the pellet by inverting the tube 6–10 times. The mixture should become viscous.
5. Add 250  $\mu$ L of Plasmid Prep Buffer C and resuspend the pellet by inverting the tube 6–10 times. The mixture should have visible white flocculation of cell debris.
6. Spin the tube at max speed in a benchtop centrifuge for 10 minutes.
7. Remove 700  $\mu$ L of supernatant to a fresh 1.5 mL tube.
8. Spin the tube at max speed in a benchtop centrifuge for 10 minutes.
9. Remove 600  $\mu$ L of supernatant to a fresh 2 mL tube.
10. Add 1.2 mL ethanol to the tube and invert the tube 6–10 times to mix.
11. Spin the tube at max speed in a benchtop centrifuge for 20 minutes.
12. Pour off the supernatant.
13. Spin the tube at max speed in a benchtop centrifuge for 1 minute.
14. Aspirate the remaining alcohol mixture, avoiding the white pellet.
15. Allow the tube to dry for 10 minutes at room temperature.
16. Add 100  $\mu$ L of Bacmid Resuspension Buffer and place the culture in a 37°C shaker for 1 hour.

### D. Screening for Correct Transposition

1. Add the following to a 0.2 mL PCR tube:

- a. 1  $\mu\text{L}$  of bacmid DNA prep
  - b. 1  $\mu\text{L}$  of 5  $\mu\text{M}$  M13 forward and reverse primers, or pooled junction oligo mix
  - c. 8  $\mu\text{L}$  deionized water
  - d. 10  $\mu\text{L}$  Phusion High-Fidelity Master Mix
2. Add the tube to a thermocycler with the following conditions: initial denaturation at 98°C for 45 seconds, followed by 30 cycles (each 98°C for 30 seconds, 55°C for 30 seconds, and 72°C for 30 seconds—or 30 seconds per kb if amplifying a full insert with M13 primers), followed by a 10-minute final elongation at 72°C and cooling to 4°C.
  3. Set up electrophoresis equipment and fill it with 1 $\times$  TAE. Insert 0.8% agarose gel into the dock.
  4. In a fresh 0.2 mL PCR tube, mix 5  $\mu\text{L}$  of the PCR product with 5  $\mu\text{L}$  of 2 $\times$  loading dye and load the full volume onto a 0.8% agarose gel. Load 1  $\mu\text{L}$  of the 1 kb Plus DNA Ladder with loading dye in a neighboring well.
  5. Electrophorese for 1 hour at 100V.
  6. Image the gel with ultraviolet light. The expected amplicons in the multiplex reaction will be 570 bp/840 bp. The expected amplicon in the M13 reaction will be the entire expression clone insert (between Tn7 left and right arms).
  7. If the amplicons are correct, make a 15% glycerol stock of the correct construct and store it at  $-80^{\circ}\text{C}$ .
  8. The prepared DNA contains *E. coli* genomic DNA, residual expression and helper plasmids, and the bacmid DNA, so concentration for the bacmid DNA cannot be accurately determined.
  9. Reserve this DNA at  $-20^{\circ}\text{C}$  until you are ready to prepare the baculovirus.

## References

1. Luckow VA, Lee SC, Barry GF, and Olins PO (1993). Efficient generation of infectious recombinant baculoviruses by site-specific transposon-mediated insertion of foreign genes into a baculovirus genome propagated in *Escherichia coli*. *J Virol.* 67:4566–4579. doi: 10.1128/JVI.67.8.4566-4579.1993.
2. Gillette W, Frank P, Perkins S, Drew M, Grose C, Esposito D (2019). Production of Farnesylated and Methylated Proteins in an Engineered Insect Cell System. *Methods Mol Biol.* 2009:259-277. doi:10.1007/978-1-4939-9532-5\_20.